
INDICATIONS FOR USE

The F&P Optiflow™+ is a nasal cannula patient interface for use with specified respiratory gas humidifiers to treat spontaneously breathing pediatric (3 years and older) to adult patients who would benefit from receiving high flow warmed and humidified respiratory gases to the upper airway.

This device is designed to be used in a hospital, sub-acute facility or long-term care facility by appropriately qualified healthcare professionals.

WARNINGS

- Appropriate patient monitoring (e.g., oxygen saturation) must be used at all times. Failure to monitor the patient (e.g., in the event of an interruption to gas flow) may result in loss of therapy, serious injury or death.
- This device is not intended for life support. Do not use on patients who cannot tolerate a brief interruption of therapy, to avoid serious injury or death.
- Nasal delivery of respiratory gases generates flow-dependent positive airway pressure (PAP). This must be taken into account where PAP could have adverse effects on a patient.
- Failure to set up and use this device as per these instructions can compromise performance and affect safety. The device has only been validated for use with the equipment described in these instructions.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

Exposure to oxygen increases the risk of fire. To avoid risks of burns, patient injury or death:

- Do not use the system near a naked flame or any ignition source, including electrosurgery, electrocautery, laser surgery instruments, or while smoking.
- Only use lotions and/or salves that are labeled as oxygen compatible.
- Do not lubricate connections or tubing.
- Do not leave the interface on bed coverings or cushions if the device is turned on but not in use. Oxygen can make the materials more flammable.

CAUTIONS

- For single patient use only. Reuse may result in transmission of infectious substances.
- Using the cannula may cause skin reddening. Monitor for signs of skin damage or irritation and adjust the cannula if necessary.
- Do not use if packaging label perforations are torn or if packaging is open, to avoid contamination.
- Achievable flow rates may depend on individual patient and/or flow source.
- Before connecting the interface, check for adequate gas flow and ensure that the system has warmed up.

Failure to comply with the following may result in loss of therapy:

- Do not crush or stretch cannula or tubing.
- Do not use if tubing damaged with holes, tears or kinks.
- Do not reprocess (e.g., soak, wash or sterilize) any part of the device as this will cause degradation of materials.
- Do not use with an air entrainer. This will reduce the amount of flow and/or oxygen reaching the patient.
- Check for condensate regularly. Drain condensate away from patient as required, to prevent blockages or condensate from reaching the patient.

- Ensure that patients in an awake prone position (i.e., on their stomach) are not lying on the interface tubing or heated breathing tube. Reposition the breathing tube clip if required. Failure to comply may result in damage to the device.
- Using the device for longer than the specified time can result in degradation of materials or infection.

SIDE EFFECTS

The following side effects may occur during nasal high flow therapy:

- Skin reddening

SPECIFICATIONS

Airvo™ series humidifier with an AirSpiral™ heated breathing tube and chamber kit (e.g., 900PT561).

Flow range: OPT942 10–60* L/min
OPT944 10–70* L/min
OPT946 10–70* L/min

* Maximum flow rate depends on cannula size and humidifier model. Refer to your Airvo series User Manual for more information.

OR

MR850 Respiratory Humidifier in invasive mode, RT series kit with 22 mm heated inspiratory tube and chamber (e.g., RT232).

Flow range: 5–60 L/min

OR

F&P 950™ Respiratory Humidifier in Optiflow mode with F&P 950 Adult Heated Circuit kit (e.g., 950A40).

Flow range: 5–70 L/min

OR

F&P 820 System in Setting 3 or 4 with 22 mm F&P 820 Circuit kit (e.g., 820A10).

Flow range: 5–70 L/min

Note: Flow rates are expressed in body temperature, pressure, saturated (BTPS).

Note: Not all products are available in all countries.

Additional information:

- This product is intended to be used for a maximum of 14 days.
- Operating temperature range: 18 °C to 28 °C (64 °F to 82 °F)
- Permissible transport and storage conditions: -10 °C to 50 °C (14 °F to 122 °F)
- Refer to the Airvo series User Manual for additional specifications and permissible operating conditions.

SET-UP INSTRUCTIONS

- Select appropriate size. Prongs must not create a seal in the nares. A clear gap must be visible around each prong.
- Adjust headstrap to fit. Do not over-tighten.
- a Ensure headstrap clip is attached to the headstrap, to prevent cannula from being pulled out of the nares.
- b The cannula can become unattached if not used with the headstrap clip.
- Attach the breathing tube clip to a secure location to prevent the cannula from pulling on the ears or off the face.
- If using an MR850 Respiratory Humidifier, attach the breathing tube clip to the heated breathing tube but ensure probe cable is not crushed by the tubing clip.

DURING USE

- As appropriate for the patient's condition, monitor for any disruptions to the patient receiving flow.
- To prevent condensate, remove/minimize the impact of anything that may cool the tubing (e.g., a fan, air conditioning).
- Notice to user:** If a serious incident has occurred while using this device please inform your local Fisher & Paykel Healthcare representative and Competent Authority in your country.

DISPOSAL

Dispose of products and packaging according to local guidelines. Facilities should discard used products according to their standard method for contaminated waste. Follow recycling guidelines for product packaging as applicable.

INTELLECTUAL PROPERTY INFORMATION

F&P, Optiflow, Airvo, AirSpiral and 950 are trademarks of Fisher & Paykel Healthcare Limited.

SYMBOL DEFINITIONS

	Catalog number
	Refer to instructions
	Medical device
	Batch code
	Manufacturer
	Date and country of manufacture CC NZ: New Zealand MX: Mexico
	Use-by date
	Storage and transport temperature range
	Single use
	European conformity - TÜV SÜD
	Prescription only
	Authorized European Representative
	Authorized Swiss Representative
	Authorized UK Representative
	Importer
	Distributor
	Do not open with sharp object
	Recommended opening method
	Small nasal cannula size
	Medium nasal cannula size
	Large nasal cannula size
	Correct fitting instruction method
	Incorrect fitting instruction method

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