# Systemes F&P Optiflow™

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**F&P Optiflow** Flow Matters Édition 11

**Optiflow™ matters**

Guideline recommendations for the use of nasal High Flow (NHf), also High Flow Nasal Cannula (HFNC), are supported by analyzed data from research investigating the effect of NHf on clinical outcomes, such as the reduced need for tracheal intubation. When selecting an NHf system, it is important to review the entire system, including design and device limits, to provide the therapy proven to deliver the expected outcomes.

**Summary**

- The National Institute of Health (NIH), the Australia and New Zealand Intensive Care Society (ANZICS) and the Surviving Sepsis Campaign (SSC) recommend NHf for use in COVID-19 related hypoxemia.<sup>1\*</sup>
- These recommendations are supported by findings from the systematic reviews with meta-analysis.<sup>2\*</sup>
- A review conducted by Fisher & Paykel Healthcare (F&P) showed that the flow rates used in the controlled published studies<sup>3\*</sup> analyzed by the first meta-analysis<sup>2\*</sup> ranged from 10 L/min to 60 L/min and 82% of the studies required flows > 45 L/min.
- When this review was repeated on the 52 acute adult NHf controlled studies (with subjects > 18), found using a systematic search of the PubMed database, it was again shown that the flow rates used ranged from 10 L/min to 60 L/min and that 85% of the studies required flows > 45 L/min.
- F&P Optiflow systems (including F&P Optiflow interface) and humidity settings of 37 °C were widely used.

**Guideline recommendations**

Recent guidelines for the clinical management of COVID-19 from organizations such as the NIH, ANZICS and SSC recommend the use of NHf as respiratory support in adults. These recommendations are supported by systematic reviews with meta-analysis which search for, review and analyze clinical data from controlled studies, such as Randomized Controlled Trials (RCTs). F&P conducted a review of the systems and settings used in studies from which analyzed data formed the basis of these recommendations.

\*The NIH is part of the U.S. Department of Health and Human Services, is the USA federal medical research agency.  
\*\*The SSC is a collaboration between the Society of Critical Care Medicine (SCCM) and the European Society of Intensive Care Medicine (ESICM).  
www.fphcare.com

**Reviews with meta-analysis**

NHF recommendations from NIH COVID-19 Treatment Guidelines, the ANZICS COVID-19 guidelines, and the SSC Guideline on the Management of COVID-19 Adults with Coronavirus Disease 2019 (COVID-19) are supported by the following systematic reviews with meta-analysis: Tran et al. 2021, Chu et al. 2021, Ni et al. 2018, Hochway et al. 2019 and Aguilera et al. 2020.<sup>4\*</sup>

**Analyzed published studies**

These five review analyzed data from 22 published studies (mostly RCTs) and one presentation.<sup>5\*</sup> The studies represent various NHf applications, including primary respiratory support, pre-oxygenation prior to intubation, and post intubation respiratory support. The studies reported the NHf systems and settings that were used.

**Systems and settings**

The reported flow rates ranged between 10 L/min and 60 L/min with the majority requiring flow at the higher end of the range.

Of the 22 published and analyzed studies, 20 (91%) used F&P Optiflow systems, including a F&P Optiflow patient interface and a F&P humidity delivery system with humidity setting of 37 °C.

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