








Optiflow™ Nasal High Flow (NHF) therapy delivers respiratory support to your spontaneously breathing patients. It provides heated, humidified air and/or oxygen at flow rates up to 70 L/min through the unique Optiflow patient interfaces.

## MECHANISMS OF ACTION



With Optiflow NHF, you can independently titrate flow and oxygen concentration ( $\text{FiO}_2$  21 – 100%) according to your patient's needs.

The mechanisms of action differ from those of conventional therapies, as do the resulting physiological effects and clinical outcomes.


## PHYSIOLOGICAL EFFECTS

- ↑ **IMPROVES** ventilation and gas exchange
- ↓ **REDUCES** respiratory rate<sup>1-7</sup>
- ↓ **REDUCES** carbon dioxide<sup>8-10</sup>
- ↑ **INCREASES** end-expiratory lung volume<sup>1</sup>
- ↑ **IMPROVES** mucus clearance<sup>11</sup>
- ↑ **IMPROVES** oxygenation<sup>1,2,4,7,12-16</sup>

## CLINICAL OUTCOMES

- ↓ **REDUCES** escalation of care when used:
  - as a first-line respiratory support<sup>14</sup>
  - post-extubation<sup>13,17-20</sup>
- ↓ **REDUCES** mortality rate<sup>14</sup>
- ↑ **IMPROVES** symptomatic relief<sup>2,3,14</sup>
- ↑ **IMPROVES** comfort and patient compliance<sup>2,3,13,17,20</sup>






MECHANISMS OF ACTION

RESPIRATORY SUPPORT



Reduction of dead space

Clearance of expired air in the upper airways<sup>8</sup>  
Reduces rebreathing of gas with high CO<sub>2</sub> and depleted O<sub>2</sub><sup>8</sup>  
Increases alveolar ventilation<sup>8</sup>




Dynamic positive airway pressure

Breath- and flow-dependent airway pressure<sup>9,21</sup>  
Promotes slow and deep breathing<sup>9</sup>  
Increases alveolar ventilation<sup>1,8</sup>



Airway hydration


Optimal Humidity  
Prevents desiccation of the airway epithelium<sup>22</sup>  
Improves mucus clearance<sup>11,22</sup>



Patient comfort

Optimal Humidity  
Open system  
No seal required  
Comfortable<sup>2,13</sup>  
and easy to use  
Patient tolerance<sup>2,14</sup>

Supplemental oxygen when required



Summary of applications for NHF therapy



Medical society

AARC.

Piraino et al. 2021<sup>23</sup>

ACP.

Qaseem et al. 2021<sup>24</sup>

SSC.

Evans et al. 2021<sup>25</sup>

ESICM.

Rochwerger et al. 2020<sup>26</sup>

ERS.

Oczkowski et al. 2021<sup>27</sup>

TSANZ.

Barnett et al. 2022<sup>28</sup>

WHO.

WHO Guideline Development Group 2022<sup>29</sup>

	MEDICAL   SURGICAL Primary support	Pre-escalation support	De-escalation support	Complementary support	Prophylactic support
AARC. Piraino et al. 2021 <sup>23</sup>	● A				●
ACP. Qaseem et al. 2021 <sup>24</sup>	●		●		
SSC. Evans et al. 2021 <sup>25</sup>	●				
ESICM. Rochwerger et al. 2020 <sup>26</sup>	● B	●	● C		
ERS. Oczkowski et al. 2021 <sup>27</sup>	●	●	●	●	
TSANZ. Barnett et al. 2022 <sup>28</sup>	●				●
WHO. WHO Guideline Development Group 2022 <sup>29</sup>	● D				

AARC: American Association for Respiratory Care. ACP: American College of Physicians. SSC: Surviving Sepsis Campaign. ESICM: European Society of Intensive Care Medicine. ERS: European Respiratory Society. TSANZ: Thoracic Society of Australia and New Zealand. WHO: World Health Organization.

A. Hypoxemia and immunocompromised patients with ARF. B. Sepsis induced hypoxemic respiratory failure  
C. Continue to use NHF if already receiving therapy during intubation. D. Acute Hypoxemic patients with severe to critical COVID-19.

# Clinical practice guidelines



## AARC CLINICAL PRACTICE GUIDELINES

Piraino T, et al. Respiratory Care. 2021.<sup>23</sup>



### Primary support – Medical

General recommendations for the delivery of supplemental oxygen for patients who require oxygen  
Aim for SpO<sub>2</sub> range of 94–98% for most of hospitalized patients (included critically ill patients).

Aim for 88–92% for patients with COPD.

Aim for 88–95% for patients with ARDS.

Consider early initiation of NHF.

LEVEL C



### De-escalation support

NHF is preferred to COT immediately post-extubation in patients who require supplemental oxygen.

LEVEL B



### Prophylactic support

NHF is preferred to COT to avoid escalation to NIV or invasive ventilation in patients who require supplemental oxygen.

LEVEL B



### Immuno-compromised

Either NHF or COT may be used in patients who require supplemental oxygen.

LEVEL B

## ACP CLINICAL GUIDELINES

Gaseem A, et al. Annals of Internal Medicine. 2021.<sup>24</sup>



### Primary support – Medical

NHF is preferred to NIV in patients with AHRF.

CONDITIONAL RECOMMENDATION



### De-escalation support

NHF is preferred to COT in patients with post-extubation AHRF.

CONDITIONAL RECOMMENDATION



### Primary support – Medical

NHF is preferred to NIV in patients with sepsis-induced hypoxemic respiratory failure.

WEAK RECOMMENDATION

ENDORSED BY:

- Society of Critical Care Medicine
- American Association of Critical Care Nurses
- American College of Chest Physicians
- American College of Emergency Physicians
- American Thoracic Society

## ESICM CLINICAL PRACTICE GUIDELINES

Rochweg B, et al. Intensive Care Medicine. 2020.<sup>26</sup>



### Primary support – Medical

NHF is preferred to COT for patients with hypoxemic respiratory failure.

STRONG RECOMMENDATION



### Primary support – Surgical

NHF is preferred to COT in high risk and/or obese patients undergoing cardiac or thoracic surgery to prevent respiratory failure in the immediate postoperative period.

Prophylactic NHF to prevent respiratory failure in other postoperative patients is not recommended.

CONDITIONAL RECOMMENDATION



### Pre-escalation support

No recommendation is made regarding use of NHF in the peri-intubation period.

NHF during intubation should be continued for patients who are already receiving NHF.

CONTINUE NHF



### De-escalation support

NHF is preferred to COT following extubation in patients with any high-risk feature who were intubated for > 24 hours.

NIPPV is preferred to NHF in patients who would normally be extubated to NIPPV.

CONDITIONAL RECOMMENDATION

## ERS CLINICAL PRACTICE GUIDELINES

Oczkowski S, et al. European Respiratory Journal. 2021.<sup>27</sup>



### Primary support – Medical

NHF is preferred to COT or NIV in patients with acute hypoxemic respiratory failure.

CONDITIONAL RECOMMENDATION



### Primary support – Medical

Trialling NIV prior to use of NHF in patients with COPD or acute hypercapnic respiratory failure.

CONDITIONAL RECOMMENDATION



### De-escalation support

NHF is preferred to COT in low-risk non-surgical patients.

NIV is preferred to NHF in non-surgical patients at high risk of extubation failure, unless NIV is contraindicated.

CONDITIONAL RECOMMENDATION



### Primary support – Surgical

Either NHF or COT can be used in post-operative patients at low risk of respiratory complications.

Either NHF or NIV can be used in post-operative patients at high risk of respiratory complications.

CONDITIONAL RECOMMENDATION



### Complementary support

NHF is preferred to COT during breaks from NIV in patients with acute hypoxemic respiratory failure.

CONDITIONAL RECOMMENDATION



Frat et al. 2015<sup>14</sup>

*The New England Journal of Medicine*

**High-flow oxygen through nasal cannula in acute hypoxemic respiratory failure.**

**Study**

A 23-center study compared NHF to use of a non-rebreather mask (standard oxygen) and NIV as a primary treatment.

The primary outcome was the number of patients intubated at day 28 (not attained).


**Method**

310 pre-intubation patients in acute hypoxemic respiratory failure ( $\text{PaO}_2:\text{FiO}_2 \leq 300$  mmHg) were randomized to receive NHF, non-rebreather mask or NIV.


**Results**

- **NHF significantly reduced ICU** ( $p = 0.047$ ) **and 90-day mortality** ( $p = 0.02$ )
- The primary outcome was not met for all patients ( $p = 0.18$ ), however, **NHF significantly reduced the need for intubation in more acute patients ( $\text{PaO}_2:\text{FiO}_2 \leq 200$  mmHg)** ( $p = 0.009$ )
- Significant increase in ventilator-free days on NHF ( $p = 0.02$ )
- NHF significantly reduced intensity of respiratory discomfort ( $p < 0.01$ ) and dyspnea ( $p < 0.001$ )

**ICU mortality**



**Reduced intubation rate\***




\*Patients with  $\text{PaO}_2:\text{FiO}_2 \leq 200$  mmHg

Ischaki et al. 2017<sup>30</sup>

*European Respiratory Review*

**Nasal high flow therapy: a novel treatment rather than a more expensive oxygen device.**



MV = mechanical ventilation; SOT = standard oxygen treatment.  
Adapted from original paper (Ischaki et al. Eur Respir Rev. 2017.); used under Creative Commons licence 4.0.  
Please note that this material is intended exclusively for healthcare practitioners and the information conveyed constitutes neither medical advice nor instructions for use. This material should not be used for training purposes or to replace individual hospital policies or practices. Before any product use, consult the appropriate user instructions.



Cortegiani et al. 2020<sup>31</sup>

*Critical Care*

**High flow nasal therapy versus noninvasive ventilation as initial ventilatory strategy in COPD exacerbation: a multicenter non-inferiority randomized trial.**

**Study**

A 9-center RCT compared NHF to NIV as an initial ventilatory strategy in hypercapnic COPD exacerbation.

**Patients**

n = 79, Mild-to-moderate AECOPD ( $\text{pH}$  7.25-7.35,  $\text{PaCO}_2 \geq 55$  mmHg before ventilator support)

**Intervention**

NHF

**Control**

NIV

**Outcome**


Primary:  $\text{PaCO}_2$  from baseline to 2 h (non-inferiority margin 10 mmHg)

Secondary: non-inferiority of NHF to NIV in reducing  $\text{PaCO}_2$  at 6 h rate of treatment changes, dyspnea, discomfort, RR, ABG, hospital LoS, mortality

**Results**

- NHF was non-inferior to NIV in reduction of  $\text{PaCO}_2$
- Both treatments had a significant effect on  $\text{PaCO}_2$  reductions over time, and trends were similar between groups.

**Mean  $\text{PaCO}_2$  reduction from baseline at 2 hours**



Pantazopoulos et al. 2020<sup>32</sup>

*COPD: Journal of Chronic Obstructive Pulmonary Disease*

**Nasal high flow use in COPD patients with hypercapnic respiratory failure: treatment algorithm & review of the literature.**

**Study**

A review of the evidence for NHF use for treatment of stable hypercapnic COPD patients and acute hypercapnic exacerbation of COPD, (21 studies: 9 AECOPD/12 stable COPD studies), with proposed evidence-based algorithm for the dinical application of NHF In patients with AECOPD.

**Conclusions**

It may well also be used in place of NIV in the least tolerant and compliant patients, or in association with NIV to reduce mask-related side effects.


NHF seems to be effective in improving clinical and gas exchange parameters in patients with moderate hypercapnic respiratory failure, with an acceptable rate of non-responders who required additional ventilatory support.

**Results**

NHF recommended for patients with:

- $\text{pH}$ : 7.25 - 7.35
- escalate to NIV if  $\text{pH} < 7.25$

**Algorithm for NHF use in acute hypercapnic exacerbation of COPD**



Adapted from original paper (Pantazopoulos et al. COPD. 2020.); used under Creative Commons licence 4.0.  
Please note that this material is intended exclusively for healthcare practitioners and the information conveyed constitutes neither medical advice nor instructions for use. This material should not be used for training purposes or to replace individual hospital policies or practices. Before any product use, consult the appropriate user instructions.







Hernández et al. (Oct) 2016<sup>18</sup>  
*Journal of the American Medical Association*

**Effect of post-extubation high-flow nasal cannula vs noninvasive ventilation on reintubation and post-extubation respiratory failure in high-risk patients: A randomized clinical trial.**

**Design**

3 center RCT

**Patients**

n = 604, Patients at high risk for reintubation


Intervention	Control
NHF	NIV

**Outcome**

Reintubation and post-extubation respiratory failure within 72 hours

**Results**

- NHF was non-inferior to NIV for preventing reintubation and post-extubation respiratory failure.
- No patients in the NHF group suffered adverse effects requiring withdrawal of the therapy, compared to 42.9% of patients in the NIV group.



Hernández et al. (Apr) 2016<sup>19</sup>  
*Journal of the American Medical Association*

**Effect of post-extubation high-flow nasal cannula vs conventional oxygen therapy on reintubation in low-risk patients.**

**Design**

7 center RCT

**Patients**

n = 527, Patients at low risk for reintubation

**Intervention**

NHF for 24 hrs post extubation


**Control**

COT for 24 hrs post extubation

**Outcome**

Primary: reintubation within 72 hours  
Secondary: post-extubation respiratory failure, adverse events, and time to reintubation, ICU and hospital LoS

**Results**








# Flow usage


Guidelines for the use of NHF for acute respiratory support in adults are supported by peer-reviewed and published evidence.

## What flow rates and ranges are used?

Systematic search of PubMed database for acute adult NHF controlled studies with subjects n > 39.



## Flow rates used in the 52 controlled studies on acute adult NHF (with subjects n >39)\*




\*Systematic search of the PubMed database: Conducted on 17 September 2020 using pre-defined search terms. Filtered using an Excel database and checked by an internal clinical team.

# Physiological indicators for stability

There is an ever-increasing body of clinical literature which may provide guidance on the day-to-day application of Optiflow NHF therapy.

## When are the effects of Optiflow NHF seen?

Sztrymf<sup>4</sup> associated Optiflow NHF therapy with sustained beneficial effects on oxygenation and physiological parameters for patients with acute respiratory failure. Similarly Rittayamai<sup>5</sup> showed significant improvement in post-extubation patients. These studies may provide guidance on patient responses to the therapy.



# Is there a way to predict the outcome of NHF?

The validated ROX index<sup>33</sup> predicts failure in adults with AHRF receiving NHF, at 4 time intervals: 2, 6, 12 and > 12 hours. It's an easy-to-use dynamic bedside tool.

## ROX index: Predicting NHF success and failure

$$\frac{\text{SpO}_2 \div \text{FiO}_2}{\text{Respiratory rate}} = \text{ROX index}$$

'Healthy' example

$$\frac{95 \div 0.21}{15} = 30.2$$

'Patient' example

$$\frac{95 \div 0.85}{37} = 3.0$$

ROX index trend over time is more important than a single measurement.

The trend of FiO<sub>2</sub> required to maintain target SpO<sub>2</sub> (i.e. 95%) and patient respiratory rate directly effect ROX trend.



NHF success

$$\frac{\text{SpO}_2 \div \text{FiO}_2 \downarrow}{\text{Respiratory rate} \downarrow} = \text{ROX index} \uparrow$$

NHF failure

$$\frac{\text{SpO}_2 \div \text{FiO}_2 \uparrow}{\text{Respiratory rate} \uparrow} = \text{ROX index} \downarrow$$

## XY plot between respiratory rate and FiO<sub>2</sub>



## EDUCATIONAL APP



The ROX Vector App proposes a model for considering the trend in ROX values over time.

# Optiflow in practice: IOWA METHODIST MEDICAL CENTER, DES MOINES, IOWA

Jackson et al. 2021<sup>34</sup>  
*Respiratory Care*

## Implementation of high-flow nasal cannula therapy outside the intensive care setting.

### Design

Single center cohort observational study (pre and post NHF implementation)

### Patients

n = 346

### Intervention

18-month after implementing NHF therapy

### Control


Prior to NHF implementation

### Outcome



Share education and implementation process. Report patient outcomes.

### Results

#### Mortality



#### Escalation to MV



After implementation:

- 53% (n = 184) of NHF patients avoided the ICU completely
- 486 ICU days were avoided

## Implementation strategy

1. Protocol of NHF written (for undifferentiated respiratory compromise and increased oxygen requirement).
2. Education of hospital staff: Key groups included respiratory therapists; internal medicine and surgery residents; internal medicine, pulmonology, trauma, cardiology, and emergency medicine physicians; and nurses on all patient floors and in the ED.
3. At least 4 hourly assessment by respiratory therapist.
4. Study team regular review of patient safety and adverse events.

NHF education topics by audience	Physicians and residents	Nurses	Respiratory therapists
Theory and physiology of NHF therapy	●	●	●
NHF protocol	●	●	●
NHF device setup and electronic medical record documentation			●
NHF device maintenance		●	●
De-escalation and weaning	●	●	●



# Nasal high flow therapy in infants and children



The body of literature helps to define the role of NHF in pediatric respiratory care and supports:

- the use of NHF early in the course of respiratory distress is associated with improved physiological outcomes compared with standard oxygen therapy, including:<sup>35-40</sup>
  - improved breathing patterns and rapid unloading of the respiratory muscles
  - significant reduction in the work of breathing
  - rapid improvement to respiratory distress
  - improved mucosal function and secretion clearance through the delivery of heated and humidified gas
- the early use of NHF in bronchiolitis outside of the PICU, either as primary support or early rescue therapy, can lead to reduced escalation of care.<sup>41-45</sup>


## What flow rates are used?

- 2 L/kg/min for patients up to 12 kg in weight has been shown to produce a rapid improvement in respiratory distress, and a reduced need for escalation of therapy.
- Flow rates for those over 12 kg have been protocolized by the PARIS<sup>46</sup> and FIRST-ABC<sup>47</sup> research groups.

Weight (kg)	≤ 12	13 – 15	16 – 30	31 – 50	> 50
Starting flow rate	2 L/min/kg	25 – 30 L/min	35 L/min	40 L/min	50 L/min

## How do I know if my patient is responding to NHF?

Schibler 2011, Intensive Care Med Mayfield 2014, J Pediatric Child Health<sup>45</sup>



An improvement to these indicators within 60 minutes may help to identify infants who are likely to respond to NHF. No improvement may help identify infants who are likely to require escalation of care.<sup>5,11</sup>

Note: Standard O<sub>2</sub> = 100% O<sub>2</sub> NHF at 2 L/kg/min = Total flow/kg/min; FiO<sub>2</sub> titrated

# PEDIATRIC

Franklin et al. 2018<sup>35</sup>  
New England Journal of Medicine



A randomized trial of high-flow oxygen therapy in infants with bronchiolitis.

## Design

17 centered RCT

## Patients

n = 1472, infants < 12 months old with bronchiolitis

## Intervention

NHF starting at 2 L/kg/min

## Comparator


Standard oxygen therapy at < 2 L/min

## Primary outcome

Therapy failure requiring therapy escalation or PICU admission

## Results

Use of NHF as a primary treatment in the ED and general care areas resulted in a significantly lower rate of therapy failure compared with standard oxygen therapy (12 vs. 23%, p < 0.001)



1 in 9 patients experienced therapy failure on NHF 2 L/kg/min



1 in 4 patients experienced therapy failure on standard O<sub>2</sub> < 2 L/min



There were no significant differences between the secondary outcomes (PICU admissions, intubation rates and adverse events).

Ramnarayan et al. 2022<sup>48</sup>  
Journal of the American Medical Association



Effect of high-flow nasal cannula therapy vs. continuous positive airway pressure therapy on liberation from respiratory support in acutely ill children admitted to pediatric critical care units: a randomized clinical trial.

## Design

24 centered non-inferiority RCT

## Patients

n = 573 (0 - 15 years, median age: 9 months), admitted to critical care requiring respiratory support

## Intervention

NHF starting at 2 L/kg/min

## Comparator


CPAP at 7-8 cm H<sub>2</sub>O

## Primary outcome


Time to liberation from respiratory support

## Results

- When used as first-line therapy, NHF met the noninferiority criteria when compared with CPAP for time on respiratory support (NHF: 52.9 hours vs. CPAP: 47.9 hours; adjusted hazard ratio: 1.03 (95% CI: 0.86 - 1.22))
- Therapy failure occurred more frequently in the CPAP group compared with the NHF group




Of the secondary outcomes, the NHF group had significantly:



Cost benefits

Use Optiflow NHF to reduce escalation<sup>14,18</sup> thereby avoiding associated costs.



Using Optiflow NHF as a first-line therapy (both pre-intubation and post-extubation) may reduce a patient's escalation 'up the acuity curve', resulting in better patient outcomes and reduced costs of care.

Apply Airvo early for stabilization and benefit the patient throughout their stay

ED

HIGH		
MEDIUM		
LOW		

- Superiority to COT
- Reduced need for therapy escalation
- Easy communication during assessment
- Physiological markers of stabilization
- ED exit to a lower acuity setting

ICU

HIGH		
MEDIUM		
LOW		

- Superiority to COT
- Reduced need for intubation/re-intubation
- Reduced ICU length of stay\*
- Non-inferiority to NIV\*
- ICU discharge to a lower acuity setting

\* For post extubation resp. support.

GENERAL CARE

HIGH		
MEDIUM		
LOW		

- Superiority to COT
- Continue patient stability outside ICU
- Airway hydration
- Hospital discharge to community

Adjust Airvo settings to suit the patient and environment.

**AIRVO STAYS WITH THE PATIENT** →





SCAN ME

## REFERENCES

\*Clinical Practice Guidelines in blue

- Corley A, Caruana LR, Barnett AG, Tronstad O, Fraser JF. Oxygen delivery through high-flow nasal cannulae increase end-expiratory lung volume and reduce respiratory rate in post-cardiac surgical patients. *Br J Anaesth*. 2011; 107(6):998-1004.
- Roca O, Riera J, Torres F, Masclans JR. High-Flow Oxygen Therapy in Acute Respiratory Failure. *Respir Care*. 2010; 55(4):408-13.
- Lenglet H, Sztrymf B, Leroy C, Brun P, Dreyfuss D, Ricard JD. Humidified High Flow Nasal Oxygen During Respiratory Failure in the Emergency Department: Feasibility and Efficacy. *Respir Care*. 2012; 57(11):1873-8.
- Sztrymf B, Messika J, Bertrand F, Hurel D, Leon R, Dreyfuss D, et al. Beneficial effects of humidified high flow nasal oxygen in critical care patients: a prospective pilot study. *Intensive Care Med*. 2011; 37(11):1780-6.
- Rittayamai N, Tscheikuna J, Rujjivit P. High-Flow Nasal Cannula Versus Conventional Oxygen Therapy After Endotracheal Extubation: A Randomized Crossover Physiologic Study. *Respir Care*. 2014; 59(4): 485-90.
- Roca O, Pérez-Terán P, Masclans JR, Pérez L, Galve E, Evangelista A, et al. Patients with New York Heart Association class III heart failure may benefit with high flow nasal cannula supportive therapy: High flow nasal cannula in heart failure. *J Crit Care*. 2013; 28(5):741-6.
- Peters SG, Holets SR, Gay PC. High-flow nasal cannula therapy in do-not-intubate patients with hypoxic respiratory distress. *Respir Care*. 2013 Apr;58(4):597-600.
- Möller W, Celik G, Feng S, Bartenstein P, Meyer G, Eickelberg O, et al. Nasal high flow clears anatomical deadspace in upper airway models. *J Appl Physiol*. 2015; 118:1525-32.
- Mündel T, Feng S, Tatkov S, Schneider H. Mechanisms of nasal high flow on ventilation during wakefulness and sleep. *J Appl Physiol*. 2013; 114:1058-65.
- Jeong JH, Kim DH, Kim SC, Kang C, Lee SH, Kang TS, et al. Changes in arterial blood gases after use of high-flow nasal cannula therapy in the ED. *Am J Emerg Med*. 2015; 3(10):1344-9.
- Hasani A, Chapman TH, McCool D, Smith RE, Dilworth JP, Agnew JE. Domiciliary humidification improves lung mucociliary clearance in patients with bronchiectasis. *Chron Respir Dis*. 2008; 5(2):81-6.
- Ritchie JE, Williams AB, Gerard C, Hockey H. Evaluation of a humidified nasal high-flow oxygen system, using oxygraphy, capnography and measurement of upper airway pressures. *Anaesthesia Intensive Care*. 2011; 39(6):1103-10.
- Maggiore SM, Idrone FA, Vaschetto R, Festa R, Cataldo A, Antonicelli F, et al. Nasal High-Flow Versus Venturi Mask Oxygen Therapy after Extubation. Effects on Oxygenation, Comfort, and Clinical Outcome. *Am J Respir Crit Care Med*. 2014; 90(3):282-8.
- Frat JP, Thille AW, Mercat A, Girault C, Ragot S, Perbet S, et al. High-Flow Oxygen through Nasal Cannula in Acute Hypoxic Respiratory Failure. *N Engl J Med*. 2015; 372(23):2185-96.
- Masclans JR, Roca O. High-Flow Oxygen Therapy in Acute Respiratory Failure. *Clin Pulm Med*. 2012; 19(3):127-30.
- Lucangelo U, Vassallo FG, Marras E, Ferluga M, Beziza E, Comuzzi L, et al. High-Flow Nasal Interface Improves Oxygenation in Patients Undergoing Bronchoscopy. *Crit Care Res Pract*. 2012; (12):1-6.
- Parke R, McGuinness S, Eccleston M. A Preliminary Randomized Controlled Trial to Assess Effectiveness of Nasal High-Flow Oxygen in Intensive Care Patients. *Respir Care*. (Mar) 2011; 56(3): 265-70.
- Hernández G, Vaquero C, González P, Subira C, Frutos-Vivar F, Rialp G, et al. Effect of Postextubation High-Flow Nasal Cannula vs Conventional Oxygen Therapy on Reintubation in Low-Risk Patients: A Randomized Clinical Trial. *JAMA*. (Apr) 2016; 315(13):1354-61.
- Hernández G, Vaquero C, Colinas L, Cuenca R, González P, Canabal A, et al. Effect of Postextubation High-Flow Nasal Cannula vs Noninvasive Ventilation on Reintubation and Postextubation Respiratory Failure in High-Risk Patients. *JAMA*. (Oct) 2016; 316(15):1565-74.
- Stéphan F, Barrucand B, Petit P, Rézaiguia-Delclaux S, Médard A, Delannoy B, et al. High-Flow Nasal Oxygen vs Noninvasive Positive Airway Pressure in Hypoxic Patients After Cardiothoracic Surgery: A Randomized Clinical Trial. *JAMA*. 2015; 13(23):2331-9.
- Parke RL, Eccleston ML, McGuinness SP. The Effects of Flow on Airway Pressure During Nasal High-Flow Oxygen Therapy. *Respir Care*. (Aug) 2011; 56(8):1151-5.
- Williams R, Rankin N, Smith T, Galler D, Seakins P. Relationship between the humidity and temperature of inspired gas and the function of the airway mucosa. *Crit Care Med*. 1996; 24(11):1920-9.
- Piraino T, Madden M, Roberts KJ, et al. AARC Clinical Practice Guideline: Management of Adult Patients With Oxygen in the Acute Care Setting. *Respir Care*. 2022 Jan;67(1):115-128.
- Qaseem A, Etzeandia-Ikobaltzeta I, Fitterman N, et al. Appropriate Use of High-Flow Nasal Oxygen in Hospitalized Patients for Initial or Postextubation Management of Acute Respiratory Failure: A Clinical Guideline From the American College of Physicians [published correction appears in *Ann Intern Med*. 2023 Apr;176(4):584].
- Evans L, Rhodes A, Alhazzani W, et al. Surviving Sepsis Campaign: International Guidelines for Management of Sepsis and Septic Shock 2021. *Crit Care Med*. 2021 Nov 1; 49(11):1063-1143.
- Rochwerf B, Elnav S, Chaudhuri D, et al. The role for high flow nasal cannula as a respiratory support strategy in adults: a clinical practice guideline. *Intensive Care Med*. 2020; 46(12):2226-2237.
- Oczkowski S, Ergon B, Bos L, et al. ERS clinical practice guidelines: high-flow nasal cannula in acute respiratory failure. *Eur Respir J*. 2022 Apr 14;59(4):2101574.
- Barnett A, Beasley R, Buchan C, et al. Thoracic Society of Australia and New Zealand Position Statement on Acute Oxygen Use in Adults: 'Swimming between the flags'. *Respirology*. 2022;27(4):262-276.
- World Health Organisation. Clinical management of COVID-19: Living guideline. Geneva: World Health Organisation; 2022 June 23 [cited 2022 Sep]. [www.who.int/publications/i/item/WHO-2019-nCoV-clinical-2022-1](https://www.who.int/publications/i/item/WHO-2019-nCoV-clinical-2022-1).
- Ischaki E, Pantazopoulos I, Zakyntinos S. Nasal high flow therapy: a novel treatment rather than a more expensive oxygen device. *Eur Respir Rev*. 2017; 26(145):170028.
- Cortegiani A, Longhini F, Madotto F, et al. High flow nasal therapy versus noninvasive ventilation as initial ventilatory strategy in COPD exacerbation: a multicenter non-inferiority randomized trial. *Crit Care*. 2020; 24(1):692.
- Pantazopoulos I, Daniil Z, Moylan M, et al. Nasal High Flow Use in COPD Patients with Hypercapnic Respiratory Failure: Treatment Algorithm & Review of the Literature. *COPD*. 2020; 17(1):101-111.
- Roca O, Caralt B, Messika J, et al. An Index Combining Respiratory Rate and Oxygenation to Predict Outcome of Nasal High-Flow Therapy. *Am J Respir Crit Care Med*. 2019;199(11):1368-1376.
- Jackson JA, Spilman SK, Kingery LK, et al. Implementation of High-Flow Nasal Cannula Therapy Outside the Intensive Care Setting. *Respir Care*. 2021; 66(3):357-365.
- Franklin D, Bahl F, Schlapbach L, Oakley E, Craig S, Neutze J, et al. A randomized trial of high-flow oxygen therapy in infants with bronchiolitis. *N Engl J Med*. 2018 Mar;378(12):1121-1131.
- Mileši C, Baleine J, Matecki S, Durand S, Combes C, Novais A, et al. Is treatment with a high-flow nasal cannula effective in acute viral bronchiolitis? A physiologic study. *Intensive Care Med*. 2013 Jun 14; 39(6):1088-1094.
- Rubin S, Ghuman A, Deakers T, Khemani R, Ross P, Newth CJ. Effort of breathing in children receiving high-flow nasal cannula. *Pediatr Crit Care Med*. 2014 Jan;15(1):1-6.
- Pham TM, O'Malley L, Mayfield S, Martin S, Schibler A. The effect of high flow nasal cannula therapy on the work of breathing in infants with bronchiolitis. *Pediatr Pulmonol*. 2015 Jul;50(7):713-20.
- Mayfield S, Bogossian F, O'Malley L, Schibler A. High-flow nasal cannula oxygen therapy for infants with bronchiolitis: pilot study. *J Paediatr Child Health*. 2014 May;50(5):373-8.
- Greenspan JS, Wolfson MR, Shaffer TH. Airway responsiveness to low inspired gas temperature in preterm neonates. *J Pediatr*. 1991 Mar 1;118(3):443-445.
- Kepreotes E, Whitehead B, Attia J, Oldmeadow C, Collison A, Searles A, et al. High-flow warm humidified oxygen versus standard low-flow nasal cannula oxygen for moderate bronchiolitis (HFWHO RCT): an open, phase 4, randomised controlled trial. *Lancet*. 2017 Mar;389(10072):930-939.
- Testa G, Iodice F, Ricci Z, Vitale V, De Razzia F, Haiberger R, et al. Comparative evaluation of high-flow nasal cannula and conventional oxygen therapy in paediatric cardiac surgical patients: a randomized controlled trial. *Interact Cardiovasc Thorac Surg*. 2014 Sep 1;19(3):456-461.
- Durand P, Guddir T, Kyheng C, Blanc F, Vignaud O, Epaul R, et al. A randomised trial of high-flow nasal cannula in infants with moderate bronchiolitis. *Eur Respir J*. 2020 Jul;56(1).
- Wing R, James C, Maranda LS, Armsby CC. Use of high-flow nasal cannula support in the emergency department reduces the need for intubation in pediatric acute respiratory insufficiency. *Pediatr Emerg Care*. 2012 Nov;28(11):1117-1123.
- Schibler A, Pham T, Dunster K, Foster K, Barlow A, Gibbons K, et al. Reduced intubation rates for infants after introduction of high-flow nasal prong oxygen delivery. *Intensive Care Med*. 2011 May 3;37(5):847-852.
- Franklin D, Bahl FE, George S, Oakley E, Borland ML, Neutze J, et al. Effect of Early High-Flow Nasal Oxygen vs Standard Oxygen Therapy on Length of Hospital Stay in Hospitalized Children with Acute Hypoxic Respiratory Failure: The PARIS-2 Randomized Clinical Trial. *JAMA*. 2023 Jan 17; 329(3):224-34.
- Richards-Belle A, Davis P, Drikite L, et al. FIRST-line Support for Assistance in Breathing in Children (FIRST-ABC): a master protocol of two randomised trials to evaluate the non-inferiority of high-flow nasal cannula (HFNC) versus continuous positive airway pressure (CPAP) for non-invasive respiratory support in paediatric critical care. *BMJ Open*. 2020; 10:e038002.
- Ramnarayan et al. Effect of High-Flow Nasal Cannula Therapy vs. Continuous Positive Airway Pressure Therapy on Liberation From Respiratory Support in Acutely Ill Children Admitted to Pediatric Critical Care Units: A Randomized Clinical Trial. *JAMA*. 2022. 328(2), 162-172.

For more information, please contact your local Fisher & Paykel Healthcare representative.

F&P, Optiflow and Airovo are trademarks of Fisher & Paykel Healthcare. For patent information, refer to [www.fphcare.com/ip](http://www.fphcare.com/ip)  
610487 REV E © 2023 Fisher & Paykel Healthcare Limited. Intended for USA only.

[www.fphcare.com](http://www.fphcare.com)

**Fisher & Paykel**  
HEALTHCARE