

F&P **Optiflow** Nasal High Flow



Fisher & Paykel
HEALTHCARE

Understand **Optiflow™** Nasal High Flow




Optiflow Nasal High Flow (NHF) delivers respiratory support to your spontaneously breathing patients, by providing heated, humidified air and oxygen at flow rates up to 60 L/min through the unique Optiflow nasal cannula.

Read on to discover more about:

- mechanisms
- physiological effects
- clinical outcomes and how using Optiflow NHF can reduce escalation, thereby avoiding its associated costs.

MECHANISMS OF ACTION



With Optiflow NHF, you can independently titrate flow and oxygen concentration (FiO₂ 21 - 100%) according to your patient's needs.

The mechanisms of action differ from those of conventional therapies, as do the resulting physiological effects and clinical outcomes.



Read more about mechanisms at:
fphcare.com/opti/mechanisms




Reduction
of dead space

Respiratory support

Dynamic
positive airway
pressure

Airway
hydration

The effects of flow rate on
clearance of rebreathing CO₂




Adapted from Möller et al.¹

Clearance of expired air
in the upper airways¹

Reduces rebreathing of gas with
high CO₂ and depleted O₂¹

Increases alveolar ventilation¹

The effects of NHF on airway pressure,
end-expiratory lung volume and tidal volume




Adapted from Corley et al.²

Breath- and flow-dependent
airway pressure^{3,4}

Promotes slow and
deep breathing³

Increases alveolar ventilation^{1,5}

The effects of high flows of warm,
humidified air on mucociliary transport



100% Humidity 90% Humidity for 15 minutes

400µm

Adapted from Tatkov et al.¹⁴

Optimal
Humidity

Prevents desiccation
of the airway epithelium⁶

Improves mucus clearance^{6,7}

Patient
comfort

Optimal
Humidity


Open system
No seal required

Comfortable^{8,9} and easy to use

Patient tolerance^{8,10}

Supplemental
oxygen
when required

Confidence in
the delivery
of blended,
humidified
oxygen^{2,11}, from
21% to 100%



Adapted from Masclans et al.¹²

PHYSIOLOGICAL EFFECTS & CLINICAL OUTCOMES

The mechanisms of respiratory support, airway hydration, patient comfort and supplemental oxygen contribute to distinct physiological effects...



↑ **IMPROVES** ventilation and gas exchange

↓ **REDUCES** respiratory rate^{5,8,11,13-16}

↓ **REDUCES** carbon dioxide^{1,3,17}

↑ **INCREASES** tidal volume⁵

↑ **INCREASES** end-expiratory lung volume⁵

↑ **IMPROVES** mucus clearance⁷

↑ **IMPROVES** oxygenation^{2,5,8-10,12,13,16,18}

... and clinical outcomes:



Read clinical studies and other evidence at:
fphcare.com/opti/evidence-library



↓ **REDUCES** escalation of care when used:

- as a first-line respiratory support¹⁰
 - post-extubation^{9,19-22}
-

↓ **REDUCES** mortality rate¹⁰

↑ **IMPROVES** symptomatic relief^{8,10,11}

↑ **IMPROVES** comfort and patient compliance^{8,9,11,19,22}

Frat 2015

The New England Journal of Medicine

STUDY



A 23-center study¹⁰ compared NHF to use of a non-rebreather mask (standard oxygen) and NIV as a primary treatment. The primary outcome was the number of patients intubated at day 28 (not attained).

METHOD

310 pre-intubation patients in acute hypoxemic respiratory failure ($\text{PaO}_2/\text{FiO}_2 \leq 300$ mmHg) were randomized to receive NHF, non-rebreather mask or NIV.

RESULTS

- ▶ **NHF significantly reduced ICU ($p=0.047$) and 90-day mortality ($p=0.02$)**
- ▶ The primary outcome was not met for all patients ($p=0.18$), however, **NHF significantly reduced the need for intubation in more acute patients ($\text{PaO}_2/\text{FiO}_2 \leq 200$ mmHg) ($p=0.009$)**
- ▶ Significant increase in ventilator-free days on NHF ($p=0.02$)
- ▶ NHF significantly reduced intensity of respiratory discomfort ($p<0.01$) and dyspnea ($p<0.001$)



Ischaki 2017

European Respiratory Review

Acute hypoxaemic respiratory failure*

Criteria for immediate or imminent intubation are present.

NO

YES

NHF initiation

- FiO_2 100%
- Flow rate 60 L·min⁻¹
- Temperature 37°C

Intubation and invasive MV

- NHF for improving pre-oxygenation and peri-laryngoscopy oxygenation
- FiO_2 100%
 - Flow rate 60 L·min⁻¹

↓ Within 1-2 h

Monitoring

Presence of prognostic factors

NO

YES

Titration

- FiO_2 based on target SpO_2 [$>88-90\%$]
- Flow rate based on $<25-30$ breaths·min⁻¹ and patient comfort
- Temperature based on patient comfort.

Noninvasive MV

Short trial [1-2 h]



Monitoring

Presence of prognostic factors within hours [maximum 48 h]

NO

YES

Weaning from NHF

Firstly decrease FiO_2 . When $\text{FiO}_2 < 0.4$ decrease flow rate by 5 L·min⁻¹.

Intubation and invasive MV

- NHF for improving pre-oxygenation and peri-laryngoscopy oxygenation
- FiO_2 100%
 - Flow rate 60 L·min⁻¹

*Adapted from original paper¹¹; used under Creative Commons licence 4.0. MV = mechanical ventilation; SOT = standard oxygen treatment.

Please note that this material is intended exclusively for healthcare practitioners and the information conveyed constitutes neither medical advice nor instructions for use. This material should not be used for training purposes or to replace individual hospital policies or practices. Before any product use, consult the appropriate user instructions.

Hernández (Apr) 2016

Journal of the American Medical Association

STUDY



A 7-center study²⁰ compared the efficacy of NHF to use of conventional oxygen therapy (COT) post-extubation. The primary outcome was reintubation within 72 hours.

METHOD

527 patients at low risk of reintubation (age < 65; APACHE score < 12; BMI < 30 etc.) were randomized to receive NHF or COT (via nasal prongs or a non-rebreather).

RESULTS

- ▶ **NHF significantly reduced reintubation** ($p=0.004$) and post-extubation respiratory failure ($p=0.03$)
- ▶ Successfully extubated patients (in both groups) had a shorter duration of mechanical ventilation ($p<0.001$), ICU stay ($p<0.001$) and hospital stay ($p=0.005$)



Hernández (Oct) 2016

Journal of the American Medical Association

STUDY

A 3-center non-inferiority study²¹ compared use of NHF to bi-level positive airway pressure (BPAP) post-extubation. The primary outcomes were reintubation and post-extubation respiratory failure within 72 hours.

METHOD

604 patients at high risk of reintubation (age > 65; APACHE score > 12; BMI > 30 etc.) were randomized to receive NHF or BPAP. The non-inferiority margin was 10%.

RESULTS

- ▶ NHF was non-inferior to BPAP for **preventing reintubation**: 22.8% (66/290) NHF group vs. 19.1% (60/314) BPAP group reintubated
- ▶ NHF was non-inferior to BPAP for **preventing post-extubation respiratory failure**: 26.9% (78/290) NHF group vs. 39.8% (125/314) BPAP group had post-extubation respiratory failure
- ▶ No patients in the NHF group suffered adverse effects requiring withdrawal of the therapy, compared to 42.9% of patients in the BPAP group ($p<0.001$)
- ▶ Median ICU length of stay was lower in the NHF group: 3 days (NHF) vs. 4 days (BPAP) ($p=0.048$)



Read clinical studies and other evidence at:
fpicare.com/opti/evidence-library



USAGE


There is an ever-increasing body of clinical literature which may provide guidance on the day-to-day application of Optiflow NHF

When are the effects of Optiflow NHF seen?

Sztrymf¹³ associated Optiflow NHF with sustained beneficial effects on oxygenation and physiological parameters for patients with acute respiratory failure.

Similarly Rittayamai¹⁴ showed significant improvement in post-extubation patients.

These studies may provide guidance on patient responses to the therapy.





View more frequently asked usage questions at:
fphcare.com/opti/usage



Is there a way to predict the outcome of NHF?

The validated ROX index²⁵ predicts failure in adults with AHRF receiving NHF, at 4 time intervals: 2, 6, 12 and > 12 hours. It's an easy-to-use dynamic bedside tool.

$$\frac{\text{SpO}_2 / \text{FiO}_2}{\text{Respiratory Rate}} = \text{ROX index}$$

Example at 6 hours

SpO₂ = 88%

FiO₂ = .70

RR = 28 breaths/minute

$$\frac{88 / .70}{28} = 4.48$$

In the example above, the resulting score of 4.48 is greater than the score for predicted failure at 6 hours (3.47 as shown in the ROX Score table right). Therefore, continued NHF treatment should be considered.

ROX score margin for failure over time

Time Point (Hours of NHF use)	ROX Score	Positive Predictive Value %
2 hours	< 2.85	98
6 hours	< 3.47	98-99
12 hours	< 3.85	99
> 12 hours	< 4.88	80



View more frequently asked usage questions at:
fphcare.com/opti/usage




What flow rates and ranges are used?

The adjacent table lists starting flows and flow ranges used in clinical studies.
5,9,10,13,16,19-22,25-31

Guidance source	Category description	Flow L/min											
		10	15	20	25	30	35	40	45	50	55	60	
RESPIRATORY DISTRESS	Macé et al 2019											●	
	Hernández et al Oct 2016											●	
	Hernández et al Apr 2016					●							
	Bell et al 2015											●	
	Frat et al 2015								■	●			
	Stéphan et al 2015											●	
	Maggiore et al 2014											●	
	Peters et al 2013								■	●			
	Sztrymf et al 2011								■	●			
	Parke et al 2011					●							
	Corley et al 2011					●							
	CHRONIC	Storgaard et al 2018				●							
Nagata et al. 2018							●						
Cirio et al 2016												■	●
Rea et al 2010						■							

Key: ■ Flow range ● Starting flow ● Mean flow




What is the approximate average dynamic pressure generated?

Average pressure increases approximately 0.5 - 1 cmH₂O per 10 L/min.^{2,4,32}

Pressure ranges are cannula and patient dependent. For illustrative purposes only.

COST BENEFITS

Use Optiflow NHF to reduce escalation^{10,20} thereby avoiding associated costs



Using Optiflow NHF as a first-line therapy (both pre-intubation and post-extubation) may reduce a patient's escalation 'up the acuity curve', resulting in better patient outcomes and reduced costs of care.

OPTIFLOW IN PRACTICE



Watch the Berkshire video at:
www.fphcare.com/us/hospital/adult-respiratory/optiflow/optiflow-in-practice/



Introducing **Optiflow** to the **Royal Berkshire** **Hospital**

This video shows the usage of AIRVO 2 & Optiflow Nasal High Flow therapy in different departments of the Royal Berkshire Hospital in Reading, UK. It shows the benefits they have found to both patients and hospital since its introduction.



Evaluate **F&P Optiflow**

Publications in the NEJM and JAMA suggest **Optiflow NHF may improve patient outcomes¹⁰ and reduce the need for higher level support^{20,21} thereby avoiding the associated costs³³.**

Fisher & Paykel Healthcare will provide training and equipment during an Optiflow NHF evaluation to help you achieve these goals in your hospital. Let us customize an evaluation to suit you.
Visit fphcare.com/opti/eval



REFERENCES

1. Möller W, Celik G, Feng S, Bartenstein P, Meyer G, Eickelberg O et al. Nasal high flow clears anatomical deadspace in upper airway models. *J Appl Physiol*. 2015; 118:1525-32.
2. Ritchie JE, Williams AB, Gerard C, Hockey H. Evaluation of a humidified nasal high-flow oxygen system, using oxygraphy, capnography and measurement of upper airway pressures. *Anaesth Intensive Care*. 2011; 39(6):1103-10.
3. Mündel T, Feng S, Tatkov S, Schneider H. Mechanisms of nasal high flow on ventilation during wakefulness and sleep. *J Appl Physiol*. 2013; 114:1058-65.
4. Parke RL, Eccleston ML, McGuinness SP. The Effects of Flow on Airway Pressure During Nasal High-Flow Oxygen Therapy. *Respir Care*. (Aug) 2011; 56(8):1151-5.
5. Corley A, Caruana LR, Barnett AG, Tronstad O, Fraser JF. Oxygen delivery through high-flow nasal cannulae increase end-expiratory lung volume and reduce respiratory rate in post-cardiac surgical patients. *Br J Anaesth*. 2011; 107(6):998-1004.
6. Williams R, Rankin N, Smith T, Galler D, Seakins P. Relationship between the humidity and temperature of inspired gas and the function of the airway mucosa. *Crit Care Med*. 1996; 24(11):1920-9.
7. Hasani A, Chapman TH, McCool D, Smith RE, Dilworth JP, Agnew JE. Domiciliary humidification improves lung mucociliary clearance in patients with bronchiectasis. *Chron Respir Dis*. 2008; 5(2):81-6.
8. Roca O, Riera J, Torres F, Masclans JR. High-Flow Oxygen Therapy in Acute Respiratory Failure. *Respir Care*. 2010; 55(4):408-13.
9. Maggiore SM, Idone FA, Vaschetto R, Festa R, Cataldo A, Antonicelli F et al. Nasal High-Flow Versus Venturi Mask Oxygen Therapy after Extubation. Effects on Oxygenation, Comfort, and Clinical Outcome. *Am J Respir Crit Care Med*. 2014; 90(3):282-8.
10. Frat JP, Thille AW, Mercat A, Girault C, Ragot S, Perbet S et al. High-Flow Oxygen through Nasal Cannula in Acute Hypoxemic Respiratory Failure. *N Engl J Med*. 2015; 372(23):2185-96.
11. Lenglet H, Sztrymf B, Leroy C, Brun P, Dreyfuss D, Ricard JD. Humidified High Flow Nasal Oxygen During Respiratory Failure in the Emergency Department: Feasibility and Efficacy. *Respir Care*. 2012; 57(11):1873-8.
12. Masclans JR, Roca O. High-Flow Oxygen Therapy in Acute Respiratory Failure. *Clin Pulm Med*. 2012; 19(3):127-30.
13. Sztrymf B, Messika J, Bertrand F, Hurel D, Leon R, Dreyfuss D et al. Beneficial effects of humidified high flow nasal oxygen in critical care patients: a prospective pilot study. *Intensive Care Med*. 2011; 37(11):1780-6.
14. Rittayamai N, Tscheikuna J, Rujiwit P. High-Flow Nasal Cannula Versus Conventional Oxygen Therapy After Endotracheal Extubation: A Randomized Crossover Physiologic Study. *Respir Care*. 2014; 59(4): 485-90.
15. Roca O, Pérez-Terán P, Masclans JR, Pérez L, Galve E, Evangelista A et al. Patients with New York Heart Association class III heart failure may benefit with high flow nasal cannula supportive therapy: High flow nasal cannula in heart failure. *J Crit Care*. 2013; 28(5):741-6.
16. Peters S, Holets S, Gay P. High-Flow Nasal Cannula Therapy in Do-Not-Intubated Patients with Hypoxemic Respiratory Distress. *Respir Care*. 2013; 58(4): 597-600.
17. Jeong JH, Kim DH, Kim SC, Kang C, Lee SH, Kang TS et al. Changes in arterial blood gases after use of high-flow nasal cannula therapy in the ED. *Am J Emerg Med*. 2015; 3(10):1344-9.
18. Lucangelo U, Vassallo FG, Marras E, Ferluga M, Beziza E, Comuzzi L et al. High-Flow Nasal Interface Improves Oxygenation in Patients Undergoing Bronchoscopy. *Crit Care Res Pract*. 2012; (12):1-6.

REFERENCES

19. Parke R, McGuinness S, Eccleston M. A Preliminary Randomized Controlled Trial to Assess Effectiveness of Nasal High-Flow Oxygen in Intensive Care Patients. *Respir Care*. (Mar) 2011; 56(3): 265-70.
20. Hernández G, Vaquero C, González P, Subira C, Frutos-Vivar F, Rialp G et al. Effect of Postextubation High-Flow Nasal Cannula vs Conventional Oxygen Therapy on Reintubation in Low-Risk Patients: A Randomized Clinical Trial. *JAMA*. (Apr) 2016; 315(13):1354-61.
21. Hernández G, Vaquero C, Colinas L, Cuenca R, González P, Canabal A et al. Effect of Postextubation High-Flow Nasal Cannula vs Noninvasive Ventilation on Reintubation and Postextubation Respiratory Failure in High-Risk Patients. *JAMA*. (Oct) 2016; 316(15):1565-74.
22. Stéphan F, Barrucand B, Petit P, Rézaiguia-Delclaux S, Médard A, Delannoy B et al. High-Flow Nasal Oxygen vs Noninvasive Positive Airway Pressure in Hypoxemic Patients After Cardiothoracic Surgery: A Randomized Clinical Trial. *JAMA*. 2015; 13(23):2331-9.
23. Ischaki E, Pantazopoulos I, Zakynthinos S. Nasal high flow therapy: a novel treatment rather than a more expensive oxygen device. *Eur Respir Rev*. 2017;26(145):170028.
24. Rittayamai N, Tscheikuna J, Praphruetkit N, Kijpinyochai S. Use of High-Flow Nasal Cannula for Acute Dyspnea and Hypoxemia in the Emergency Department. *Respir Care*. 2015; 60(10):1377-82.
25. Roca O, Caralt B, Messika J, Samper M, Sztrymf B, Hernandez G, et al. An Index Combining Respiratory Rate and Oxygenation to Predict Outcome of Nasal High Flow Therapy. *Am J Respir Crit Care Med*. 2018 Dec 21. [Epub ahead of print]
26. Macé J, Marjanovic N, Faranpour F, Mimoz O, Frerebeau M, Violeau M, et al. Early high-flow nasal cannula oxygen therapy in adults with acute hypoxemic respiratory failure in the ED: A before-after study. *Am J Emerg Med*. 2019 Mar 5. [Epub ahead of print]
27. Bell N, Hutchinson CL, Green TC, Rogan E, Bein KJ, Dinh MM. Randomised control trial of humidified high flow nasal cannulae versus standard oxygen in the emergency department. *Emerg Med Australas*. 2015 Dec; 27(6):537-41.
28. Storgaard LH, Hockey HU, Laursen BS, Weinreich UM. Long-term effects of oxygen-enriched high-flow nasal cannula treatment in COPD patients with chronic hypoxemic respiratory failure. *Int J Chron Obstruct Pulmon Dis*. 2018; 13:1195-205.
29. Nagata K, Kikuchi T, Horie T, Shiraki A, Kitajima T, Kadowaki T et al. Domiciliary High-Flow Nasal Cannula Oxygen Therapy for Patients with Stable Hypercapnic Chronic Obstructive Pulmonary Disease. A Multicenter Randomized Crossover Trial. *Ann Am Thorac Soc*. 2018;15(4):432-9.
30. Cirio S, Piran M, Vitacca M, Piaggi G, Ceriana P, Prazzoli M et al. Effects of heated and humidified high flow gases during high-intensity constant-load exercise on severe COPD patients with ventilatory limitation. *Respir Med*. 2016;118:128-32.
31. Rea H, McAuley S, Jayaram L, Garrett J, Hockey H, Storey L et al. The clinical utility of long-term humidification therapy in chronic airway disease. *Respir Med*. 2010; 104(4): 525-33.
32. Groves N, Tobin A. High flow nasal oxygen generates positive airway pressure in adult volunteers. *Aust Crit Care*. 2007; 20(4):126-31.
33. Eaton Turner E, Jenks M. Cost-effectiveness analysis of the use of high-flow oxygen through nasal cannula in intensive care units in NHS England. *Expert Rev Pharmacoecon Outcomes Res*. 2018; 18(3):331-7.
34. Takov S. *Mucociliary Transport Video Microscopy*. https://www.youtube.com/watch?v=HMdrhwEnY6M&list=PLonAnS_1BEgrIjk745MrBNpRoZFwN114e



Optiflow is a trademark of Fisher & Paykel Healthcare.
For patent information, refer to www.fphcare.com/ip

610487 REV C US © 2019 Fisher & Paykel Healthcare Limited

www.fphcare.com



Fisher & Paykel
HEALTHCARE