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
SSC INTERNATIONALE LEITLINIEN

Evans L, et al. Critical Care Medicine. 2021.²⁵

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- Society of Critical Care Medicine
- American Association of Critical Care Nurses
- American College of Chest Physicians
- American College of Emergency Physicians
- American Thoracic Society

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Gaseem A, et al. Annals of Internal Medicine. 2021.²⁷

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Clinical management of COVID-19: living guideline. 2022.²⁹

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Rochwerg et al. 2019 ³¹

Intensive Care Medicine

High flow nasal cannula compared with conventional oxygen therapy for acute hypoxemic respiratory failure: A SYSTEMATIC REVIEW AND META-ANALYSIS.

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

KI = Konfidenzintervall; RR = relatives Risiko

Frat et al. 2015 ¹⁴

The New England Journal of Medicine

High-flow oxygen through nasal cannula in acute hypoxemic respiratory failure.

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


Ischaki et al. 2017 ³²

European Respiratory Review

Nasal high flow therapy: a novel treatment rather than a more expensive oxygen device.

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Cortegiani et al. 2020 ³³

Critical Care

High flow nasal therapy versus noninvasive ventilation as initial ventilatory strategy in COPD exacerbation: a multicenter non-inferiority randomized trial.

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

Chaudhuri et al. 2020 ³⁵

Chest

High-flow nasal cannula in the immediate postoperative period: a systematic review and meta-analysis.

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
Stephan et al. 2015 ²⁰

Journal of the American Medical Association

High-flow nasal oxygen vs noninvasive positive airway pressure in hypoxemic patients after cardiothoracic surgery: a randomized clinical trial.

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
Yu et al. 2017 ³⁶

Canadian Respiratory Journal

Effect of high-flow nasal cannula versus conventional oxygen therapy for patients with thoracoscopic lobectomy after extubation.

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

Granton et al. 2020 ³⁷

Critical Care Medicine

High-flow nasal cannula compared with conventional oxygen therapy or noninvasive ventilation immediately postextubation: a systematic review and meta-analysis.

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
Hernández et al. (Okt) 2016 ¹⁹

Journal of the American Medical Association

Effect of post-extubation high-flow nasal cannula vs noninvasive ventilation on reintubation and post-extubation respiratory failure in high-risk patients: A randomized clinical trial.

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
Hernández et al. (Apr) 2016 ¹⁸

Journal of the American Medical Association

Effect of post-extubation high-flow nasal cannula vs conventional oxygen therapy on reintubation in low-risk patients.

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
Ergebnisse

Thille et al. 2019 ³⁸*Journal of the American Medical Association*

Effect of post-extubation high-flow nasal oxygen with noninvasive ventilation vs high-flow nasal oxygen alone on reintubation among patients at high risk of extubation failure: a randomized clinical trial.

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

Spoletini et al. 2018 ³⁹*Journal of Critical Care*

High-flow nasal therapy vs standard oxygen during breaks off noninvasive ventilation for acute respiratory failure.

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


Pirret et al. 2017 ⁴⁰*Intensive & Critical Care Nursing*

Nasal high flow oxygen therapy in the ward setting: A prospective observational study.

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- Corley A, Caruana LR, Barnett AG, et al. Oxygen delivery through high-flow nasal cannulae increase end-expiratory lung volume and reduce respiratory rate in post-cardiac surgical patients. *Br J Anaesth*. 2011; 107(6):998-1004.
- Roca O, Riera J, Torres F, Masclans JR. High-Flow Oxygen Therapy in Acute Respiratory Failure. *Respir Care*. 2010; 55(4):408-13.
- Lenglet H, Sztymf B, Leroy C, et al. Humidified High Flow Nasal Oxygen During Respiratory Failure in the Emergency Department: Feasibility and Efficacy. *Respir Care*. 2012; 57(11):1873-8.
- Sztymf B, Messika J, Bertrand F, et al. Beneficial effects of humidified high flow nasal oxygen in critical care patients: a prospective pilot study. *Intensive Care Med*. 2011; 37(11):1780-6.
- Rittayamai N, Tscheikuna J, Rujiwit P. High-Flow Nasal Cannula Versus Conventional Oxygen Therapy After Endotracheal Extubation: A Randomized Crossover Physiologic Study. *Respir Care*. 2014; 59(4): 485-90.
- Roca O, Pérez-Terán P, Masclans JR, et al. Patients with New York Heart Association class III heart failure may benefit with high flow nasal cannula supportive therapy: High flow nasal cannula in heart failure. *J Crit Care*. 2013; 28(5):741-6.
- Peters S, Holets S, Gay P. High-Flow Nasal Cannula Therapy in Do-Not-Intubated Patients with Hypoxemic Respiratory Distress. *Respir Care*. 2013; 58(4): 597-600.
- Möller W, Celik G, Feng S, et al. Nasal high flow clears anatomical deadspace in upper airway models. *J Appl Physiol*. 2015; 118:1525-32.
- Mündel T, Feng S, Tatkov S, Schneider H. Mechanisms of nasal high flow on ventilation during wakefulness and sleep. *J Appl Physiol*. 2013; 114:1058-65.
- Jeong JH, Kim DH, Kim SC, et al. Changes in arterial blood gases after use of high-flow nasal cannula therapy in the ED. *Am J Emerg Med*. 2015; 3(10):1344-9.
- Hasani A, Chapman TH, McCool D, et al. Domiciliary humidification improves lung mucociliary clearance in patients with bronchiectasis. *Chron Respir Dis*. 2008; 5(2):81-6.
- Ritchie JE, Williams AB, Gerard C, Hockey H. Evaluation of a humidified high-flow nasal oxygen system, using oxymography, capnography and measurement of upper airway pressures. *Anaesth Intensive Care*. 2011; 39(6):1103-10.
- Maggiore SM, Idone FA, Vascchetto R, et al. Nasal High-Flow Versus Venturi Mask Oxygen Therapy after Extubation. Effects on Oxygenation, Comfort, and Clinical Outcome. *Am J Respir Crit Care Med*. 2014; 90(3):282-8.
- Frat JP, Thille AW, Mercat A, et al. High-Flow Oxygen through Nasal Cannula in Acute Hypoxemic Respiratory Failure. *N Engl J Med*. 2015; 372(23):2185-96.
- Masclans JR, Roca O. High-Flow Oxygen Therapy in Acute Respiratory Failure. *Clin Pulm Med*. 2012; 19(3):127-30.
- Lucangelo U, Vassallo FG, Marras E, et al. High-Flow Nasal Interface Improves Oxygenation in Patients Undergoing Bronchoscopy. *Crit Care Res Pract*. 2012; 12(1):1-6.
- Parke R, McGuinness S, Eccleston M. A Preliminary Randomized Controlled Trial to Assess Effectiveness of Nasal High-Flow Oxygen in Intensive Care Patients. *Respir Care*. (Mar) 2011; 56(3): 265-70.
- Hernández G, Vaquero C, Colinas L, et al. Effect of Post-extubation High-Flow Nasal Cannula vs Noninvasive Ventilation on Reintubation and Postextubation Respiratory Failure in High-Risk Patients. *JAMA*. (Oct) 2016; 316(15):1565-74.
- Hernández G, Vaquero C, González P, et al. Effect of Postextubation High-Flow Nasal Cannula vs Conventional Oxygen Therapy on Reintubation in Low-Risk Patients: A Randomized Clinical Trial. *JAMA*. (Apr) 2016; 315(13):1354-61.
- Stéphan F, Barrucand B, Petit P, et al. High-Flow Nasal Oxygen vs Noninvasive Positive Airway Pressure in Hypoxemic Patients After Cardiothoracic Surgery: A Randomized Clinical Trial. *JAMA*. 2015; 13(23):2331-9.
- Parke RL, Eccleston ML, McGuinness SP. The Effects of Flow on Airway Pressure During Nasal High-Flow Oxygen Therapy. *Respir Care*. (Aug) 2011; 56(8):1151-5.
- Williams R, Rankin N, Smith T, et al. Relationship between the humidity and temperature of inspired gas and the function of the airway mucosa. *Crit Care Med*. 1996; 24(11):1920-9.
- Rochweg B, Einav S, Chaudhuri D, et al. The role for high flow nasal cannula as a respiratory support strategy in adults: a clinical practice guideline. *Intensive Care Med*. 2020; 46(12):2226-2237.
- Oczkowski S, Ergun B, Bos L, et al. ERS Clinical Practice Guidelines: High-flow nasal cannula in acute respiratory failure. *Eur Respir J*. 2021 Oct 28 [Epub ahead of print].
- Evans L, Rhodes A, Alhazzani W, et al. Surviving Sepsis Campaign: International Guidelines for Management of Sepsis and Septic Shock 2021. *Crit Care Med*. 2021 Nov 1; 49(11):1063-1143.
- Piraino T, Madden M, Roberts KJ, et al. AARC Clinical Practice Guideline: Management of Adult Patients With Oxygen in the Acute Care Setting. *Respir Care*. 2022 Jan;67(1):115-128.
- Qaseem A, Etcheandia-Ikbalzeta I, Fitterman N, et al. Appropriate Use of High-Flow Nasal Oxygen in Hospitalized Patients for Initial or Postextubation Management of Acute Respiratory Failure: A Clinical Guideline From the American College of Physicians [published correction appears in *Ann Intern Med*. 2023 Apr;176(4):584].
- Barnett A, Beasley R, Buchan C, et al. Thoracic Society of Australia and New Zealand Position Statement on Acute Oxygen Use in Adults: 'Swimming between the flags'. *Respirology*. 2022;27(4):262-276.
- World Health Organisation. Clinical management of COVID-19: Living guideline. Geneva: World Health Organisation; 2022 June 23 [cited 2022 Sep]. www.who.int/publications/i/item/WHO-2019-nCoV-clinical-2022-1.
- Tasaka S, Ohshimo S, Takeuchi M, Yasuda H, Ichikado K, Tsushima K, et al. ARDS Clinical Practice Guideline 2021. *J Intensive Care*. 2022; 10(1):32.
- Rochweg B, Granton D, Wang DX, et al. High flow nasal cannula compared with conventional oxygen therapy for acute hypoxemic respiratory failure: a systematic review and meta-analysis. *Intensive Care Med*. 2019; 45(5):563-572.
- Ischaki E, Pantazopoulos I, Zakyntinos S. Nasal high flow therapy: a novel treatment rather than a more expensive oxygen device. *Eur Respir Rev*. 2017; 26(145):170028.
- Cortegiani A, Longhini F, Madotto F, et al. High flow nasal therapy versus noninvasive ventilation as initial ventilatory strategy in COPD exacerbation: a multicenter non-inferiority randomized trial. *Crit Care*. 2020; 24(1):692.
- Pantazopoulos I, Daniil Z, Moylan M, et al. Nasal High Flow Use in COPD Patients with Hypercapnic Respiratory Failure: Treatment Algorithm & Review of the Literature. *COPD*. 2020; 17(1):101-111.
- Chaudhuri D, Granton D, Wang DX, et al. High-Flow Nasal Cannula in the Immediate Postoperative Period: A Systematic Review and Meta-analysis. *Chest*. 2020; 158(5):1934-1946.
- Yu Y, Qian X, Liu C, Zhu C. Effect of High-Flow Nasal Cannula versus Conventional Oxygen Therapy for Patients with Thoracoscopic Lobectomy after Extubation. *Can Respir J*. 2017; 7894631. [Epub].
- Granton D, Chaudhuri D, Wang D, et al. High-Flow Nasal Cannula Compared With Conventional Oxygen Therapy or Noninvasive Ventilation Immediately Postextubation: A Systematic Review and Meta-Analysis. *Crit Care Med*. 2020; 48(11):e1129-e1136.
- Thille AW, Muller G, Gacouin A, et al. Effect of Postextubation High-Flow Nasal Oxygen With Noninvasive Ventilation vs High-Flow Nasal Oxygen Alone on Reintubation Among Patients at High Risk of Extubation Failure: A Randomized Clinical Trial. *JAMA*. 2019; 322(15):1465-1475.
- Spoletini G, Mega C, Pisani L, et al. High-flow nasal therapy vs standard oxygen during breaks off noninvasive ventilation for acute respiratory failure: A pilot randomized controlled trial. *J Crit Care*. 2018; 48:418-425.
- Pirret AM, Takerei SF, Matheson CL, et al. Nasal high flow oxygen therapy in the ward setting: A prospective observational study. *Intensive Crit Care Nurs*. 2017 Oct;42:127-134.
- Roca O, Caralt B, Messika J, et al. An Index Combining Respiratory Rate and Oxygenation to Predict Outcome of Nasal High-Flow Therapy. *Am J Respir Crit Care Med*. 2019;199(11):1368-1376.
- Jackson JA, Spilman SK, Kingery LK, et al. Implementation of High-Flow Nasal Cannula Therapy Outside the Intensive Care Setting. *Respir Care*. 2020; 66(3):357-365.

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610925 REV C © 2023 Fisher & Paykel Healthcare Limited

www.fphcare.com

Fisher & Paykel
HEALTHCARE