F&P **Optiflow** Alto fluxo nasal

Fisher & Paykel
HEALTHCARE

Entenda o Alto fluxo nasal **Optiflow™**

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1. Möller W, Celik G, Feng S, Bartenstein P, Meyer G, Eickelberg O et al. Nasal high flow clears anatomical deadspace in upper airway models. *J Appl Physiol*. 2015; 118:1525-32.
2. Ritchie JE, Williams AB, Gerard C, Hockey H. Evaluation of a humidified nasal high-flow oxygen system, using oxygraphy, capnography and measurement of upper airway pressures. *Anaesthesia Intensive Care*. 2011; 39(6):1103-10.
3. Mündel T, Feng S, Tatkov S, Schneider H. Mechanisms of nasal high flow on ventilation during wakefulness and sleep. *J Appl Physiol*. 2013; 114:1058-65.
4. Parke RL, Eccleston ML, McGuinness SP. The Effects of Flow on Airway Pressure During Nasal High-Flow Oxygen Therapy. *Respir Care*. (Aug) 2011; 56(8):1151-5.
5. Corley A, Caruana LR, Barnett AG, Tronstad O, Fraser JF. Oxygen delivery through high-flow nasal cannulae increase end-expiratory lung volume and reduce respiratory rate in post-cardiac surgical patients. *Br J Anaesth*. 2011; 107(6):998-1004.
6. Williams R, Rankin N, Smith T, Galler D, Seakins P. Relationship between the humidity and temperature of inspired gas and the function of the airway mucosa. *Crit Care Med*. 1996; 24(11):1920-9.
7. Hasani A, Chapman TH, McCool D, Smith RE, Dilworth JP, Agnew JE. Domiciliary humidification improves lung mucociliary clearance in patients with bronchiectasis. *Chron Respir Dis*. 2008; 5(2):81-6.
8. Roca O, Riera J, Torres F, Masclans JR. High-Flow Oxygen Therapy in Acute Respiratory Failure. *Respir Care*. 2010; 55(4):408-13.
9. Maggiore SM, Idone FA, Vaschetto R, Festa R, Cataldo A, Antonicelli F et al. Nasal High-Flow Versus Venturi Mask Oxygen Therapy after Extubation. Effects on Oxygenation, Comfort, and Clinical Outcome. *Am J Respir Crit Care Med*. 2014; 90(3):282-8.
10. Frat JP, Thille AW, Mercat A, Girault C, Ragot S, Perbet S et al. High-Flow Oxygen through Nasal Cannula in Acute Hypoxemic Respiratory Failure. *N Engl J Med*. 2015; 372(23):2185-96.
11. Lenglet H, Sztymf B, Leroy C, Brun P, Dreyfuss D, Ricard JD. Humidified High Flow Nasal Oxygen During Respiratory Failure in the Emergency Department: Feasibility and Efficacy. *Respir Care*. 2012; 57(11):1873-8.
12. Masclans JR, Roca O. High-Flow Oxygen Therapy in Acute Respiratory Failure. *Clin Pulm Med*. 2012; 19(3):127-30.
13. Sztymf B, Messika J, Bertrand F, Hurel D, Leon R, Dreyfuss D et al. Beneficial effects of humidified high flow nasal oxygen in critical care patients: a prospective pilot study. *Intensive Care Med*. 2011; 37(11):1780-6.
14. Rittayamai N, Tscheikuna J, Rujiwit P. High-Flow Nasal Cannula Versus Conventional Oxygen Therapy After Endotracheal Extubation: A Randomized Crossover Physiologic Study. *Respir Care*. 2014; 59(4): 485-90.
15. Roca O, Pérez-Terán P, Masclans JR, Pérez L, Galve E, Evangelista A et al. Patients with New York Heart Association class III heart failure may benefit with high flow nasal cannula supportive therapy: High flow nasal cannula in heart failure. *J Crit Care*. 2013; 28(5):741-6.
16. Peters S, Holets S, Gay P. High-Flow Nasal Cannula Therapy in Do-Not-Intubated Patients with Hypoxemic Respiratory Distress. *Respir Care*. 2013; 58(4): 597-600.
17. Jeong JH, Kim DH, Kim SC, Kang C, Lee SH, Kang TS et al. Changes in arterial blood gases after use of high-flow nasal cannula therapy in the ED. *Am J Emerg Med*. 2015; 3(10):1344-9.
18. Lucangelo U, Vassallo FG, Marras E, Ferluga M, Beziza E, Comuzzi L et al. High-Flow Nasal Interface Improves Oxygenation in Patients Undergoing Bronchoscopy. *Crit Care Res Pract*. 2012; (12):1-6.
19. Parke R, McGuinness S, Eccleston M. A Preliminary Randomized Controlled Trial to Assess Effectiveness of Nasal High-Flow Oxygen in Intensive Care Patients. *Respir Care*. (Mar) 2011; 56(3): 265-70.
20. Hernández G, Vaquero C, González P, Subira C, Frutos-Vivar F, Rialp G et al. Effect of Postextubation High-Flow Nasal Cannula vs Conventional Oxygen Therapy on Reintubation in Low-Risk Patients: A Randomized Clinical Trial. *JAMA*. (Apr) 2016; 315(13):1354-61.
21. Hernández G, Vaquero C, Colinas L, Cuenca R, González P, Canabal A et al. Effect of Postextubation High-Flow Nasal Cannula vs Noninvasive Ventilation on Reintubation and Postextubation Respiratory Failure in High-Risk Patients. *JAMA*. (Oct) 2016; 316(15):1565-74.
22. Stéphan F, Barrucand B, Petit P, Rézaiguia-Delclaux S, Médard A, Delannoy B et al. High-Flow Nasal Oxygen vs Noninvasive Positive Airway Pressure in Hypoxemic Patients After Cardiothoracic Surgery: A Randomized Clinical Trial. *JAMA*. 2015; 13(23):2331-9.
23. Ischaki E, Pantazopoulos I, Zakynthinos S. Nasal high flow therapy: a novel treatment rather than a more expensive oxygen device. *Eur Respir Rev*. 2017;26(145):170028.
24. Rittayamai N, Tscheikuna J, Praphruetkit N, Kijpinyochai S. Use of High-Flow Nasal Cannula for Acute Dyspnea and Hypoxemia in the Emergency Department. *Respir Care*. 2015; 60(10):1377-82.
25. Bell N, Hutchinson CL, Green TC, Rogan E, Bein KJ, Dinh MM. Randomised control trial of humidified high flow nasal cannulae versus standard oxygen in the emergency department. *Emerg Med Australas*. 2015 Dec; 27(6):537-41.
26. Storgaard LH, Hockey HU, Laursen BS, Weinreich UM. Long-term effects of oxygen-enriched high-flow nasal cannula treatment in COPD patients with chronic hypoxemic respiratory failure. *Int J Chron Obstruct Pulmon Dis*. 2018; 13:1195-205.
27. Nagata K, Kikuchi T, Horie T, Shiraki A, Kitajima T, Kadowaki T et al. Domiciliary High-Flow Nasal Cannula Oxygen Therapy for Patients with Stable Hypercapnic Chronic Obstructive Pulmonary Disease. A Multicenter Randomized Crossover Trial. *Ann Am Thorac Soc*. 2018;15(4):432-9.
28. Cirio S, Piran M, Vitacca M, Piaggi G, Ceriana P, Prazzoli M et al. Effects of heated and humidified high flow gases during high-intensity constant-load exercise on severe COPD patients with ventilatory limitation. *Respir Med*. 2016;118:128-32.
29. Rea H, McAuley S, Jayaram L, Garrett J, Hockey H, Storey L et al. The clinical utility of long-term humidification therapy in chronic airway disease. *Respir Med*. 2010; 104(4): 525-33.
30. Groves N, Tobin A. High flow nasal oxygen generates positive airway pressure in adult volunteers. *Aust Crit Care*. 2007; 20(4):126-31.
31. Eaton Turner E, Jenks M. Cost-effectiveness analysis of the use of high-flow oxygen through nasal cannula in intensive care units in NHS England. *Expert Rev Pharmacoecon Outcomes Res*. 2018; 18(3):331-7.