

English (en)

ADULT CIRCUIT INSPIRATORY HEATED WITH MR290 AUTOFEED CHAMBER

For the delivery of breathing gases.

TECHNICAL SPECIFICATIONS

Compatible with Fisher & Paykel Healthcare MR850 humidifiers. Refer to humidifier user instructions.

INTERFACE CONNECTIONS	ISO 5356-1 Conical Connectors
COMPRESSIBLE VOLUME	1.00 L
COMPLIANCE @ 60 cmH₂O	1.64 ± 0.08 mL/cmH ₂ O (including 0.02 mL/cmH ₂ O measurement uncertainty)
RESISTANCE TO FLOW @ 30 L/min	0.74 ± 0.06 cmH ₂ O (including 0.03 cmH ₂ O measurement uncertainty)
MAXIMUM OPERATING PRESSURE	8 kPa
BREATHING TUBE LENGTH	1.6 m
MINIMUM TUBE INTERNAL DIAMETER	20 mm
GAS LEAKAGE @ 60 cmH₂O	<40 mL/min

OVERALL PERFORMANCE AT 20 °C TO 26 °C AMBIENT TEMPERATURE	
	NONINVASIVE MODE
Humidification Output	>12 mg/L
Flow Rate	5 – 120 L/min

SYMBOL DEFINITIONS

	Incorrect water level, replace MR290 chamber	Rx only	Prescription only
	Correct water level in the MR290 chamber	CE 0123	CE Marking 93/42/EEC
	Consult operating instructions		Date of manufacture
	Type BF applied part		Manufacturer
	Not made with phthalates		Use-by date
	Not made with natural rubber latex		Caution/Consult instructions for use
	Single use		7 Days maximum use
	Lot number		Transportation and storage temperature limits
	Reference number		European Union authorised representative

WARNINGS, CAUTIONS, AND NOTES

⚠ WARNINGS

- DO NOT reuse this product. Reuse may result in transmission of infectious substances, interruption to treatment, serious harm or death.
- The use of breathing circuits, chambers, accessories or combinations which are not recommended by Fisher & Paykel Healthcare may result in poor humidification system performance, ventilator malfunction, and harm to the patient/user.
- Appropriate patient monitoring (e.g. oxygen saturation) must be used at all times. Failure to monitor the patient (e.g. in the event of an interruption to gas flow) may result in serious harm or death.
- DO NOT touch the heater-plate or chamber base. Surfaces may exceed 85 °C. Failure to comply may result in a skin burn.

- Failure to comply with the following warnings may impair performance of the device or compromise safety (including potentially cause serious harm):**
- When mounting a humidifier adjacent to a patient, ensure that the humidifier is always positioned lower than the patient.
 -  DO NOT use beyond 7 days maximum duration of use.
 - DO NOT cover the circuit with materials such as blankets, towels or bed linen.
 - DO NOT stretch or milk the tubing.
 - DO NOT soak, wash, or sterilize this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.
 - DO NOT use the chamber if the water level rises above the maximum water level line.
 - DO NOT use the chamber if the seals are not intact when received, or if it has been dropped.
 - DO NOT operate the chamber at an angle in excess of 10°.
 - DO NOT spike the water source until the blue caps have been removed. Should the primary float fail, splashing into the circuit may occur if the chamber is being operated in excess of 80 L/min.
 - Visually inspect breathing sets for damage (e.g. a crushed tube or cracked connector) before use, and replace if damaged.
 - DO NOT crush tube.
 - The water source must be at least 50 cm higher than the chamber.
 - DO NOT use CPAP valve with air entrainer.
 - DO NOT use heated wire breathing circuits without gas flow. If gas flow is interrupted, turn the humidifier off.

- Ensure there is a water supply connected to the chamber and that water is present within the chamber.
- DO NOT fill the chamber with water in excess of 37 °C.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

Attention

- Ensure appropriate ventilator or flow source alarms are set before connecting breathing set to patient.
- Avoid prolonged contact of heated tube with patient's skin.
- Perform a pressure and leak test on the breathing system and check for occlusions before connecting to a patient.
- Use USP sterile water for inhalation or equivalent for humidification. DO NOT add other substances to the water.
- Regularly monitor and drain condensate build-up in the circuit.
- Check all connections are tight before use.
- Check that the heater wire is evenly distributed along the circuit and not bunched or kinked.

NOTES

- For use under the supervision of trained medical personnel.
- Dispose of product according to hospital protocol. User may be exposed to breathing tract fluids during disposal.
- The responsible organization is accountable for the compatibility of the humidifier and all of the parts and accessories used to connect to the patient and other equipment before use.

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