
F&P 850 System

F&P Optiflow

Adult Optiflow™ Circuit Kit
with MicroCell Technology™



USER INSTRUCTIONS

REF RT232

⊗ -10 °C 50 °C 14 CE 0123 MD Rx only

www.fphcare.com Fisher & Paykel HEALTHCARE

Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Email: info@fphcare.co.nz Web: www.fphcare.com

Australia (AU) (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 **Brazil (BR)** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo - SP, Brazil Tel: +55 11 2548 7002 **China (CN)** 代理人/售后服务机构:费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 **France (FR)** [EC REP] Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Email: c.s@fphcare.fr **Germany (DE)** Fisher & Paykel Healthcare GmbH, Wiesenstrasse 49, 73614 Schorndorf, Germany Tel: +49 7181 98599 0 **Hong Kong (HK)** Tel: +852 2116 0032 **India (IN)** Tel: +91 80 2309 6400 **Japan (JP)** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea (KR)** Tel: +82 2 6205 6900 **Mexico (MX)** Tel: +52 55 9130 1626 **Poland (PL)** Fisher & Paykel Healthcare Poland Sp. z o.o., Pl. Andersa 7, 61-894 Poznań, Poland Tel: +48 664 846 464 **Russia (RU)** Tel: +7 495 782 21 50 **Switzerland (CH)** [CH REP] Fisher & Paykel Healthcare GmbH, Säntisstrasse 2, 9501 Wil / SG, Switzerland Tel: 0800 83 47 63 **Taiwan (TW)** Tel: +886 2 8751 1739 **Turkey (TR)** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Östim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 **UK (GB)** [UK REP] Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 **USA (US)/Canada (CA)** Tel: 1800 446 3908 or +1 949 453 4000

Austria (AT) Tel: 0800 29 31 23 **Benelux (BE NL LU)** Tel: +31 40 216 3555 **Denmark (DK)** Tel: +45 70 26 37 70 **Finland (FI)** Tel: +358 9 251 66 123 **Ireland (IE)** Tel: 1800 409 011 **Italy (IT)** Tel: +39 06 7839 2939 **Norway (NO)** Tel: +47 21 60 13 53 **Spain (ES)** Tel: +34 902 013 346 **Sweden (SE)** Tel: +46 8 564 76 680

For patent information, see www.fphcare.com/ip
F&P 850 System, F&P Optiflow and MicroCell Technology are trademarks of Fisher & Paykel Healthcare Ltd.

REF 900165 REV C 2023-08 © 2023 Fisher & Paykel Healthcare Limited

English en

Symbols

	This product is intended to be used for a maximum of 14 days		Storage temperature range
	Single Use Only		Use by
	Incorrect water level, replace MR290 Chamber		Correct water level in the MR290 chamber
	Authorised representative in the European community		Consult instructions for use.
	European Conformity		Date and country of manufacture CC NZ: New Zealand MX: Mexico
	Manufacturer		Catalogue number
	Batch code	Rx only	Prescription only
	Medical device		European Union authorised representative
	Switzerland authorised representative		Medical Device
	UK responsible person		Caution/Consult instructions for use
	Distributor		Importer

Intended Use

The Fisher & Paykel Healthcare RT232 is a single use consumable kit intended for use with the MR850 to deliver heated and humidified nasal high flow therapy to spontaneously breathing adult patients who require breathing support. The product is designed for use in the hospital environment and must be prescribed by a physician.

Technical Specifications

Unless otherwise stated, flow rates are expressed at BTPS (body temperature and pressure, saturated) conditions.

Inspiratory Limb Length	1.5m
Compressible volume	1 L
Flow Range	<60 L/min (breathing circuit)
Connection to Flow Source	ISO 5356-1 Conical Connector
Minimum Internal Diameter	14mm
Flow Resistance @ 30 L/min BTPS (including 0.02 cmH ₂ O measurement uncertainty)	0.85 ± 0.07 cmH ₂ O
Maximum tube surface temperature	44°C
Maximum chamber operating pressure	8 kPa
Gas Leakage @ 25 cmH ₂ O	<0.6 L/min

OVERALL PERFORMANCE AT 20°C TO 26°C AMBIENT TEMPERATURE	
Humidifier Mode	INVASIVE MODE
Humidification Output	>33 mg/L
Flow Rate	5–60 L/min

Setup

- MR850 Humidifier in invasive mode. RT232 kit with 22mm heated inspiratory tube and chamber, ventilator in high flow mode.
- Refer to interface user instructions for maximum flow rates (e.g. OPT944, OPT1044)

Attention

- Check all connections are tight before use.
- Check that the heater wire is evenly distributed along the circuit and not bunched or kinked.
- Avoid prolonged contact with patient’s skin.
- Set appropriate ventilator alarms.
- Check breathing circuit and interface for condensation every 6 hours. Drain as required to prevent loss of therapy.
- Before connecting the interface, check for adequate gas flow and ensure that the system has warmed up.
- Use USP sterile water for inhalation or equivalent

Warnings

Failure to comply with the following warnings may impair performance of the device or compromise safety (including potentially cause serious harm):

- There are residual risks associated with use of this product, even if used as intended. Following all instructions and warnings provided, the risks of hypoxic injury, skin burns, airway burns, airway injury, lung injury, electric shock injury, musculoskeletal injury, and hypothermia remain. These risks may result in serious injury or death.
- Visually inspect components and accessories for damage before use and replace if damaged. Use of damaged components or accessories may impair performance of this device or compromise safety.
- The use of breathing circuit/chamber/interface combinations not recommended by Fisher & Paykel Healthcare may result in poor humidification performance, ventilator malfunction and harm to the patient/user.
- For single patient use only. Reuse may result in transmission of infectious substances, interruption to treatment, or serious harm. Attempting to reprocess will result in degradation of materials and render the product defective.
- Nasal delivery of respiratory gases generates flow-dependent positive airway pressure (PAP). This must be taken into account where PAP could have adverse effects on a patient.
- Appropriate patient monitoring must be used at all times. Failure to monitor the patient may result in loss of therapy, serious injury or death.
- When mounting a humidifier adjacent to a patient ensure that the humidifier is always positioned lower than the patient.
- Ensure there is a water supply connected to the chamber and that water is present within the chamber.
- Do not use if packaging is not sealed.
- Do not use heated wire breathing circuits without gas flow. If gas flow is interrupted turn the humidifier off.
- Do not use the chamber if the seals are not intact when received, or if it has been dropped.
- Do not spike the water source until the blue caps have been removed.
- Do not reinsert the blue cap into the chamber after removal. This could damage components, impair the performance of this device or compromise safety.
- Do not use the chamber if the water level rises above the maximum water level line.
- Do not operate the chamber at an angle in excess of 10°.
- Do not fill the chamber with water in excess of 37 °C.
- Do not touch the heater plate or chamber base. Surfaces may exceed 85 °C.
- Do not use near a naked flame, to avoid fires.
- Do not stretch or milk the tubing.
- Do not cover circuit with materials such as blankets, towels or bed linen as this can affect the quality of the therapy or injure the patient.
- Do not soak, wash, sterilize, or reuse this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65.

Notes

- Discard according to appropriate hospital protocol. User may be exposed to breathing tract fluids during disposal.
- This breathing set is validated for use with MR850-system humidifiers only.
- The responsible organisation is accountable for the compatibility of the humidifier and all of the parts and accessories used to connect to the patient and other equipment before use.
- Notice to User: If a serious incident has occurred while using this device, please inform your local Fisher & Paykel Healthcare representative and Competent Authority in your country.

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