



USER INSTRUCTIONS REF RT319

CE 0123 Rx only

MD

-10 °C 50 °C

For patent information, see www.fphcare.com/ip

Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Email: info@fphcare.co.nz Web: www.fphcare.com
Australia (AU) (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 **Brazil (BR)** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo – SP, Brazil Tel: +55 11 2548 7002 **China (CN)**
 代理人/售后服务机构: 费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 **France (FR)** EC REP Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France
 Tel: +33 1 6446 5201 Email: c.s@fphcare.fr **Germany (DE)** Fisher & Paykel Healthcare GmbH, Wiesenstrasse 49, 73614 Schorndorf, Germany Tel: +49 7181 98599 0 **Hong Kong (HK)** Tel: +852 2116 0032 **India (IN)** Tel: +91 80 2309 6400 **Japan (JP)** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea (KR)**
 Tel: +82 2 6205 6900 **Mexico (MX)** Tel: +52 55 9130 1626 **Poland (PL)** Fisher & Paykel Healthcare Poland Sp. z o.o., Pl. Andersa 7, 61-894 Poznań, Poland Tel: +48 664 846 464 **Russia (RU)** Tel: +7 495 782 21 50 **Switzerland (CH)** CH REP Fisher & Paykel Healthcare GmbH, Säntisstrasse 2,
 9501 Wil / SG, Switzerland Tel: 0800 83 47 63 **Taiwan (TW)** Tel: +886 2 8751 1739 **Turkey (TR)** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 **UK (GB)** UK REP
 Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 **USA (US)/Canada (CA)** Tel: 1800 446 3908 or +1 949 453 4000

Austria (AT) Tel: 0800 29 31 23 **Benelux (BE NL LU)** Tel: +31 40 216 3555 **Denmark (DK)** Tel: +45 70 26 37 70 **Finland (FI)** Tel: +358 9 251 66 123 **Ireland (IE)** Tel: 1800 409 011 **Italy (IT)** Tel: +39 06 7839 2939 **Norway (NO)** Tel: +47 21 60 13 53 **Spain (ES)** Tel: +34 902 013 346 **Sweden (SE)**
 Tel: +46 8 564 76 680

English en[ⓘ]

BI-LEVEL /CPAP CIRCUIT

The RT319 circuit is indicated for use with spontaneously breathing adult patients requiring respiratory support, Bi-level or CPAP therapy. Therapy is to be set up and maintained by trained medical professionals in a hospital or institutional environment with patient monitoring in place.

CONTRAINDICATIONS

When using the exhalation port, it should not be used on patients who:

- Are unconscious, unable to breathe spontaneously, uncooperative or unresponsive.
- Have copious secretions, at risk of nausea/vomiting, or at high risk of aspiration of emesis.

If symptoms of these conditions occur, discontinue treatment immediately.

TECHNICAL SPECIFICATIONS

Compatible with Fisher & Paykel Healthcare MR850 Humidifiers. Refer to the MR850 Respiratory Humidifier instructions for use for additional safety information, clinical claims and instructions for setup and use.

Unless otherwise stated, flow rates are expressed at BTPS (body temperature and pressure, saturated) conditions.

Inspiratory Limb with Extension, Filter, and MR290 Chamber


INTERFACE CONNECTIONS	ISO 5356-1 Conical Connectors
COMPRESSIBLE VOLUME	1,7 L
COMPLIANCE @ 60 cmH₂O	4.28 ± 0.44 mL/cmH ₂ O
	(including 0.08 mL/cmH ₂ O measurement uncertainty)
RESISTANCE TO FLOW @ 30 L/min	
With Accessories	0.64 ± 0.07 cmH ₂ O
Without Accessories	0.26 ± 0.04 cmH ₂ O
	(including 0.03 cmH ₂ O measurement uncertainty)
MINIMUM OPERATING PRESSURE (BI-LEVEL/CPAP)	4 cmH ₂ O
MAXIMUM OPERATING PRESSURE	80 cmH ₂ O
CIRCUIT LENGTH	1.8 m
GAS LEAKAGE @ 60 cmH₂O	<40 mL/min

OVERALL PERFORMANCE AT 20 °C TO 26 °C AMBIENT TEMPERATURE		
Humidifier Mode	NONINVASIVE MODE	INVASIVE MODE
Humidification Output	>12 mg/L	>33 mg/L
Flow Rate	3 – 120 L/min	16 – 60 L/min

FILTRATION EFFICIENCY

Viral Organism	>99.94%
Bacterial Organism	ΦX174 Bacteriophage
Mean particle size	>99.996%
NaCl	Staphylococcus aureus
	3 µm
	96.73%

SYMBOL DEFINITIONS

	Incorrect water level, replace MR290 chamber		Correct water level in the MR290 chamber
Rx only	Prescription only		European Conformity
	Consult operating instructions		Date and country of manufacture CC NZ: New Zealand MX: Mexico
	Use-by date		Manufacturer
	Lot number		Single use
	Reference number		Caution/Consult instructions for use
	Medical device		7 Days maximum use
	European Union authorised representative		Transportation and storage temperature limits
	Distributor		Importer
	UK responsible person		Switzerland authorised representative

WARNINGS, CAUTIONS, AND NOTES

WARNINGS

- There are residual risks associated with use of this product, even if used as intended. Following all instructions and warnings provided, the risks of hypoxic injury, skin burns, airway burns, airway injury, lung injury, electric shock injury, musculoskeletal injury, and hypothermia remain. These risks may result in serious injury or death.
- DO NOT reuse this product. Reuse may result in transmission of infectious substances, interruption to treatment, serious harm or death.

- The use of breathing circuits, chambers, accessories, or combinations which are not approved by Fisher & Paykel Healthcare may result in poor humidification system performance, ventilator malfunction and harm to the patient/user.
- Appropriate patient monitoring (e.g. oxygen saturation) must be used at all times. Failure to monitor the patient (e.g. in the event of an interruption to gas flow) may result in serious harm or death.
- DO NOT touch the heater-plate or chamber base. Surfaces may exceed 85 °C. Failure to comply may result in a skin burn.
- Exhalation port must be used with a non-vented interface on a single limb system. Failure to comply may lead to patient inhaling excess carbon dioxide resulting in hypercapnia.
- Monitor the filter when in use. Replace the filter following standard hospital procedures if blockage occurs due to:
 - Excessive condensation
 - Nebulized medication.

- Failure to do so may result in excessive carbon dioxide inhalation.
- DO NOT use this circuit to provide pressure support therapy via a tracheostomy tube. Copious secretions produced by tracheostomy patients may block the disposable exhalation port and lead to excessive carbon dioxide rebreathing.

Failure to comply with the following warnings may impair performance of the device or compromise safety (including potentially cause serious harm):

- Ensure that the humidifier is always positioned lower than the patient’s airway.
-  DO NOT use beyond 7 days maximum duration of use.
- DO NOT cover circuit with materials such as blankets, towels or bed linen.
- DO NOT cover circuit with materials such as blankets, towels or bed linen.
- DO NOT soak, wash, or sterilize this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.
- DO NOT use the chamber if the water level rises above the maximum water level line.
- DO NOT use the chamber if the seals are not intact when received, or if it has been dropped.
- DO NOT operate the chamber at an angle in excess of 10°.
- DO NOT spike the water source until the blue caps have been removed. Should the primary float fail, splashing into the circuit may occur if the chamber is being operated in excess of 80 L/min.
- Visually inspect breathing sets for damage (e.g. a crushed tube or cracked connector) before use, and replace if damaged.
- DO NOT crush tube.
- The water source must be at least 50 cm higher than the chamber.
- DO NOT use heated wire breathing circuits without gas flow. If gas flow is interrupted, turn the humidifier off.
- Ensure there is a water supply connected to the chamber and that water is present within the chamber.
- Change filter every 24 hours, or sooner if noticeable deterioration occurs, following standard hospital procedure.
- DO NOT fill the chamber with water in excess of 37 °C.
- DO NOT block or seal the vent holes on the exhalation port.
- DO NOT use with any expiratory limbs.
- Monitor the filter when in use. Replace the filter following standard hospital procedures if blockage occurs due to:
 - Excessive condensation
 - Nebulized medication

- Failure to do so may result in excessive carbon dioxide inhalation.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit www.fhpcare.com/prop65

Attention

- Ensure appropriate ventilator or flow source alarms are set before connecting breathing set to patient.
- Avoid prolonged contact of heated tubes with patient’s skin.
- Before connecting to patient, ensure that flow and pressure testing applicable to the ventilator has been completed.
- Perform a pressure and leak test on the breathing system and check for occlusions before connecting to a patient.
- Use USP sterile water for inhalation or equivalent for humidification. DO NOT add other substances to the water.
- Regularly monitor and drain condensate build-up in the circuit.
- Check all connections are tight before use.
- Check that the heater wire is evenly distributed along the circuit and not bunched or kinked.
- When nebulized drugs are used, resistance to flow should be monitored, and the filter replaced following standard hospital procedure.
- Check for occlusions in both inspiratory and pressure line before connecting to patient.

NOTES

- For use under the supervision of trained medical personnel.
- Refer to Hospital protocol for disposal methods. User may be exposed to breathing tract fluids during disposal.
- The responsible organization is accountable for the compatibility of the humidifier and all of the parts and accessories used to connect to the patient and other equipment before use.
- Notice to User: If a serious incident has occurred while using this device, please inform your local Fisher & Paykel Healthcare representative and Competent Authority in your country.

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