

Infant continuous flow circuit (0.3 to 10 L/min) single heated with MR290 autofeed chamber and pressure relief valve

**USER INSTRUCTIONS REF RT329**

**CE 0123 Rx only**

For patent information, see [www.fphcare.com/ip](http://www.fphcare.com/ip)

**Manufacturer** Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: [info@fphcare.co.nz](mailto:info@fphcare.co.nz) Web: [www.fphcare.com](http://www.fphcare.com) **Importer/ Distributor Australia** (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 Fax: +61 3 9871 4998 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **Brazil** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraisópolis, 04004-000, São Paulo - SP, Brazil Tel: +55 11 2548 7002 **China** 代理人/售后服务机构 费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 传真: +86 20 32052132 **Denmark** Tel: +45 70 26 37 70 Fax: +46 83 66 310 **Finland** Tel: +358 9 251 66 123 Fax: +46 83 66 310 **France** Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 Email: [c.s@fphcare.fr](mailto:c.s@fphcare.fr) **Germany** Fisher & Paykel Healthcare GmbH, Deutschland, Österreich, Schweiz, Wiesenstrasse 49, D 73614 Schorndorf, Germany Tel: +49 7181 98599 0 Fax: +49 7181 98599 66 **Hong Kong** Tel: +852 2116 0032 Fax: +852 2116 0085 **India** Tel: +91 80 2309 6400 **Ireland** Tel: 1800 409 011 Fax: +44 1628 626 146 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea** Tel: +82 2 6205 6900 Fax: +82 2 6309 6901 **Mexico** Tel: + 52 55 9130 1626 **Norway** Tel: +47 21 60 13 53 Fax: +47 22 99 60 10 **Poland** Tel: 00800 49 119 77 Fax: 00800 49 119 78 **Russia** and Fax: +7 495 782 21 50 **Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **Taiwan** Tel: +886 2 8751 1739 Fax: +886 2 8751 5625 **Turkey** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 Fax: +90 312 354 31 01 **UK** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 Fax: +44 1628 626 146 **USA/Canada** Tel: 1800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001

English <sup>en</sup>

**INFANT CONTINUOUS FLOW CIRCUIT (0,3 TO 10 L/MIN) SINGLE HEATED WITH MR290 AUTOFEED CHAMBER AND PRESSURE RELIEF VALVE**

For delivery of breathing gases.

**TECHNICAL SPECIFICATIONS**

Compatible with Fisher & Paykel Healthcare MR850 humidifiers. Refer to humidifier user instructions.

<b>INTERFACE CONNECTIONS</b>	ISO 5356-1 Conical Connectors
<b>RESISTANCE TO FLOW @ 2.5 L/min</b> (including 0.03 cmH <sub>2</sub> O measurement uncertainty)	0.10 ± 0.05 cmH <sub>2</sub> O
<b>MAXIMUM OPERATING PRESSURE</b>	8 kPa
<b>BREATHING TUBE LENGTH</b>	1.6 m
<b>GAS LEAKAGE @ 60 cmH<sub>2</sub>O</b>	<30 mL/min
<b>MANIFOLD RELIEF PRESSURE</b>	40 cmH <sub>2</sub> O

**OVERALL PERFORMANCE AT 20 °C TO 26 °C AMBIENT TEMPERATURE**








	INVASIVE MODE
Humidification Output	>12 mg/L
Flow Rate	0.3-10 L/min

**APPROVED FOR USE WITH CANNULA**

PART NUMBER	CANNULA SIZE	MAXIMUM DELIVERED FLOW (L/MIN)*
BC2425	Premature	6.0
BC2435	Neonatal	6.0
BC2745	Infant	7.0
BC2755	Intermediate Infant	7.0
BC3780	Pediatric	8.0

\*Flows are to be used as a guide only.

**SYMBOL DEFINITIONS**

	Incorrect water level, replace MR290 chamber	<b>Rx only</b>	Prescription only
	Correct water level in the MR290 chamber		CE Marking 93/42/EEC
	Consult operating instructions		Date of manufacture
	Type BF applied part		Manufacturer
	Not made with phthalates		Use-by date
	Not made with natural rubber latex		Caution/Consult instructions for use
	Single use		7 Days maximum use
	Lot number		Transportation and storage temperature limits
	Reference number		European Union authorised representative

**WARNINGS, CAUTIONS, AND NOTES**

** WARNINGS**

- DO NOT reuse this product. Reuse may result in transmission of infectious substances, interruption to treatment, serious harm or death.
- The use of breathing circuits, chambers, accessories, or combinations which are not approved by Fisher & Paykel Healthcare may result in poor humidification system performance, ventilator malfunction, and harm to the patient/user.
- Appropriate patient monitoring (e.g. oxygen saturation) must be used at all times. Failure to monitor the patient (e.g. in the event of an interruption to gas flow) may result in serious harm or death.
- DO NOT touch the heater plate or chamber base. Surfaces may exceed 85 °C. Failure to comply may result in a skin burn.
- DO NOT bypass the pressure relief valve as this may cause damage/injury.

**Failure to comply with any of the following warnings may impair performance of the device or compromise safety (including potentially cause serious harm):**

- When mounting a humidifier adjacent to a patient, ensure that the humidifier is always positioned lower than the patient.
-  DO NOT use beyond 7 days maximum duration of use.
- DO NOT cover the circuit with materials such as blankets, towels or bed linen.
- DO NOT stretch or milk the tubing.
- DO NOT soak, wash, or sterilize this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.

- DO NOT use the chamber if the water level rises above the maximum water level line.
- DO NOT use the chamber if the seals are not intact when received, or if it has been dropped.
- DO NOT operate the chamber at an angle in excess of 10°.
- DO NOT spike the water source until the blue caps have been removed. Should the primary float fail, splashing into the circuit may occur if the chamber is being operated in excess of 80 L/min.
- Visually inspect breathing sets for damage (e.g. a crushed tube or cracked connector) before use, and replace if damaged.
- The water source must be at least 50 cm higher than the chamber.
- DO NOT use heated wire breathing circuits without gas flow. If gas flow is interrupted, turn the humidifier off.
- Heated tubes should be kept out of direct contact with the patient’s skin. At flow rates below 0.8 L/min, the surface temperature can reach 50 °C.
- Ensure there is a water supply connected to the chamber and that water is present within the chamber.
- Fill the chamber before connecting to the flow source.
- DO NOT fill the chamber with water in excess of 37 °C.
- DO NOT connect the pressure relief valve to the outlet port of the chamber, as this may impair the operation of the valve.
- Only use pressure relief valve with blue cover with Nasal High Flow Therapy.
- Discard product if the blue cover on the pressure relief valve has been damaged or is missing.
- Discard product if the pressure relief valve is damaged or damage is suspected.
- Ensure any unused ports have their caps and/or plugs in place before use.

- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65
- Ensure appropriate ventilator or flow source alarms are set before connecting breathing set to patient.
- Avoid prolonged contact of heated tubes with patient’s skin.
- Perform a pressure and leak test on the breathing system and check for occlusions before connecting to a patient.
- Use USP sterile water or equivalent for humidification. DO NOT add other substances to the water.

- Regularly monitor and drain condensate build-up in the circuit.
- Check that all connections, caps and/or plugs are tight before use.
- Check that the heater wire is evenly distributed along the circuit and not bunched or kinked.
- DO NOT use with MR810 humidifier.
- Delivered flow to the patient may differ from the set flow depending on cannula type and size.
- Check that the cannula fits the patient correctly. Ensure that the patient’s nares and mouth remain clear at all times.
- DO NOT use an in-line nebulizer between the pressure relief manifold and the dry side of the humidification chamber.

**NOTES**

- For use under the supervision of trained medical personnel.
- Dispose of product according to hospital protocol. User may be exposed to breathing tract fluids during disposal.
- Ensure clinical judgement is used when prescribing flow rates.
- The responsible organization is accountable for the compatibility of the humidifier and all of the parts and accessories used to connect to the patient and other equipment before use.

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