


USER INSTRUCTIONS

REF

RT330

Optiflow™ Tubing Kit

CE 0123

Rx only

MD

-10 °C

50 °C

F&P and Optiflow are trademarks of Fisher & Paykel Healthcare Limited.
For patent information, see www.fphcare.com/ip

Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Email: info@fphcare.co.nz Web: www.fphcare.com
Australia (AU) (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 **Brazil (BR)** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo – SP, Brazil Tel: +55 11 2548 7002 **China (CN)** 代理人/售后服务机构:费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 **France (FR)** **[EC REP]** Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Email: c.s.@fphcare.fr **Germany (DE)** Fisher & Paykel Healthcare GmbH, Wiesenstrasse 49, 73614 Schorndorf, Germany Tel: +49 7181 98599 0 **Hong Kong (HK)** Tel: +852 2116 0032 **India (IN)** Tel: +91 80 2309 6400 **Japan (JP)** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea (KR)** Tel: +82 2 6205 6900 **Mexico (MX)** Tel: +52 55 9130 1626 **Poland (PL)** Fisher & Paykel Healthcare Poland Sp. z o.o., Pl. Andersa 7, 61-894 Poznań, Poland Tel: +48 664 846 464 **Russia (RU)** Tel: +7 495 782 21 50 **Switzerland (CH)** **[CH REP]** Fisher & Paykel Healthcare GmbH, Säntisstrasse 2, 9501 Wil / SG, Switzerland Tel: 0800 83 47 63 **Taiwan (TW)** Tel: +886 2 8751 1739 **Turkey (TR)** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 **UK (GB)** **[UK REP]** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 **USA (US)/Canada (CA)** Tel: 1800 446 3908 or +1 949 453 4000

Austria (AT) Tel: 0800 29 31 23 **Benelux (BE NL LU)** Tel: +31 40 216 3555 **Denmark (DK)** Tel: +45 70 26 37 70 **Finland (FI)** Tel: +358 9 251 66 123 **Ireland (IE)** Tel: 1800 409 011 **Italy (IT)** Tel: +39 06 7839 2939 **Norway (NO)** Tel: +47 21 60 13 53 **Spain (ES)** Tel: +34 902 013 346 **Sweden (SE)** Tel: +46 8 564 76 680

English (en)

OPTIFLOW TUBING KIT

For use with High Flow Nasal Therapy using a Fisher & Paykel Healthcare approved nasal cannula.

TECHNICAL SPECIFICATIONS

Compatible with Fisher & Paykel Healthcare MR850 humidifiers. Refer to the MR850 Respiratory Humidifier instructions for use for additional safety information, clinical claims and instructions for setup and use.

Approved for use with Fisher & Paykel Optiflow Junior and Optiflow Junior 2/2+ nasal cannula. Refer to cannula user instructions for maximum flow rates*.

*Flows are to be used as a guide only.

Unless otherwise stated, flow rates are expressed at BTPS (body temperature and pressure, saturated) conditions.

INTERFACE CONNECTIONS	ISO 5356-1 Conical Connectors
RESISTANCE TO FLOW (WITHOUT CANNULA) @ 2.5 L/min:	0.10 ± 0.03 cmH ₂ O
	(including 0.03 cmH ₂ O measurement uncertainty)
MAXIMUM OPERATING PRESSURE	8 kPa
BREATHING TUBE LENGTH	1.75 m / 5' 9"
GAS LEAKAGE @ 60 cmH₂O	<30 mL/min
MANIFOLD RELIEF PRESSURE	40 cmH ₂ O (0.57 psi/3.9 kPa/ 0.04 bar)

OVERALL PERFORMANCE AT 20 °C TO 26 °C AMBIENT TEMPERATURE	
Humidifier Mode	INVASIVE MODE
Humidification Output	>12 mg/L
Flow Rate	1.1 - 36 L/min

SYMBOL DEFINITIONS

 	Incorrect water level, replace MR290 chamber	Rx only	Prescription only
 	Correct water level in the MR290 chamber	CE 0123	European Conformity
 	Consult operating instructions	 	Date and country of manufacture CC NZ: New Zealand MX: Mexico
	Single use		Manufacturer
	Lot number		Use-by date
	Reference number		Caution/Consult instructions for use
	Medical device		7 Days maximum use
	Importer		Transportation and storage temperature limits
	Distributor		European Union authorised representative
	UK responsible person		Switzerland authorised representative

WARNINGS, CAUTIONS, AND NOTES

WARNINGS

- There are residual risks associated with use of this product, even if used as intended. Following all instructions and warnings provided, the risks of hypoxic injury, skin burns, airway burns, airway injury, lung injury, electric shock injury, musculoskeletal injury, and hypothermia remain. These risks may result in serious injury or death.
- DO NOT reuse this product. Reuse may result in transmission of infectious substances, interruption to treatment, serious harm or death.
- The use of breathing circuits, chambers, accessories, or combinations which are not approved by Fisher & Paykel Healthcare may result in poor humidification system performance, ventilator malfunction, and harm to the patient/ user.
- Appropriate patient monitoring (e.g. oxygen saturation) must be used at all times. Failure to monitor the patient (e.g. in the event of an interruption to gas flow) may result in serious harm or death.
- DO NOT touch the heater-plate or chamber base. Surfaces may exceed 85 °C. Failure to comply may result in a skin burn.
- DO NOT bypass the pressure relief valve as this may cause damage/injury.

Failure to comply with the following warnings may impair performance of the device or compromise safety (including potentially cause serious harm):

- When mounting a humidifier adjacent to a patient, ensure that the humidifier is always positioned lower than the patient.
-  DO NOT use beyond 7 days maximum duration of use.
- DO NOT cover the circuit with materials such as blankets, towels or bed linen.

- DO NOT stretch or milk the tubing.
- DO NOT soak, wash, or sterilize this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.
- DO NOT use the chamber if the water level rises above the maximum water level line.
- DO NOT use the chamber if the seals are not intact when received, or if it has been dropped.
- DO NOT operate the chamber at an angle in excess of 10°.
- DO NOT spike the water source until the blue caps have been removed.
- Visually inspect breathing sets for damage (e.g. a crushed tube or cracked connector) before use, and replace if damaged.
- The water source must be at least 50 cm higher than the chamber.
- DO NOT use heated wire breathing circuits without gas flow. If gas flow is interrupted, turn the humidifier off.
- Ensure there is a water supply connected to the chamber and that water is present within the chamber.
- DO NOT fill the chamber with water in excess of 37 °C.
- DO NOT connect the pressure relief valve to the outlet port of the chamber as this may impair the operation of the valve.
- Only use pressure relief valve with blue cover with Nasal High Flow Therapy.
- Discard product if the blue cover on the pressure relief valve has been removed or is missing.
- Ensure any unused ports have their caps and/or plugs in place before use.
- Before attaching oxygen analyser to pressure relief valve, remove cap from oxygen analyser port.
- Discard product if the pressure relief valve is damaged or damage is suspected.
- DO NOT use an in-line nebulizer between the pressure relief manifold and the dry side of the humidification chamber.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

Attention

- Avoid prolonged contact of heated tubes with patient's skin.
- Use USP sterile water for inhalation or equivalent for humidification. DO NOT add other substances to the water.
- Check all connections are tight before use.
- DO NOT use with Fisher & Paykel Healthcare MR810 humidifier.
- Fill the chamber before connecting to flow source.

NOTES

- Ensure appropriate ventilator or flow source alarms are set before connecting breathing set to patient.
- For use under the supervision of trained medical personnel.
- Refer to hospital protocol for disposal methods. User may be exposed to breathing tract fluids during disposal.
- The responsible organization is accountable for the compatibility of the humidifier and all of the parts and accessories used to connect to the patient and other equipment before use.
- Notice to User: If a serious incident has occurred while using this device please inform your local Fisher & Paykel Healthcare representative and Competent Authority in your country.

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