






F&P 850 F&P Evaqua 2

The F&P Respiratory Care Continuum

Adult Evaqua™ 2 dual heated ventilator circuit & Optiflow™+ nasal cannula

USER INSTRUCTIONS REF RT481


CE 0123 Rx only

F&P, Evaqua and Optiflow are trademarks of Fisher & Paykel Healthcare Limited. For patent information, see www.fphcare.com/ip

Manufacturer Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.co.nz Web: www.fphcare.com **Importer/ Distributor Australia** (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 Fax: +61 3 9871 4998 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **Brazil** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo - SP, Brazil Tel: +55 11 2548 7002 **China** 代理人/售后服务机构 费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 传真: +86 20 32052132 **Denmark** Tel: +45 70 26 37 70 Fax: +46 83 66 310 **Finland** Tel: +358 9 251 66 123 Fax: +46 83 66 310 **France** EC REP Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 Email: c.s@fphcare.fr **Germany** Fisher & Paykel Healthcare GmbH, Deutschland, Österreich, Schweiz, Wiesenstrasse 49, D 73614 Schorndorf, Germany Tel: +49 7181 98599 0 Fax: +49 7181 98599 66 **Hong Kong** Tel: +852 2116 0032 Fax: +852 2116 0085 **India** Tel: +91 80 2309 6400 **Ireland** Tel: 1800 409 011 Fax: +44 1628 626 146 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea** Tel: +82 2 6205 6900 Fax: +82 2 6309 6901 **Mexico** Tel: +52 55 9130 1626 **Norway** Tel: +47 21 60 13 53 Fax: +47 22 99 60 10 **Poland** Tel: 00800 49 119 77 Fax: 00800 49 119 78 **Russia** Tel. and Fax: +7 495 782 21 50 **Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **Taiwan** Tel: +886 2 8751 1739 Fax: +886 2 8751 5625 **Turkey** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 Fax: +90 312 354 31 01 **UK** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 Fax: +44 1628 626 146 **USA/Canada** Tel: 1800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001

POST EXTUBATION KIT FOR RT481

ADULT EVAQUA 2 DUAL HEATED VENTILATOR CIRCUIT AND OPTIFLOW+ NASAL CANNULA

For delivery of breathing gases to adult patients requiring respiratory support. For use in an Intensive Care environment in conjunction with life support equipment.

TECHNICAL SPECIFICATIONS

Compatible with Fisher & Paykel Healthcare MR850 humidifiers.

Refer to humidifier user instructions.

INTERFACE CONNECTIONS	ISO 5356-1 Conical connectors
COMPRESSIBLE VOLUME	1.60 L
COMPLIANCE @ 60 cmH₂O	
Without Accessories	2.18 ± 0.11 mL/cmH ₂ O
With Accessories	2.6 ± 0.56 mL/cmH ₂ O
	(including 0.02 mL/cmH ₂ O measurement uncertainty)

RESISTANCE TO FLOW @ 30 L/min	
Without Accessories	
Inspiratory limb	0.77 ± 0.05 cmH ₂ O
Expiratory limb	0.20 ± 0.04 cmH ₂ O
With Accessories	
Expiratory limb	0.99 ± 0.11 cmH ₂ O
	(including 0.03 cmH ₂ O measurement uncertainty)




MAXIMUM OPERATING PRESSURE	8 kPa
BREATHING TUBE LENGTH	Inspiratory: 1.6 m Expiratory: 1.6 m
MINIMUM TUBE INTERNAL DIAMETER	20 mm
GAS LEAKAGE @ 60 cmH₂O	<40 mL/min
TIDAL VOLUME	>120 mL

OVERALL PERFORMANCE AT 20 °C TO 26 °C AMBIENT TEMPERATURE		
	INVASIVE MODE	NONINVASIVE MODE
Humidification Output	>33 mg/L	>12 mg/L
Flow Rate	5 – 60 L/min	5 – 120 L/min

FILTRATION EFFICIENCY

Viral Organism	>99.99%
Bacterial Organism	ΦX174 Bacteriophage >99.999%
Organism	Bacillus subtilis
Mean particle size	3 µm
NaCl	98.04%

SYMBOL DEFINITIONS

	Water level too high	Rx only	Prescription only
	Maximum water level	CE 0123	CE Marking 93/42/EEC
	Consult operating instructions		Date of manufacture
	Type BF applied part		Manufacturer
	Not made with phthalates		Use-by date
	Not made with natural rubber latex	14	14 Days maximum use
	Single use		Transportation and storage temperature limits
LOT	Lot number	EC REP	European Union authorised representative
REF	Reference number		Caution/Consult instructions for use

POST EXTUBATION KIT FOR USE WITH RT481
Nasal cannula patient interface for delivery of humidified respiratory gases.

CONTRAINDICATIONS

This therapy may not be suitable where continuous positive airway pressure (CPAP) is contraindicated. This includes non-spontaneous breathing.

Setup
MR850 Humidifier in invasive mode, RT481 kit with 22mm heated inspiratory tube and chamber, ventilator in high flow mode

Flow Range: <60 L/min

- Select appropriate size. Prongs must not create a seal in the nares. A clear gap must be visible around each prong.
- Adjust headstrap to fit. Do not over-tighten.
- Ensure headstrap clip is attached, to prevent cannula from being pulled out of the nares.

- Cannula can become detached if not used with the headstrap clip.
- Attach tubing clip to clothing/bedding to prevent cannula from pulling off face.
- Attach the tubing clip to the breathing circuit but ensure probe cable is not crushed by tubing clip.

WARNINGS, CAUTIONS AND NOTES

WARNINGS

- DO NOT reuse this product. Reuse may result in transmission of infectious substances, interruption to treatment, serious harm or death.
- The use of breathing circuits, chambers, accessories, or combinations which are not approved by Fisher & Paykel Healthcare may result in poor humidification system performance, ventilator malfunction and harm to the patient/ user.
- Appropriate patient monitoring (e.g. oxygen saturation) must be used at all times. Failure to monitor the patient (e.g. in the event of an interruption to gas flow) may result in serious harm or death.
- DO NOT touch heater plate or chamber base. Surfaces may exceed 85 °C. Failure to comply may result in a skin burn.
- Failure to comply with any of the following warnings may impair performance of the device or compromise safety (including potentially causing serious harm):**
- When mounting a humidifier adjacent to a patient, ensure that the humidifier is always positioned lower than the patient.

- 14** DO NOT use beyond 14 days maximum duration of use.
- DO NOT cover the circuit with materials such as blankets, towels or bed linen.
- DO NOT stretch or milk the tubing.
- DO NOT soak, wash, or sterilize this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.
- DO NOT use the chamber if the water level rises above the maximum water level line.
- DO NOT use the chamber if the seals are not intact when received, or if it has been dropped.
- DO NOT operate the chamber at an angle in excess of 10°.
- DO NOT spike the water source until the blue caps have been removed. Should the primary float fail, splashing into the circuit may occur if the chamber is being operated in excess of 80 L/min.

- Visually inspect breathing sets for damage (e.g. a crushed tube or cracked connector) before use and replace if damaged.
- The water source must be at least 50 cm higher than the chamber.
- DO NOT use heated wire breathing circuits without gas flow. If gas flow is interrupted, turn the humidifier off.
- Ensure there is a water supply connected to the chamber and that water is present within the chamber.
- DO NOT fill the chamber with water in excess of 37 °C.
- DO NOT use medications containing Tyloxapol (such as Tacholiquin) as this may damage the tubing and lead to a loss of ventilation pressure.
- Remove any sources of ignition; such as cigarettes or open flames.
- When nebulized drugs are used, resistance to flow should be monitored and the filter replaced, following standard hospital procedure.

- Change filter every 24 hours or sooner if noticeable deterioration/build-up occurs, following standard hospital procedure.
- Ensure appropriate ventilator or flow source alarms are set before connecting breathing set to patient.
- Avoid prolonged contact of heated tubes with patient’s skin.
- Perform a pressure and leak test on the breathing system and check for occlusions before connection to a patient.
- Use USP Sterile Water For Inhalation or equivalent for humidification. DO NOT add other substances to the water.
- Check breathing circuits for condensation every 6 hours and drain if required.
- Check all connections are tight before use.
- Check that the heater wire is evenly distributed along the circuit and not bunched or kinked.
- DO NOT use if packaging is not sealed.
- Use only the MR290HFV with the SensorMedics 3100B.
- Before connecting the interface, check for adequate gas flow and ensure that the system has warmed up.
- Failure to use the set-up described above can compromise performance and affect patient safety.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65.

NOTES

- For use under the supervision of trained medical personnel.
- Dispose of product according to appropriate hospital protocol. User may be exposed to breathing tract fluids during disposal.
- The responsible organization is accountable for the compatibility of the humidifier and all of the parts and accessories used to connect to the patient and other equipment before use.

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