# Surgical

## Humidification System

### SH870

#### INSTRUCTIONS FOR USE

REF SH870AEA

REF SH870AEK

REF SH870AFU

REF SH870JHU

# Surgical

## Humidification System

|      |                   |    |
|------|-------------------|----|
| en   | ENGLISH           | 01 |
| de   | DEUTSCH           | 06 |
| es   | ESPAÑOL           | 11 |
| fi   | SUOMI             | 16 |
| fr   | FRANÇAIS          | 21 |
| frca | FRANÇAIS CANADIEN | 26 |
| it   | ITALIANO          | 31 |
| ko   | 한국어               | 36 |
| nl   | NEDERLANDS        | 41 |
| pt   | PORTUGUÊS         | 46 |
| sv   | SVENSKA           | 51 |
| zht  | 繁體中文              | 56 |

## INTENDED USE

To warm and humidify carbon dioxide (CO<sub>2</sub>) from a gas control device prior to entry to the surgical cavity.

### Patient Profile

Patients of all ages undergoing laparoscopic, open abdominal or open cardiothoracic surgical procedures where (CO<sub>2</sub>) insufflation can be used.

### Operator Profile

To be operated by a healthcare professional with experience in the set-up of surgical equipment.


### General Precautions

- Only use devices and accessories approved by Fisher & Paykel Healthcare. Unauthorized devices and accessories may impair performance or compromise patient safety.
- Monitor patient temperature throughout the procedure and adjust patient warming protocols to prevent hyperthermia or hypothermia.

### Warning

- Do not immerse in liquid, to avoid the risk of electric shock.
- Connect to a supply mains with a protective earth, to avoid risk of electric shock.
- To avoid the risk of Electromagnetic interference, avoid stacking humidifier on other equipment. See Page 4 for appropriate distances from other equipment.
- Portable RF communications equipment including antennas, can effect medical electrical equipment and should therefore be used no closer than 12 cm to any part of the SH870 System including cables.
- Other cables and accessories may negatively affect Electromagnetic Compatibility performance.
- The operation of high frequency surgical apparatus, short wave or microwave equipment in the vicinity of the humidifier may adversely affect its function. If this occurs, the humidifier should be removed from the vicinity of such devices.

## CONTROLS AND INDICATORS

- 1. Standby/On.** Pressing this button turns the humidifier on or stand-by
- 2. On Indicator.** A solid light indicates that the humidifier is on.
- 3. Heater-wire Adaptor Plug Indicator.** A flashing light indicates that the heater-wire plug is disconnected.
- 4. Flow Out/Water Out Indicator.** A flashing light indicates that the chamber requires refilling and/or the gas supply is running low. Check both the chamber water level and gas supply.
- 5. Hardware Fault Indicator.** A flashing or solid light indicates a fault with the humidifier and it therefore requires servicing. Please refer to the SH870 Product Technical Manual (Ref 185048422). Humidified gas will not be delivered when this indicator is on.

## TO SETUP

1. Remove the humidifier from its packaging. Plug the power cord into the power supply and ensure the plug is accessible at all times.
2. Ensure that the humidifier is level and securely mounted. A range of mounting solutions are available from Fisher & Paykel Healthcare. Please contact your Fisher & Paykel Healthcare representative or authorized distributor for assistance in selecting the most suitable mounting bracket.
3. At the end of the procedure ensure the humidifier is on standby by pressing the "Standby/On" button.

## ROUTINE MAINTENANCE AND SERVICING

In order to keep your humidifier in good working order, it is necessary for qualified service personnel to perform maintenance at regular intervals. A full technical description including maintenance and service data, is contained in the SH870 Product Technical Manual (Ref 185048422) available from your Fisher & Paykel Healthcare office or authorized distributor.

### Regular Maintenance Period

12 months

### Maintenance Aspects\*

Check for damage. Carry out a performance test

\* Refer to the SH870 Product Technical Manual for detailed instructions (Ref 185048422).

## GENERAL CLEANING INSTRUCTIONS

### Surgical Humidifier

1. Disconnect the humidifier from any electrical outlet
2. Using a damp cloth, clean the humidifier with one of the following solutions:
  - Isopropyl alcohol
  - Dishwashing detergent
3. Wipe the humidifier clean of any cleaning residues before use.

## SPECIFICATIONS

### General Specifications

| Specifications for the F&P HumiGard™ Surgical Humidifier |                                     |
|--|-------------------------------------|
| Mode of operation  | Continuous                          |
| Dimensions (without chamber fitted)                      | 94 x 154 x 135 mm                   |
| Weight (without chamber fitted)                          | 1.7 kg                              |
| Supply voltage   | 230 V- (SH870A--) 115 V- (SH870J--) |
| Supply frequency   | 50/60 Hz                            |
| Power input  | 200 VA                              |
| Heater-wire  | 22 V ~ at 30 W                      |
| Heater-plate   | 150 W                               |
| Heater-plate high temperature                            | 93 ± 6 °C cut-out                   |
| Operating environment                                    | 18 °C to 26 °C                      |
| Humidity performance (0 to 10 L/min)                     | >33 mg/L                            |
| Transport and storage conditions                         | -10 °C to 50 °C                     |
| Classification   | Class 01: Type BF Applied Part      |

| Specifications for Laparoscopic Surgery<br>with ST320 Humidified Insufflation Kit | Specifications for Open Surgery<br>with ST320 Humidified Insufflation Kit |
|---|---|
| Input set flow rate: ≤ 45 L/min   | Recommended flow rate: 10 L/min   |

Refer to the ST320 user instructions (Ref 185048421) for setup. **NOTE:** Not made with natural rubber latex.

## GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC EMISSIONS

The SH870 Surgical Humidification system is intended for use in the electromagnetic environment specified below. The customer or the user of the SH870 Surgical Humidification system should ensure that it is used in such an environment.

| EMISSIONS TEST  | COMPLIANCE | ELECTROMAGNETIC ENVIRONMENT – GUIDANCE   |
|---|------------|--|
| RF EMISSIONS<br>CISPR 11                                | GROUP 1    | The SH870 Humidifier uses RF energy only for its internal function. Therefore, its RF emissions are very low and not likely to cause any Interference in nearby electronic equipment.                              |
| RF EMISSIONS<br>CISPR 11                                | CLASS A    |  |
| Harmonic emissions<br>IEC 61000-3-2                     | CLASS A    | The SH870 Humidifier is suitable for use in all establishments, including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Voltage fluctuations/flicker<br>emissions IEC 61000-3-3 | Complies   |  |

Note: The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). If it used in a residential environment (for which CISPR 11 Class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orientating the equipment.

## GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY

The SH870 Surgical Humidifier is intended for use in the electromagnetic environment specified below. The customer or the user of the SH870 Surgical Humidifier should ensure that it is used in such an environment.

| IMMUNITY TEST  | IEC 60601<br>TEST LEVEL   | COMPLIANCE LEVEL  | ELECTROMAGNETIC ENVIRONMENT – GUIDANCE   |
|--|---|---|--|
| Electrostatic<br>discharge (ESD)<br>IEC 61000-4-2  | ± 8 kV contact<br>± 15 kV air   | ± 8 kV contact<br>± 15 kV air   | Floors should be made of wood, concrete or ceramic tile. If floors are covered with synthetic material, the Relative Humidity should be at least 30%.  |
| Electrical fast<br>transient/burst<br>IEC 61000-4-4  | ± 2 kV for power<br>supply lines<br>± 1 kV for input/output lines   | ± 2 kV for power<br>supply lines<br>± 1 kV for input/output lines   | Mains-power quality should be that of a typical commercial or hospital environment.  |
| Surge<br>IEC 61000-4-5   | ± 1 kV line(s) to line(s)<br>± 2 kV line(s) to earth  | ± 1 kV line(s) to line(s)<br>± 2 kV line(s) to earth  | Mains-power quality should be that of a typical commercial or hospital environment.  |
| Voltage dips, short<br>interruptions and voltage<br>variations on power supply<br>input lines IEC 61000-4-11 | <5% $U_t$ (>95% dip in $U_t$<br>for 0.5 cycle)<br>40% $U_t$ (60% dip in $U_t$<br>for 5 cycles)<br>70% $U_t$ (30% dip in $U_t$<br>for 25 cycles)<br><5% $U_t$ (>95% dip in $U_t$<br>for 5 seconds) | <5% $U_t$ (>95% dip in $U_t$<br>for 0.5 cycle)<br>40% $U_t$ (60% dip in $U_t$<br>for 5 cycles)<br>70% $U_t$ (30% dip in $U_t$<br>for 25 cycles)<br><5% $U_t$ (>95% dip in $U_t$<br>for 5 seconds) | Mains-power quality should be that of a typical commercial or hospital environment. If the user of the SH870 Surgical Humidifier requires continued operation during power mains usage, it is recommended that the SH870 Surgical Humidifier be powered from an uninterruptible power supply or battery. |
| Power frequency<br>(50/60 Hz) magnetic field<br>IEC 61000-4-8  | 30 A/m  | 30 A/m  | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.  |

Note:  $U_t$  is the A.C. mains voltage prior to application of the test level.

## GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY

The SH870 Surgical Humidifier is intended for use in the electromagnetic environment specified below. The customer or the user of the SH870 Surgical Humidifier should ensure that it is used in such an environment.

| IMMUNITY TEST                 | IEC 60601 TEST LEVEL        | COMPLIANCE LEVEL | ELECTROMAGNETIC ENVIRONMENT – GUIDANCE   |
|-------------------------------|-----------------------------|------------------|--|
| Conducted RF<br>IEC 61000-4-6 | 3 Vrms<br>150 kHz to 80 MHz | 3 V/m            | Portable and mobile RF communications equipment should be used no closer to any part of the SH870 Surgical Humidifier, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.<br><b>Recommended separation distance</b><br>$d = 1.2\sqrt{P}$<br>$d = 1.2\sqrt{P} \text{ 80 MHz to 800 MHz}$<br>$d = 2.3\sqrt{P} \text{ 800 MHz to 2.5 GHz}$   |
| Radiated RF<br>IEC 61000-4-3  | 10 V/m<br>80 MHz to 2.7 GHz | 10 V/m           | Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site surveyor, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol.<br> |

**Note 1.** At 80 MHz and 800 MHz, the higher-frequency range applies. **Note 2.** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast, cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the SH870 Surgical Humidifier is used exceeds the applicable RF compliance level above, the SH870 Surgical Humidifier should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the SH870 Surgical Humidifier.
- b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

## RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE SH870 HUMIDIFIER

The SH870 Surgical Humidification System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the SH870 Humidifier can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the SH870 Humidifier as recommended below, according to the maximum output power of the communications equipment.

| RATED MAXIMUM OUTPUT POWER OF TRANSMITTER (W) | SEPARATION DISTANCE ACCORDING TO THE FREQUENCY OF TRANSMITTER (m) |  |   |
|---|---|--|---|
|   | 150 MHZ TO 80 MHZ<br>$d = 1.2\sqrt{P}$                            | 80 MHZ TO 800 MHZ<br>$d = 1.2\sqrt{P}$ | 800 MHZ TO 2.5 GHZ<br>$d = 2.3\sqrt{P}$ |
| 0.01  | 0.12  | 0.12                                   | 0.23                                    |
| 0.1   | 0.38  | 0.38                                   | 0.73                                    |
| 1   | 1.2   | 1.2                                    | 2.3                                     |
| 10  | 3.8   | 3.8                                    | 7.3                                     |
| 100   | 12  | 12                                     | 23                                      |

## SYMBOL

**IPX1** Resistant to vertical falling drips

STANDBY/ON Button

Fragile. Handle With Care

On Indicator

This Way Up

Heater-wire Adaptor Plug Indicator

Regulatory Compliance Mark

Flow Out/Water Out Indicator

Type BF Applied Part

Hardware Fault Indicator

Follow Instructions For Use

Electrostatic Sensitive Devices

## DEVICE DISPOSAL INSTRUCTIONS

The SH870 Humidifier is intended for repeat use. If disposal is required, please do not discard as regular rubbish. Dispose of electronics according to local guidelines.

## APPROVED ACCESSORIES

### Laparoscopic Surgery Accessories

- ST320 Humidified Insufflation Kit

### Open Surgery Accessories

- ST320 Humidified Insufflation Kit
- ST300DF VITA-diffuser™ manufactured by Cardia Innovation AB

**NOTE:** Accessories may not be available in all countries. Contact your local Fisher & Paykel Healthcare office or authorized distributor for more information.

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For more information, please contact your local  
Fisher & Paykel Healthcare representative

**Manufacturer**

Fisher & Paykel Healthcare  
Ltd,  
15 Maurice Paykel Place,  
East Tamaki, Auckland 2013  
PO Box 14 348 Panmure,  
Auckland 1741,  
New Zealand

Tel: +64 9 574 0100  
Fax: +64 9 574 0158  
Email: info@fphcare.co.nz  
Web: www.fphcare.com  
Importer/ Distributor  
**Australia (Sponsor)**  
Fisher & Paykel Healthcare Pty  
Ltd,  
19-31 King Street, Nunawading,  
Melbourne, Victoria 3131.  
Tel: +61 3 9871 4900  
Fax: +61 3 9871 4998

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Tel: +90 312 354 34 12  
Fax: +90 312 354 31 01

**UK**

Fisher & Paykel Healthcare Ltd,  
Unit 16, Cordwallis Park,  
Clivemont Road, Maidenhead,  
Berkshire SL6 7BU, UK

Tel: 0800 132 189  
Fax: +44 1628 626 146

**USA/Canada**

Tel: 1800 446 3908  
or +1 949 453 4000  
Fax: +1 949 453 4001

€ 0123

REF 185048419 REV F 2020-08 © 2020 Fisher &amp; Paykel Healthcare Limited

For patent information, see [www.fphcare.com/ip](http://www.fphcare.com/ip).

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İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi,