

This guide is intended for multi-patient use of the F&P Solo™ mask in the hospital, or other clinical settings, where proper disinfection or sterilization of the device must occur between patients. When using the mask for single-patient use in the home, please refer to your F&P Solo Use and Care Guide under ‘Cleaning Your Mask at Home’.

Mask Components	Thermal Disinfection		
	In line with EN ISO 15883, the uncertainty range for disinfection temperature is -0 °C and +5 °C		
	1 min at 90 °C (194 °F)	10 mins at 80 °C (176 °F)	30 mins at 75 °C (167 °F)
F&P Solo Mask Headgear	Replace these parts with new parts between patients		
<ul style="list-style-type: none"> <li>- Frame</li> <li>- Tube</li> <li>- Swivel</li> <li>- Nasal Cushion*</li> <li>- Pillows Cushion*</li> </ul>	-	Maximum 25 Cycles	Maximum 25 Cycles

\*Replace cushion if it discolors.  
Each cycle is limited to one day of use.

### ⚠ General Warnings and Cautions

The guidelines provided in this guide have been validated as being capable of preparing the masks for re-use. Fisher & Paykel Healthcare does not recommend any deviations from the recommended method of reprocessing.	
Warning	Deviations from the recommended methods of reprocessing may result in inadequate disinfection/sterilization and adverse health effects.
Warning	Frequency of cleaning, methods of cleaning or the use of cleaning agents, other than those specified in the accompanying documents or exceeding the number of processing cycles, can have an adverse effect on the mask and consequently the safety or the quality of the therapy.
Caution	Only clean the mask with cleaning products specified in the cleaning guidelines provided in these instructions, as any other products may result in an adverse reaction if cleaning residue is left behind.
Caution	Do not clean the mask with products containing alcohol, antiseptic, bleach, chlorine, or moisturizer as these may compromise the function of the mask.
Caution	Do not clean the mask in a dishwasher as this may compromise the function of the mask.

## Validated Disinfection and Sterilization Procedures

<b>Point of Use</b>	Remove any gross soiling on the mask by wiping or rinsing in potable water.
<b>Disassembly</b>	Disassemble the mask according to the Disassembling your Mask guidelines provided in the F&P Solo Use and Care Guide. <b>⚠ Caution:</b> Failure to fully disassemble the mask may compromise the effectiveness of cleaning and the subsequent disinfection steps which may result in a range of harms, some of which are severe may compromise the effectiveness of cleaning and the subsequent disinfection steps which may result in a range of harms, some of which are severe.
<b>Cleaning</b>	Masks should always be properly cleaned before further reprocessing. <b>⚠ WARNING:</b> Failure to properly clean this device may result in ineffective disinfection/sterilization and adverse health effects.
	1. Soak all reprocessable parts in a solution of mild alkaline anionic detergent according to the manufacturer's instructions (e.g. Alconox™ <sup>1</sup> at 1% concentration and 50 °C) for 5 minutes, ensuring that no air bubbles are present on the device surface.
	2. Using a soft non-metallic brush (e.g. medium-hard toothbrush), vigorously brush the mask parts until visibly clean but not for less than 30 seconds. An additional narrow bottle brush is required to access the Mask Tube. Brush the components while submerged and out of the solution. Pay close attention to crevices, cavities and corners, such as the join between the silicone cushion and plastic surrounds and the internal surfaces of the frame.
	3. Vigorously brush the joint between the frame and tube alternating between the inside and outside surfaces for at least 30 seconds.
	4. Manually rotate the swivel for at least 10 seconds, while submerged in the detergent solution.
	5. Repeat steps 1 to 4 (total soaking time of 10 minutes).
	6. Rinse by submerging in 5 liters of fresh demineralised water and agitating for at least 10 seconds. Repeat rinse.
	7. Manually rotate the swivel for at least 10 seconds, while submerged in the demineralised water.
	8. Ensure the mask is visibly clean before continuing. If not, repeat steps 1-7.
<b>Drying</b>	Make sure mask components are thoroughly dry before continuing by leaving them to air dry out of direct sunlight. Failure to do this may result in ineffective disinfection.
<b>Disinfection or Sterilization</b>	<b>Thermal Disinfection</b>
	<b>For all reprocessable components:</b>
	Run a disinfection cycle in a legally marketed automatic washer-disinfector using one of the time/temperature combinations shown in the table above. These time/temperature combinations are consistent with EN ISO 15883 and are equivalent in terms of known thermal inactivation kinetics of vegetative microorganisms, meeting an A <sub>0</sub> of 600
<b>Caution</b>	The use of STERRAD™ <sup>2</sup> or Autoclave as a reprocessing method may lead to mask degradation and adverse health effects.
<b>Drying</b>	Allow all components to air dry out of direct sunlight.
<b>Inspection</b>	Before each use, inspect the mask for deterioration. If any visible deterioration of a mask component is apparent, the mask component should be discarded and replaced.
<b>Reassembly</b>	Reassemble the mask according to the Assembling your Mask guidelines provided in the F&P Solo Use and Care Guide. To ensure traceability, mask components belonging to an individual mask must be reassembled back into the same mask after cleaning, excluding components requiring replacement as marked in the table above.
<b>Storage</b>	Ensure the mask is completely dry before storing. Store in clean conditions out of direct sunlight. Storage temperature: -20 °C to 50 °C (-4 °F to 122 °F).

1. Alconox is a trademark of Alconox, Inc.

2. Sterrad is a registered trademark of ASP Global Manufacturing GmbH.

F&P and F&P Solo are trademarks of Fisher & Paykel Healthcare Ltd. For patent information, see [www.fphcare.com/ip](http://www.fphcare.com/ip).

SUI 631076 REV A 2023-10  
© 2023 Fisher & Paykel Healthcare Limited



Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Email: [info@fphcare.co.nz](mailto:info@fphcare.co.nz) Web: [www.fphcare.com](http://www.fphcare.com)

**Australia (AU)** (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 **Brazil (BR)** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo – SP, Brazil Tel: +55 11 2548 7002 **China (CN)** 代理人/售后服务机构: 费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 **France (FR)** Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Email: [c.s@fphcare.fr](mailto:c.s@fphcare.fr)  
**Germany (DE)** Fisher & Paykel Healthcare GmbH, Wiesenstrasse 49, 73614 Schorndorf, Germany Tel: +49 7181 98599 0 **Hong Kong (HK)** Tel: +852 2116 0032 **India (IN)** Tel: +91 80 2309 6400 **Japan (JP)** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea (KR)** Tel: +82 2 6205 6900 **Mexico (MX)** Tel: +52 55 9130 1626  
**Poland (PL)** Fisher & Paykel Healthcare Poland Sp. z o.o., Pl. Andersa 7, 61-894 Poznań, Poland Tel: +48 664 846 464 **Russia (RU)** Tel: +7 495 782 21 50 **Switzerland (CH)** Fisher & Paykel Healthcare GmbH, Säntisstrasse 2, 9501 Wil / SG, Switzerland Tel: 0800 83 47 63 **Taiwan (TW)** Tel: +886 2 8751 1739  
**Turkey (TR)** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 3412 **UK (GB)** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 **USA (US)/Canada (CA)** Tel: 1800 446 3908 or +1 949 453 4000

**Austria (AT)** Tel: 0800 29 31 23 **Benelux (BE NL LU)** Tel: +31 40 216 3555 **Denmark (DK)** Tel: +45 70 26 37 70 **Finland (FI)** Tel: +358 9 251 66 123 **Ireland (IE)** Tel: 1800 409 011 **Italy (IT)** Tel: +39 06 7839 2939 **Norway (NO)** Tel: +47 21 60 13 53 **Spain (ES)** Tel: +34 902 013 346 **Sweden (SE)** Tel: +46 8 564 76 680