


# Choosing the right F&P Solo Pillows cushion size

## Start with the Medium cushion

The Medium cushion fits the majority of users. Follow the fitting tips, turn CPAP\*-therapy pressure on and lie down in bed.



**If you experience persistent leaks,**  
try the Large cushion.



**If you experience discomfort**  
in your nostrils, try the Small cushion.



\*We refer to CPAP in this instance, but this may also include other positive airway pressure device modes such as APAP or BiLevel. Masks must be used with a PAP device to deliver therapy.

Always read the label and follow the directions for use. Prescription only. This product may not be right for you. Your healthcare provider will advise you whether this product is suitable for you/your condition. If symptoms persist, talk to your healthcare provider. F&P, F&P Solo and Solo are trademarks of Fisher & Paykel Healthcare Limited. For patent information, please see [www.fphcare.com/ip](http://www.fphcare.com/ip). 633086 REV A © 2024 Fisher & Paykel Healthcare Limited.