Humidified Insufflation Kit

ST320

INSTRUCTIONS FOR USE
LAPAROSCOPIC AND OPEN SURGERY

REF ST320

ST320

INSTRUCTIONS FOR USE LAPAROSCOPIC AND OPEN SURGERY

en	ENGLISH	1
de	DEUTSCH	6
es	ESPAÑOL	11
fi	SUOMI	16
fr	FRANÇAIS	21
frca	FRANÇAIS CANADIEN	26
it	ITALIANO	31
ko	한국어	36
nl	NEDERLANDS	41
pt	PORTUGUÊS	46
sv	SVENSKA	51
zht	繁體中文	56

ADDITIONAL COMPONENTS

BARBED CONNECTOR
GERIFFELTER KONNEKTOR
CONECTOR DE LENGÜETAS
VÄKÄPÄINEN LIITIN
RACCORD CANNELÉ
RACCORD CANNELÉ
CONNETTORE UNCINATO
바브 커넥터
GERIBBELDE CONNECTOR
CONECTOR FARPADO
HULLINGFÖRSEDD ANSLUTNING
插管式連接頭


BYPASSCONNECTOR
BYPASS-KONNEKTOR
CONECTOR DE DERIVACIÓN
OHITUSLIITIN
RACCORD DE DÉRIVATION
RACCORD DE DÉRIVATION
CONNETTORE DI BYPASS
바이패스커넥터
BYPASS-CONNECTOR
CONECTOR DE DERIVAÇÃO
FÖRBIKOPPLINGSANSLUTNING
旁通連接頭


Humidified Insufflation Kit

ST320 HUMIDIFIED INSUFFLATION KIT

SETUP FOR LAPAROSCOPIC & OPEN SURGERY

TO SET UP

- 01.** Remove the pre-assembled, non-sterile components from the tray. Depress the finger guard on the humidifier and slide on pre-assembled components (chamber and dry line).
- 02.** Remove the filter end from the chamber.
- 03.** Connect filter end to gas source, using the barbed connector if required;
 - a. laparoscopic surgery: insufflator
 - b. open surgery: (carbon dioxide) CO₂ gas supply stand.
- 04.** Fill chamber with sterile water to level indicator (~180 mL) using funnel.
- 05.** Remove the sterile insufflation tube from the pouch.
- 06.** Connect the heater wire adapter to the insufflation tube elbow.
- 07.** Push the insufflation tube elbow onto the chamber.
- 08.** Connect the luer lock to the patient interface;
 - a. laparoscopic surgery: cannula
 - b. open surgery: VITA-diffuser™.
- 09.** Turn on the humidifier.

INTENDED USE

The ST320 Humidified Insufflation Kit is an accessory to the SH870 Surgical Humidifier, intended to warm and humidify carbon dioxide (CO₂) from a gas control device prior to entry to the surgical cavity.

PATIENT PROFILE

Patients of all ages undergoing laparoscopic, open abdominal or open cardiothoracic surgical procedures where CO₂ insufflation can be used.

OPERATOR PROFILE

To be operated by a healthcare provider with experience in the set-up of surgical equipment.

INSTRUCTIONS FOR USE

WARNINGS

General Warnings

- Ensure only CO₂ gas is used. Misconnection of flammable gases could result in serious harm.
- Monitor for systemic acidosis. CO₂ insufflation may increase acidosis.
- Intended for single use only. Reuse may result in transmission of infectious substances, interruption of treatment, serious harm or death.

Warnings for Laparoscopic Procedures

- Purge the system of air using CO₂ before insufflation, to prevent air embolism.

Warnings for Open Surgery Procedures

- Insufficient insufflation may result in inefficient de-airing; monitor for gas embolism in cardiac surgery.
- In-line gas monitoring is required to avoid systemic acidosis in cardiac surgery.

PRECAUTIONS

General Precautions

- Inspect sterile packaging for damage before use. Discard if damaged.
- Monitor patient temperature throughout the procedure and adjust patient warming protocols to prevent hyperthermia or hypothermia.
- If using implantables, please consult the manufacturer for suitability for use with humidified CO₂.
- Only use devices and accessories approved by Fisher & Paykel Healthcare. For correct use of approved accessories and devices, follow associated instructions for use.
- Monitor to ensure there is always water in the chamber.

Precautions for Open Surgery Procedures

- Ensure appropriate ventilation in the operating area to reduce risk of CO₂ related fatigue.

TROUBLESHOOTING: VERESS NEEDLE USE

The use of flow rates >1 L/min when using a Veress needle may result in unstable and inaccurate insufflator pressure readings.

To ensure stable and accurate pressure readings on the insufflator when using the Veress needle for flow rates >1 L/min, the following guidance applies:

- Ensure that the humidification chamber is filled to the recommended water level of ~180 mL, or bypass the chamber using the supplied by pass connector while the Veress needle is in use.
- Make sure insufflators are regularly serviced and/or calibrated in line with the manufacturer's instructions.

APPROVED DEVICES AND ACCESSORIES

- SH870 F&P HumiGard Surgical Humidifier

Open Surgery Accessory only

- ST300DF VITA-diffuser™ manufactured by Cardia Innovation AB

NOTE: Accessories may not be available in all countries. Contact your local F&P Healthcare office or authorized distributor for more information.

SPECIFICATIONS


General Specifications

- Humidification chamber reservoir: ~180 mL sterile water
- Luer fittings conform to ISO 594-2:1998
- Not made with natural rubber latex

F&P Healthcare has a policy of continuous product improvement and reserves the right to alter specifications without notice.

Specifications for Laparoscopic Surgery	Specifications for Open Surgery
Maximum input set flow rate: ≤ 45 L/min	Recommended flow rate: 10 L/min
Maximum achieved flow rate: ≤ 15 L/min	
Maximum input set pressure: ≤ 30 mmHg	

SYMBOLS

	Sterile using irradiation		Packaging pouch is recyclable
	Non-sterile		Single use
	Excess water in the chamber		Tray is recyclable
	Correct water level in the chamber		Use-by date
	Inspect sterile packaging for damage before use. Discard if damaged		Date of manufacture
	Peel sterile bag open in the direction indicated		Phthalates Free (PHT, DEHP, BBP, DBP Free)

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For more information, please contact your
local Fisher & Paykel Healthcare representative

Manufacturer

Fisher & Paykel Healthcare Ltd,
15 Maurice Paykel Place,
East Tamaki, Auckland 2013

PO Box 14 348 Panmure,
Auckland 1741,
New Zealand

Tel: +64 9 574 0100
Fax: +64 9 574 0158
Email: info@fphcare.co.nz
Web: www.fphcare.com

Importer/ Distributor

Australia (Sponsor)

Fisher & Paykel Healthcare Pty
Ltd,
19-31 King Street, Nunawading,
Melbourne, Victoria 3131.

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Tel: +90 312 354 34 12
Fax: +90 312 354 31 01

UK

Fisher & Paykel Healthcare Ltd,
Unit 16, Cordwallis Park,
Clivemont Road, Maidenhead,
Berkshire SL6 7BU, UK

Tel: 0800 132 189
Fax: +44 1628 626 146

USA/Canada

Tel: 1800 446 3908
or +1 949 453 4000
Fax: +1 949 453 4001

CE0123

REF 185048421 REV D 2020-02 © 2020 Fisher & Paykel Healthcare Limited

For patent information, see www.fphcare.com/ip

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Ürünleri Ticaret Limited Şirketi,