



USER INSTRUCTIONS

REF RT023 Patient Ventilator Weaning Kit

CE 0123 Rx only MD -10 °C 50 °C

F&P is a trademark of Fisher & Paykel Healthcare Limited.

Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Email: info@fphcare.co.nz Web: www.fphcare.com
Australia (AU) (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 Brazil (BR) Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo – SP, Brazil Tel: +55 11 2548 7002 China (CN)
代理人/售后服务机构:费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 France (FR) [EC REP] Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France
Tel: +33 1 6446 5201 Email: c.s.@fphcare.fr Germany (DE) Fisher & Paykel Healthcare GmbH, Wiesenstrasse 49, 73614 Schorndorf, Germany Tel: +49 7181 98599 0 Hong Kong (HK) Tel: +852 2116 0032 India (IN) Tel: +91 80 2309 6400 Japan (JP) Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 Korea (KR)
Tel: +82 2 6205 6900 Mexico (MX) Tel: +52 55 9130 1626 Poland (PL) Fisher & Paykel Healthcare Poland Sp. z o.o., Pl. Andersa 7, 61-894 Poznań, Poland Tel: +48 664 846 464 Russia (RU) Tel: +7 495 782 21 50 Switzerland (CH) [CH REP] Fisher & Paykel Healthcare GmbH, Säntisstrasse 2,
9501 Wil / SG, Switzerland Tel: 0800 83 47 63 Taiwan (TW) Tel: +886 2 8751 1739 Turkey (TR) İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Caddesi No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 UK (GB) [UK REP]
Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 USA (US)/Canada (CA) Tel: 1800 446 3908 or +1 949 453 4000

Austria (AT) Tel: 0800 29 31 23 Benelux (BE NL LU) Tel: +31 40 216 3555 Denmark (DK) Tel: +45 70 26 37 70 Finland (FI) Tel: +358 9 251 66 123 Ireland (IE) Tel: 1800 409 011 Italy (IT) Tel: +39 06 7839 2939 Norway (NO) Tel: +47 21 60 13 53 Spain (ES) Tel: +34 902 013 346 Sweden (SE)
Tel: +46 8 564 76 680

English

PATIENT VENTILATOR WEANING KIT

For delivery of breathing gases

TECHNICAL SPECIFICATIONS

Compatible with Fisher & Paykel Healthcare MR850 humidifiers and Adult RT series breathing circuits.

Refer to the MR850 Respiratory Humidifier instructions for use for additional safety information, clinical claims and instructions for setup and use.

Unless otherwise stated, flow rates are expressed at BTPS (body temperature and pressure, saturated) conditions.

INTERFACE CONNECTIONS	ISO 5356-1 Conical connectors
COMPLIANCE @ 60 cmH₂O	
Without Accessories	0.68 ± 0.21 mL/cmH ₂ O (including 0.02 mL/cmH ₂ O measurement uncertainty)
RESISTANCE TO FLOW @ 30 L/min	
Without Accessories	
Inspiratory limb:	0.11 ± 0.02 cmH ₂ O (including 0.03 cmH ₂ O measurement uncertainty)
MAXIMUM OPERATING PRESSURE	Refer to chamber instructions for use
BREATHING TUBE LENGTH	Inspiratory: 0.6 m
GAS LEAKAGE @ 60 cmH₂O	<1.2 L/min

SYMBOL DEFINITIONS			
	Consult operating instructions	Rx only	Prescription only
	Single use	CE 0123	European Conformity
	Lot number		Date of manufacture and ISO country code of manufacture (CC)
	Reference number		Manufacturer
	7 Days maximum use		Use-by date
	Medical device		Transportation and storage temperature limits
	UK responsible person		European Union authorised representative
	Switzerland authorised representative		Distributor
	Importer		Caution/consult instruction for use

WARNINGS, CAUTIONS AND NOTES

WARNINGS

- There are residual risks associated with use of this product, even if used as intended. Following all instructions and warnings provided, the risks of hypoxic injury, skin burns, airway burns, airway injury, lung injury, electric shock injury, musculoskeletal injury, and hypothermia remain. These risks may result in serious injury or death.
- DO NOT reuse this product. Reuse may result in transmission of infectious substances, interruption to treatment, serious harm or death.
- The use of breathing circuits, chambers, accessories, or combinations which are not approved by Fisher & Paykel Healthcare may result in poor humidification system performance, ventilator malfunction and harm to the patient/user.

Failure to comply with any of the following warnings may impair performance of the device or compromise safety (including potentially causing serious harm):

- When mounting a humidifier adjacent to a patient, ensure that the humidifier is always positioned lower than the patient.
-  DO NOT use beyond 7 days maximum duration of use.
- DO NOT cover with materials such as blankets, towels or bed linen.
- DO NOT stretch or milk the tubing.
- DO NOT soak, wash, or sterilize this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.
- Check all connections are tight before use.
- DO NOT use CPAP valve with an air entrainer.
- Make sure CPAP generator can generate pressure required for CPAP valve operation.
- DO NOT use in a closed loop breathing circuit.
- Inspiratory swivel can disconnect under excess strain.
- Check all components are intact.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

NOTES

- For use under the supervision of trained medical personnel.
- Dispose of product according to appropriate hospital protocol. User may be exposed to breathing tract fluids during disposal.
- The responsible organization is accountable for the compatibility of the humidifier and all of the parts and accessories used to connect to the patient and other equipment before use.
- Notice to User: If a serious incident has occurred while using this device, please inform your local Fisher & Paykel Healthcare representative and Competent Authority in your country.

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