
USER INSTRUCTIONS

REF 850A61 F&P 850 AirSpiral™ Adult NIV and NHF circuit kit

CE 0123 Rx only

F&P and AirSpiral are trademarks of Fisher & Paykel Healthcare Limited.
For patent information, see www.fphcare.com/ip

www.fphcare.com

Fisher & Paykel
HEALTHCARE

REF 901218 REV C 2021-09 © 2021 Fisher & Paykel Healthcare Limited

Manufacturer ■ Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.co.nz Web: www.fphcare.com Importer/ Distributor **Australia** (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131, Tel: +61 3 9871 4900 Fax: +61 3 9871 4998 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **Brazil** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraiso, 04004-000, São Paulo - SP, Brazil Tel: +55 11 2548 7002 **China** 代理人/售后服务机构:费雪派克医疗保健(广州)有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 传真: +86 20 32052132 **Denmark** Tel: +45 70 26 37 70 Fax: +46 83 66 310 **Finland** Tel: +358 9 251 66 123 Fax: +46 83 66 310 **France** Fisher & Paykel Healthcare SAS, 10 Av du Québec, Bât. F5, BP 512, Villebon-sur-Yvette, 91946 Courtabœuf Cedex, France Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 Email: c.s@fphcare.fr **Germany** Fisher & Paykel Healthcare GmbH, Deutschland, Österreich, Schweiz, Wiesenstrasse 49, D 73614 Schorndorf, Germany Tel: +49 7181 98599 0 Fax: +49 7181 98599 66 **Hong Kong** Tel: +852 2116 0032 Fax: +852 2116 0085 **India** Tel: +91 80 2309 6400 **Ireland** Tel: 1800 409 011 Fax: +44 1628 626 146 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea** Tel: +82 2 6205 6900 Fax: +82 2 6309 6901 **Mexico** Tel: +52 55 9130 1626 **Norway** Tel: +47 21 60 13 53 Fax: +47 22 99 60 10 **Poland** Tel: 00800 49 119 78 **Russia** Tel. and Fax: +7 495 782 21 50 **Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **Taiwan** Tel: +886 2 8751 1739 Fax: +886 2 8751 5625 **Turkey** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249, Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 Fax: +90 312 354 31 01 **UK** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 Fax: +44 1628 626 146 **USA/Canada** Tel: 1800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001

F&P 850 AIRSPIRAL™ ADULT NIV AND NHF CIRCUIT KIT

For the delivery of heated, humidified breathing gases to spontaneously breathing adult patients requiring respiratory support. This breathing set is suitable for use with Fisher & Paykel Healthcare MR850 Humidifiers in hospital and long-term care environments.

CONTRAINDICATIONS

When using the exhalation port, it should not be used on patients who:

- Are unconscious, are unable to breathe spontaneously, are uncooperative, or are unresponsive.
- Have copious secretions, are at risk of nausea or vomiting, or are at risk of aspiration of vomitus.

If symptoms of these conditions occur, discontinue treatment immediately.

TECHNICAL SPECIFICATIONS

Compatible with Fisher & Paykel Healthcare MR850 Humidifiers. Refer to humidifier user instructions.

INTERFACE CONNECTIONS

ISO 5356-1 Conical Connectors

COMPLIANCE @ 60 cmH₂O

With Accessories 1.46 ± 0.13 mL/cmH₂O

Without Accessories 1.25 ± 0.12 mL/cmH₂O

(including 0.06 mL/cmH₂O measurement uncertainty)

RESISTANCE TO FLOW @ 30 L/min

With Accessories 1.23 ± 0.13 cmH₂O

Without Accessories 1.11 ± 0.13 cmH₂O

(including 0.03 cmH₂O measurement uncertainty)

MAXIMUM OPERATING PRESSURE

8 kPa

CIRCUIT LENGTH

1.6 m

MINIMUM INTERNAL DIAMETER

17 mm

COMPRESSIBLE VOLUME

1.0 L

GAS LEAKAGE @ 60 cmH₂O

<40 mL/min

OVERALL PERFORMANCE AT 20 °C TO 26 °C AMBIENT TEMPERATURE

	NONINVASIVE MODE	INVASIVE MODE
Humidification Output	>12 mg/L	>33 mg/L
Flow Rate	10-120 L/min	10-60 L/min

FILTRATION EFFICIENCY

Viral Organism	>99.99%
Bacterial Organism	>99.99%
Mean particle size	Bacillus subtilis
NaCl	3 µm
	98.04%

	Incorrect water level, replace MR290 chamber	Rx only	Prescription only
	Correct water level in the MR290 chamber	CE 0123	CE Marking 93/42/EEC
	Consult operating instructions		Date of manufacture
	Type BF applied part		Manufacturer
	Not made with phthalates		Use-by date
	Not made with natural rubber latex		Caution/Consult instructions for use
	Single use		14 Days maximum use
	Lot number		Transportation and storage temperature limits
	Reference number		European Union authorised representative
	To Be Disposed		Do not open with blade
	Patient end		Circuit end
	Medical Device		

WARNINGS, CAUTIONS AND NOTES

⚠️ WARNINGS

- DO NOT reuse this product. Reuse may result in transmission of infectious substances, interruption to treatment, serious harm or death.
- The use of breathing circuits, chambers, accessories, or combinations which are not approved by Fisher & Paykel Healthcare may result in poor humidification system performance, ventilator malfunction and harm to the patient/user.
- Appropriate patient monitoring (e.g. oxygen saturation) must be used at all times. Failure to monitor the patient (e.g. in the event of an interruption to gas flow) may result in serious harm or death.
- DO NOT touch the heater-plate or chamber base. Surfaces may exceed 85 °C. Failure to comply may result in a skin burn.
- Exhalation port must be used with a non-vented interface on a single limb system. Failure to comply may lead to patient inhaling excess carbon dioxide resulting in hypercapnia.

Failure to comply with the following warnings may impair performance of the device or compromise safety (including potentially cause serious harm):

- DO NOT use beyond 14 days maximum duration of use.
- DO NOT soak, wash, or sterilize this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.
- DO NOT use the chamber if the water level rises above the maximum water level line.
- DO NOT use the chamber if the seals are not intact when received, or if it has been dropped.
- When mounting a humidifier adjacent to a patient, ensure that humidifier is always positioned lower than the patient.
- DO NOT operate the chamber at an angle in excess of 10°.
- DO NOT spike the water source until the blue caps have been removed. Should the primary float fail, splashing into the circuit may occur if the chamber is being operated in excess of 80 L/min.
- DO NOT use heated wire breathing circuits without gas flow. If gas flow is interrupted, turn the humidifier off.
- Ensure there is a water supply connected to the chamber and that water is present within the chamber.
- Use USP Sterile Water for Inhalation or equivalent for humidification. DO NOT add other substances to the water.
- The water source must be at least 50 cm higher than the chamber.
- DO NOT fill the chamber with water in excess of 37 °C.
- DO NOT cover circuit with materials such as blankets, towels or bed linens.
- Avoid prolonged contact of heated tubes with patient's skin.
- DO NOT stretch or milk the tubing.
- Change filter every 24 hours, or sooner if noticeable deterioration occurs, following standard hospital procedure.
- When nebulized drugs are used, resistance to flow should be monitored, and the filter replaced following standard hospital procedure.
- DO NOT block or seal the vent holes on the exhalation port.
- Ensure appropriate ventilator or flow source alarms are set before connecting breathing set to patient.
- Before connecting to patient, ensure that flow and pressure testing applicable to the ventilator has been completed.
- Visually inspect breathing sets for damage (e.g. a crushed tube or cracked connector) before use, and replace if damaged.
- Regularly monitor and drain condensate build-up in the circuit.
- Check all connections are tight before use.
- Remove any sources of ignition such as cigarettes or open flames.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit www.fphcare.com/prop65

NOTES

- For use under the supervision of trained medical personnel.
- Dispose of product according to Hospital protocol. User may be exposed to breathing tract fluids during disposal.
- The responsible organization is accountable for the compatibility of the humidifier and all of the parts and accessories used to connect to the patient and other equipment before use.
- If a serious incident has occurred while using this device please notify your local Fisher & Paykel Healthcare representative and Competent Authority.

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