













F&P Visairo

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OPERATING INSTRUCTIONS

- This mask is intended for use for a maximum of 14 days.
- This mask is for single patient use.
- For use with ventilators and devices delivering NIV therapy using a single-limb circuit system.
- For use with pressurised gases provided by an external flow source or ventilator.
- Recommended for use as part of the Fisher & Paykel Healthcare noninvasive humidification systems, F&P 850 or F&P 950™, to deliver humidified gases to the patient.
- Note:** F&P 950 System may not be available in all countries.
- Dispose of mask according to hospital protocol.
- Notice to user: Any serious incident that has occurred in relation to the device should be reported to the competent authority in your country and your local Fisher & Paykel Healthcare representative.

WARNINGS

- Appropriate patient monitoring, e.g., oxygen saturation, must be adhered to at all times. Failure to monitor the patient, e.g., in the event of flow disruption, may result in serious harm or death.
- Do not use this mask outside of a hospital/institutional environment in which the patient is adequately monitored by trained medical staff as this may result in serious harm or death.
- Verify that the therapy device, i.e., ventilator or flow source, including all alarms and safety systems are functioning correctly and that it is supplying the correct pressure(s). Failure to do so may result in excessive carbon dioxide rebreathing.
- Always use with an exhalation port. Failure to do so may result in excessive carbon dioxide rebreathing.
- Connecting components inline between the mask and exhalation port may result in excessive carbon dioxide rebreathing. Ensure appropriate patient monitoring is in place.
- This mask is fitted with an anti-asphyxiation valve and is not for use with dual-limb ventilation. Use with dual-limb ventilation may result in compromised therapy.
- Do not block the anti-asphyxiation valve openings . Doing so may result in excessive carbon dioxide rebreathing in the event of flow disruption.
- Replace the mask every 14 days or sooner if noticeable damage occurs. Failure to do so may result in compromised therapy.
- Do not reuse, soak, wash or sterilise this mask. Avoid contact with chemicals, cleaning agents or hand sanitizers. Reusing or reprocessing the mask may result in transmission of infectious substances, inhalation of particulates or compromised therapy.
- There is a risk of fire associated with oxygen equipment and therapy. Do not use near sparks or open flames.
- Smoking during oxygen therapy is dangerous and is likely to result in serious injury from fire.
- Do not use lotions or salves that are labeled as oxygen compatible to avoid the risk of fire and burns.
- Do not lubricate fittings, connections, tubing or other accessories to avoid the risk of fire and burns.

FITTING INSTRUCTIONS

- Using the SizeMe guide, size the patient by positioning the guide under the patient’s nose. The three colored nasal profiles should be facing upward. The correct mask size is determined by the colored nasal profile that the patient’s nose best fits inside: A, B or C.
- Connect the mask to the flow source. Ensure the flow source is turned on.
- If the patient has an orogastric or nasogastric tube, ensure the tube is placed under the integrated TubeFit™ zone to achieve an effective seal around the tube. Do not pass the tube through the nasal opening of the mask seal. The tape holding the orogastric/nasogastric tube on the nose should be as smooth as possible to ensure a good seal.
- Undo the lower clip(s) on the headgear. Hold the mask to the patient’s face with the small opening underneath the patient’s nose and slide the loose headgear over the patient’s head.
- Connect the headgear clip(s) back onto the mask.
- Adjust the upper headgear side straps and the crown straps as needed, so that the mask seal sits comfortably under the patient’s nose. Ensure that:
 - the crown headgear strap is positioned in the middle of the patient’s head
 - there is approximately a finger-width gap between the patient’s ear and the headgear
 - the headgear straps are not too close to the patient’s eyes.
- Adjust the lower headgear side straps so that the mask seal rests securely on the patient’s chin.
- Gently pull the mask away from the patient’s face, allowing the seal to inflate and minimize leak. **Note:** Readjust the headgear and mask as required to achieve a good seal. Remember to gently pull the mask forward to inflate the seal .

WARNINGS

- Only appropriately trained medical practitioners or care providers are to fit the mask and establish therapy. Fitment by an untrained person may result in a range of patient harms, some of which are serious.
- Remove all packaging prior to use. Packaging can occlude gas pathway resulting in compromised therapy.
- Do not use the mask if the mask is damaged. Using a damaged mask may result in compromised therapy.
- Do not use the mask if the packaging is damaged or unintentionally opened before use. Doing so may result in contamination leading to infection.
- Verify that the mask is the correct size for the patient. Incorrect sizing may result in compromised therapy, pressure injuries or eye infection.
- Do not pass an orogastric or nasogastric tube through the nasal opening of the mask seal. Doing so may result in the accidental removal of the tube leading to patient discomfort, nasal trauma or liquid aspiration.
- Do not stretch the seal over the patient’s nose. Doing so may result in compromised therapy, pressure injuries, or eye infection.
- Do not overtighten any of the headgear straps . Overtightening may result in pressure injuries.

NOTES

- Facial hair, missing teeth and dentures and facial structure irregularities may compromise mask fit or seal.
- The F&P Visairo mask is available in two headgear sizes. The standard size is designed to fit most patients. If your patient population has a smaller head size, you may wish to use the Visairo mask with Small Headgear for better strap adjustment and accessibility. "Small Headgear" is stated on the packaging and headgear. For more information, please contact your local F&P representative.

English (en)

Non-vented Hospital Under Nose Mask, Anti-asphyxiation Valve Version

INTENDED USE

The Fisher & Paykel Healthcare single patient use masks are intended for use as an accessory to ventilators to enable non-invasive positive pressure ventilation (NPPV) therapy (CPAP or bi-level) to be delivered to spontaneously breathing adult patients (>30 kg) with respiratory insufficiency or respiratory failure who have been prescribed NPPV. The masks are to be fitted and therapy maintained, by trained medical practitioners in a hospital/ institutional environment with patient monitoring in place.

CONTRAINDICATIONS

This mask should not be used on patients who:

- are having cardiac or respiratory arrest, or severe hemodynamic instability
- are unconscious, unable to breathe spontaneously, uncooperative, unresponsive or unable to remove the mask
- have an upper airway obstruction, or an inability to clear secretions (impaired cough or swallow reflexes, excessive reflux, epistaxis, hiatal hernia)
- have copious secretions, are at risk of nausea/vomiting or at high risk of aspiration of emesis
- have had head or facial surgery, trauma or burns
- have severe upper gastro-intestinal bleeding, or barotrauma (undrained pneumothorax).

If symptoms of these conditions occur, discontinue treatment immediately.

SIDE EFFECTS

After using this mask you may experience:

- redness, irritation, or pressure injury
- discomfort (eye irritation, claustrophobia, gastric insufflation, noise, dryness, nasal congestion, desynchrony)
- aspiration of fluid collected in the mask.

Please consult the appropriate medical practitioner if any of the above symptoms occur.

Operating Pressure Range	4 – 40 cmH₂O
Interface Connections	ISO 5356-1 Conical Connectors
Mask Dead Space	<200 cm ³
Anti-asphyxiation valve open to atmospheric pressure:	0.41 cmH ₂ O
Anti-asphyxiation valve closed to atmospheric pressure:	0.92 cmH ₂ O
Resistance to flow through mask	0.23 cmH ₂ O @ 50 L/min, 0.51 cmH ₂ O @ 100 L/min

- This mask was not made with natural rubber latex.
- This mask was not made with phthalates.

SYMBOL DEFINITIONS

	Non-vented	REF	Catalog Number	EC REP	European Union Authorised Representative
	Single-limb Compatible	LOT	Batch Number	CE 0123	European Conformity – TÜV SÜD
size A	Size A		Date and Country of Manufacture	UK REP	UK Authorized Representative
size B	Size B		Manufacturer		Packaging Materials
size C	Size C	MD	Medical Device	CH REP	Switzerland Authorized Representative
size A	Size A Small Headgear		Expiry Date		Warning
size B	Size B Small Headgear	Rx only	Prescription Only	14	14 Days' Maximum Use
size C	Size C Small Headgear		Not Made with Phthalates (DEHP, BBP, DBP)		Temperature Limits
	Follow Instructions for Use		Importer		Do Not Use If Package is Damaged
	Consult Instructions for Use		Keep Dry		Distributor
	Not Made with Natural Rubber Latex		Single Use		

^[1] REF RT075A/B/C F&P Visairo™ RT075 Non-vented Hospital Under Nose Mask, Anti-asphyxiation Valve Version

^[2] REF RT075A-SH / RT075B-SH / RT075C-SH F&P Visairo™ RT075 Non-vented Hospital Under Nose Mask, Anti-asphyxiation Valve Version – Small Headgear

Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Email: info@fphcare.co.nz Web: www.fphcare.com

Australia (AU) (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 **Brazil (BR)** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo - SP, Brazil Tel: +55 11 2548 7002 **China (CN)** 代理人/售后服务机构:费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486

France (FR) Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Email: c.s@fphcare.fr **Germany (DE)** Fisher & Paykel

Healthcare GmbH, Wiesenstrasse 49, 73614 Schorndorf, Germany Tel: +49 7181 98599 0 **Hong Kong (HK)** Tel: +852 2116 0032 **India (IN)** Tel: +91 80 2309 6400 **Japan (JP)** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115

Korea (KR) Tel: +82 2 6205 6900 **Mexico (MX)** Tel: +52 55 9130 1626 **Poland (PL)** Fisher & Paykel Healthcare Poland Sp. z o.o., Pl. Andersa 7, 61-894 Poznań, Poland Tel: +48 664 846 464 **Russia (RU)** Tel: +7 495 782 21 50

Switzerland (CH) **CH REP** Fisher & Paykel Healthcare GmbH, Säntisstrasse 2, 9501 Wil / SG, Switzerland Tel: 0800 83 47 63 **Taiwan (TW)** Tel: +886 2 8751 1739 **Turkey (TR)** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 **UK (GB)** **UK REP** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 **USA (US)/Canada (CA)** Tel: 1800 446 3908 or +1 949 453 4000

Austria (AT) Tel: 0800 29 31 23 **Benelux (BE NL LU)** Tel: +31 40 216 3555 **Denmark (DK)** Tel: +45 70 26 37 70 **Finland (FI)** Tel: +358 9 251 66 123 **Ireland (IE)** Tel: 1800 409 011 **Italy (IT)** Tel: +39 06 7839 2939 **Norway (NO)** Tel: +47 21 60 13 53 **Spain (ES)** Tel: +34 902 013 346 **Sweden (SE)** Tel: +46 8 564 76 680

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REF 619435 REV F 2023-12 © 2023 Fisher & Paykel Healthcare Limited
F&P, Tubefit, Visairo and F&P 950 are trademarks of Fisher and Paykel Healthcare Limited.