











**Non-vented Hospital Under Nose Mask**  
Anti-asphyxiation Valve Version

**When using NG/OG Tubes**

**CE 0123 Rx only**

**[14]**

F&P, Tubefit, Visairo and F&P 950 are trademarks of Fisher and Paykel Healthcare Limited.

## English (en)

Non-vented Hospital Under Nose Mask, Anti-asphyxiation Valve Version

### INTENDED USE

The Fisher & Paykel Healthcare single-patient-use masks are intended for use as an accessory to ventilators to enable noninvasive positive pressure ventilation (NPPV) therapy (CPAP or bi-level) to be delivered to spontaneously breathing adult patients (>66 lb, >30 kg) with respiratory insufficiency or respiratory failure who have been prescribed NPPV. The masks are to be fitted and therapy maintained by trained medical practitioners in a hospital/institutional environment with patient monitoring in place.

### TECHNICAL SPECIFICATIONS

OPERATING PRESSURE RANGE	4 - 40 cmH <sub>2</sub> O
INTERFACE CONNECTIONS	ISO 5356-1 Conical Connectors
MASK DEAD SPACE	<200 cm <sup>3</sup>
ANTI-ASPHYXIA VALVE FLAP OPEN TO ATMOSPHERIC PRESSURE	0.41 cmH <sub>2</sub> O
ANTI-ASPHYXIA VALVE FLAP CLOSED TO ATMOSPHERIC PRESSURE	0.92 cmH <sub>2</sub> O
RESISTANCE TO FLOW THROUGH MASK	0.23 cmH <sub>2</sub> O @ 50 L/min, 0.51 cmH <sub>2</sub> O @ 100 L/min

- [14]** This product is intended for use for a maximum of 14 days.
- For use with pressurized gases provided by an external flow source or ventilator.
- For use with ventilators/devices delivering 'NIV' therapy using a vented single limb circuit system.
- Can be used with a heated pass over humidifier to deliver humidified gases to the patient.
- When humidifying gases, use as part of Fisher & Paykel Healthcare noninvasive humidification systems F&P MR850 or F&P 950™.

**Note:** F&P 950 System may not be available in all countries.  
F&P 950 System is currently not cleared for sale in the USA.

- This product was not made with natural rubber latex.
- Dispose of mask according to hospital protocol.
- If a serious incident has occurred while using this device, please contact your local Fisher & Paykel Healthcare representative and Competent Authority.

### LABEL - MASK SYMBOLS

	Non-vented	Seal Size A	Seal Size B	Seal Size C
	Standard Headgear			
	Small Headgear			

### SYMBOL DEFINITIONS

	Rx only	Prescription only	<b>[14]</b>	14 Days maximum use
	<b>CE 0123</b>	CE Marking 93/42/EEC		European Union authorised representative
	Reference number		Do not use if package is damaged	

### BEFORE USE

- Ensure you have read and understood the instructions completely.
- Remove all packaging prior to use.
- Do not use if package is damaged.
- Inspect the mask for damage. Discard the mask if any parts are broken or if the seal is torn.
- Ensure that the correct mask is being used for the therapy device. Refer to setup guide on the front most page
- Verify that the mask is the correct size for the patient.
- Verify operation of the anti-asphyxiation valve before use.
- Verify that the therapy device (i.e. ventilator or flow source), including all alarms and safety systems, are functioning correctly and that it is supplying the correct pressure(s).
- Clean the patient's face and then inspect for signs of redness, irritation or discomfort. Do not use if these such signs are present.
- Ensure adequate patient monitoring is in place.

### WARNINGS, CAUTIONS AND CONTRAINDICATIONS

#### WARNINGS

- For CPAP/Bi-Level use with a single limb circuit **an exhalation port must be fitted**.
- Do not block or try to seal the exhalation port.**
- Verify that the therapy device, including alarms and safety systems, are functioning correctly prior to use.
- This mask is fitted with an anti-asphyxiation valve and is not for use with dual limb ventilation systems. With no system flow the external openings in the valve should allow room air to entrain into the mask. With system flow the flap should close the external openings and allow system air to flow into the mask.
- Do not block anti-asphyxiation valve vents.**

- Replace the mask if the anti-asphyxiation valve does not operate or becomes fouled with secretions.
- The mask must only be worn when the therapy is being delivered.
- The mask must be fitted and therapy established by an appropriately trained medical practitioner or care provider.
- This mask may only be used in a hospital or clinical setting where the patient is adequately monitored by trained medical staff. Failure to monitor the patient may result in loss of therapy, serious injury or death.
- This mask is for single patient use. **Do not reuse.** Reuse may result in transmission of infectious substances, interruption to treatment, serious harm or death. Do not soak, wash, sterilize, or reuse this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.
- At low CPAP or EPAP pressures the flow through the exhalation port may be too low to clear all exhaled gas from the mask. Some rebreathing may occur.
- Facial hair, missing teeth/dentures and facial structure irregularities may compromise mask seal.

#### CAUTIONS

- Do not overtighten any of the headgear straps. Overtightening may cause patient discomfort and/or leaks.
- If the patient experiences skin redness, irritation or discomfort, discontinue use and contact a physician.

#### CONTRAINDICATIONS

This mask should not be used on patients who:

- Have cardiac or respiratory arrest, or severe hemodynamic instability
- Are unconscious, unable to breathe spontaneously, uncooperative, unresponsive or unable to remove the mask
- Have an upper airway obstruction, or an inability to clear secretions (impaired cough or swallow reflexes, excessive reflux, epistaxis, or hiatal hernia)
- Have copious secretions, at risk of nausea/vomiting, or at high risk of aspiration of emesis
- Have had head or facial surgery, trauma or burns
- Have severe upper gastro-intestinal bleeding, or barotrauma (un-drained pneumothorax).

If symptoms of these conditions occur, discontinue treatment immediately.

#### FITTING INSTRUCTIONS

- Using the SizeMe guide, size the patient by positioning the guide under the patient's nose. The three colored nasal profiles should be facing upward. The correct mask size is determined by the colored nasal profile that the patient's nose best fits inside: A, B or C.
- Connect the mask to the flow source. Ensure the flow source is turned on.
- If the patient has an orogastric or nasogastric tube, ensure the tube is placed under the integrated TubeFit™ zone to achieve an effective seal around the tube. Do not pass the tube through the nasal opening of the mask seal. The tape holding the orogastric/nasogastric tube on the nose should be as smooth as possible to ensure a good seal.
- Undo the lower clip(s) on the headgear. Hold the mask to the patient's face with the small opening underneath the patient's nose and slide the loose headgear over the patient's head.
- Connect the headgear clip(s) back onto the mask.
- Adjust the upper headgear side straps and the crown straps as needed, so that the mask seal sits comfortably under the patient's nose. Ensure that:
  - the crown headgear strap is positioned in the middle of the patient's head
  - there is approximately a finger-width gap between the patient's ear and the headgear
  - the headgear straps are not too close to the patient's eyes.
- Adjust the lower headgear side straps so that the mask seal rests securely on the patient's chin.
- Gently pull the mask away from the patient's face, allowing the seal to inflate and minimize leak.

**Note:** Readjust headgear/mask as required to achieve an effective seal. Do not overtighten the mask. When fitted properly, you should be able to slide a finger underneath the upper and lower side straps. The F&P Visairo mask is available in two headgear sizes. The standard size is designed to fit most patients. If your patient population has a smaller head size, you may wish to use the Visairo mask with Small Headgear for better strap adjustment and accessibility. 'Small Headgear' is stated on the packaging and headgear. For more information, please contact your local F&P representative.

For patent information, refer to [www.fphcare.com/ip](http://www.fphcare.com/ip)

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