Next generation care


Instructions for use

[REF] WJR210, WJR212, WJR214

For use with

F&P Optiflow™ Junior 2

F&P Optiflow Junior 2+

www.fphcare.com



Fisher & Paykel
HEALTHCARE

I



F&P WigglewiNG™ fitting



A. RIP into two

R**L**

B. STICK to face




C. FOLD to secure NG tube

F&P WigglewiNG removal

F&P Optiflow Junior 2 ¹ Cannula Size	F&P WigglewiNG
OJR410 XS 	WJR210
OJR412 S 	WJR212
OJR414 M 	
OJR416 L 	
OJR418 XL 	

F&P Optiflow Junior 2+ ¹ Cannula Size	F&P WigglewiNG
OJR520 XXL 	WJR214

¹ Refer to the F&P Optiflow Junior 2 and F&P Optiflow Junior 2+ nasal cannula user instructions for appropriate sizing.

English	4
العربية (Arabic)	7
български (Bulgarian)	9
Česky (Czech)	12
Dansk (Danish)	15
Deutsch (German)	18
Ελληνικά (Greek)	21
Español (Spanish)	24
Eesti keel (Estonian)	27
Suomi (Finnish)	30
Français (French)	33
עברית (Hebrew)	36
Hrvatski (Croatian)	38
Magyar (Hungarian)	41
Bahasa Indonesia (Indonesian)	44
Italiano (Italian)	47
日本語 (Japanese)	50
한국어 (Korean)	53
Lietuvių (Lithuanian)	56
Latviešu (Latvian)	59
Nederlands (Dutch)	62
Norsk (Norwegian)	65
Polski (Polish)	68
Português (Portuguese)	71
Português (Brasil) (Brazilian Portuguese)	74
Română (Romanian)	77
Русский (Russian)	80
Slovenčina (Slovak)	83
Slovenščina (Slovenian)	86
Српски (Serbian)	89
Svenska (Swedish)	92
ไทย (Thai)	95
Türkçe (Turkish)	98
Tiếng Việt (Vietnamese)	101
简体中文版 (Simplified Chinese)	103
繁體中文版 (Traditional Chinese)	105

F&P, Optiflow, Wigglepads, WigglewiNG, and the Seahorse, Crab, Starfish, Octopus, Turtle, Dolphin and Stingray images shown in this user instruction are trade marks of Fisher & Paykel Healthcare Limited.

For patent information, see www.fphcare.com/ip

California residents please be advised of the following, pursuant to Proposition 65:

WARNING

This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm.

For more information, please visit: www.fphcare.com/prop65

Intended use

The F&P WigglewiNG is intended to secure nasogastric tubing to the face in conjunction with the F&P Optiflow Junior 2 and F&P Optiflow Junior 2+ nasal cannula ranges.

The F&P WigglewiNG is a single use accessory intended for use on neonates, infants and children in a hospital environment.

Contraindications/side effects

For contraindications and side effects refer to Interface User Instruction.

GENERAL WARNINGS

- This product is only designed and verified for use with equipment, accessories and spare parts approved by F&P. Unauthorized equipment, accessories or spare parts which are used with this product may impair performance of this product or compromise safety (including potentially causing serious patient injury).
- Appropriate patient monitoring (e.g. oxygen saturation) must be used at all times. Failure to monitor the patient (e.g. in the event of an interruption to gas flow) may result in serious injury or death.
- If using supplemental oxygen, keep ignition sources away from the patient.
- Pre-existing craniofacial abnormalities, deformities, malformations or trauma may be exacerbated by NHF interfaces and/or the retention mechanism and may not permit therapy to be delivered as intended, leading to further injury or death.
- The use of this product is not without risk, even if used as intended. Following all instructions and warnings provided, risks of hypoxic injury and skin damage remain. These risks may result in serious injury or death.
- Application of NHF is known to generate positive airway pressure which may exacerbate preexisting, untreated serious air leak syndrome which may lead to further serious injury or death.
- Any serious incident that has occurred in relation to this device should be reported to your Fisher & Paykel Healthcare representative and the local competent authority.

GENERAL CAUTIONS

- This product is intended to be used for a maximum of 7 days. Using this product beyond 7 days may impair performance of this product or compromise safety (including potentially causing serious patient injury).
- Regularly monitor the patient to ensure skin integrity and that the skin underneath the F&P WigglewiNG remains dry. A barrier film may be used between the NG tubing and the patient's upper lip to prevent irritation.
- DO NOT soak, sterilize or reuse this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.
- Reuse may result in the transmission of infectious substances, interruption to treatment, serious injury or death.
- DO NOT place the F&P WigglewiNG on the patient's eyes, ears or injured skin.
- DO NOT stretch the cannula on application; this may cause increased pressure to the patient's skin. If necessary, the cannula may be repositioned.
- Ensure the cannula is placed directly on the F&P WigglewiNG. Direct skin contact caused by cannula misalignment may result in skin breakdown.
- DO NOT use if the product or its packaging has been tampered with.
- Failure to apply and use this product within the directions, transport, storage and operating conditions specified in the labeling and user instructions may impair performance of this product or compromise safety (including potentially causing serious patient injury).

I F&P WigglewiNG fitting

For the correct F&P WigglewiNG positioning, refer to the placement of the F&P Wigglepads™ 2 in the F&P Optiflow Junior 2 and F&P Optiflow Junior 2+ nasal cannula user instructions.

A. RIP into two

R Patient's right nostril L Patient's left nostril

- Prepare the patient's skin according to hospital protocol.
- 1. Determine which nostril you would like to place the NG tube in.
 - Hold the F&P WigglewiNG over the patient's cheek on the same side as this nostril.
- 2. Face the white tabs (marked **2**) towards you and the hook pads away from you.

- Place the straighter edge towards the patient's nose.
- 3.** Rip off the entire upper third of the F&P WigglewiNG (distal to the patient's chin) by holding the middle section and ripping off the upper section.
- 4.** This will leave you with the lower 2/3 of the F&P WigglewiNG.

B. Stick to face

- 5.** Peel off the blue tab (marked **①**).
- 6.** Stick the F&P WigglewiNG to the patient's cheek
- 7.** Peel off the white tabs (marked **②** and F&P WigglewiNG).
 - Insert the NG tube according to hospital protocol.
- 8.** Place the NG tube along the perforations.

C. Fold to secure the NG tube

- 9.** Fold the F&P WigglewiNG over to secure the NG tube.
- 10.** Apply the section removed in Step 3 as you would a Wigglepad on the opposite cheek.
 - Remove the existing Wigglepad from the nasal cannula.
- 11.** Secure the cannula in place so the cannula bridge rests close to the nose without touching the septum. Do not stretch the cannula during application.

II F&P WigglewiNG removal

- Place fingertip on the outside edge of the F&P WigglewiNG and gently peel nasal cannula away. Starting from the outside, peel towards the nose.
- 1.** Lift the edge of the F&P WigglewiNG. Use a damp cloth to wipe the patient's skin and the underside of the F&P WigglewiNG while gently peeling away from the patient's face.
- 2.** Holding the NG tube securely between the patient's nose and the F&P WigglewiNG, detach the NG tube from the F&P WigglewiNG by tearing along the perforations.

Disposal

There are no specific disposal requirements for this device. Disposal protocols appropriate to single use devices that may be contaminated with potentially infectious substances should be followed.

Operating conditions

Ambient temperature range: 18 to 26 °C.

Note: This product may be used in an incubator.

Approved F&P respiratory interfaces

F&P Optiflow Junior 2 nasal cannula
F&P Optiflow Junior 2+ nasal cannula

Refer to the F&P Optiflow Junior 2 and F&P Optiflow Junior 2+ nasal cannula user instructions for appropriate sizing.

Symbol definitions

	Do not reuse
	Not made with phthalates (DEHP, DBP, BBP)
LOT	Batch code
	Consult instructions for use. fphcare.com/oj2-ifu
EC REP	Authorized representative in the European Community
	Use by date YYYY-MM-DD
	Date and Country of Manufacture NZ: New Zealand YYYY-MM-DD
	Maximum 7 days use
	This product is not made with natural rubber latex
	Storage temperature range -10 °C (14 °F) / +50 °C (122 °F)
REF	Catalogue number
Rx only	Prescription only
	European Conformity
MD	Medical device
	Manufacturer
	Importer
	Distributor
CH REP	Switzerland authorized representative
UK REP	UK responsible person

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REF 607831 REV D 2024-03 © 2024 Fisher & Paykel Healthcare Ltd