Am J Respir Crit Care Med. 2018 Dec 21. doi: 10.1164/rccm.201803-0589OC. [Epub ahead of print]

An Index Combining Respiratory Rate and Oxygenation to Predict Outcome of Nasal High Flow Therapy.

Roca O^{1,2}, Caralt B³, Messika J⁴, Samper M⁵, Sztrymf B⁶, Hernández G⁷, García-de-Acilu M³, Frat JP⁸ Masclans JR^{10,5}, Ricard JD¹¹.

Roca et al. 2018

American Journal of Respiratory and Critical Care Medicine

An Index Combining Respiratory Rate and Oxygenation to Predict Outcome of Nasal High Flow Therapy.



Introducing the ROX Index

Roca et al. 2018

An index combining respiratory rate and oxygenation to predict outcome of Nasal High Flow (NHF) therapy.

$$SpO_2/FiO_2$$
-RR



Key points

Roca et al. 2018

- The ROX index can help to predict NHF therapy outcomes in patients with acute hypoxic respiratory failure (AHRF) due to pneumonia.
- The index can be calculated easily and repeatedly at the patient's bedside.
- The index is dynamic.



Background

Roca et al. 2018

In a pilot study, Roca (2016) previously demonstrated that the ROX index:

- outperformed the diagnostic accuracy of its individual variables
- participants with an index ≥ 4.88 and > 12 hours of ongoing NHF therapy were less likely to be intubated, even after adjusting for potential covariates.

It was anticipated that the ratio of the ROX index may predict failure in adult AHRF patients on NHF therapy.

Further independent validation of this index was required.

Aim

- To validate the ROX index value of ≥ 4.88 as a determinant of NHF success at 2, 6, or 12 hours.
- To identify a series of ROX score margins for NHF failure at four time points.
- To perform external validations using the FLORALI study cohort (Frat, 2015).



Background

Roca et al. 2018

Setting

Five French and Spanish ICUs 2016 - 2017

Population

N = 191, adults with AHRF* followed until death or discharge

*AHRF defined according to the IDSA / ATS 2007 guidelines.



A multi-center prospective non-randomized observational cohort study.

INCLUSION CRITERIA	EXCLUSION CRITERIA
ICU patients with AHRF due to pneumonia, on NHF therapy	< 18 years old
	Indication for immediate intubation
	Limitation of therapy (DNR)
	Elective intubation for diagnostic or therapeutic procedures

Primary outcome

NHF failure* requiring invasive ventilation within study completion time.

*Failure defined by set of predetermined criteria for intubation.



Roca et al. 2018

Pre-determined criteria for intubation

- GCS < 12
- Cardiac arrest
- Arrhythmias
- Severe hemodynamic instability
- Two or more of the following
 - PaO₂ < 60 mmHg or SpO₂ < 90% on NHF ≥ 30 L/min FiO₂ 100%
 - Respiratory acidosis PaCO₂ > 50 mmHg
 or Pv CO₂ > 55 mmHg, pH < 7.25
 - RR > 30 BPM
 - Inability to clear secretions



Treatment regimen

- Participants were on FM O₂ ≥ 10 L/min with SpO₂ < 92% and RR ≥ 25 breaths/min. NHF was initiated at 30 L/min, FiO₂ 100%
- FiO₂ titrated with a target $SpO_2 > 92\%$
- Flow rate adjusted according to the maximum tolerated rate which was achieved within the first 10 mins of NHF treatment - for all participants.



Outcome

- Accuracy of the ROX index for predicting NHF outcomes (need for intubation or not)
 - Determined through Cox proportional hazards modelling
 - ROX index was calculated as the ratio of SpO_2/FiO_2 to RR
 - Variables with a positive association with NHF success were in the numerator field, and the variables with an inverse relation to NHF success were in the denominator field.

NHF failure definition

Requiring intubation and invasive mechanical ventilation.

NHF success definition

Not requiring intubation or invasive mechanical ventilation.



Baseline characteristics

CHARACTERISTIC				
Age (years)	63 (52 – 72)			
Pneumonia severity index	116 (80 – 152)			
APACHE II	16 (12 – 21)			
SOFA ICU admission	5 (3 – 7)			
Presence of shock (%)	20.4			

At baseline these study participants were sicker than those enrolled in the pilot (Roca, 2016).



Roca et al. 2018

ROX index threshold ≥ 4.88 (determinant for NHF success) was validated with three time points:

- 2 h (HR 0.434 [95% CI 0.264 0.715]; P= 0.001)
- 6 h (HR 0.304 [95% CI 0.182 0.509]; P<0.001)
- 12 h (HR 0.291 [95% CI 0.161 0.524]; P<0.001)

Three time points were validated for this threshold as opposed to one as previously described ≥ 12hrs (Roca, 2016).



Respiratory variables during NHF therapy comparing success vs. failure (median)

VARIABLE*	TIME	SUCCESS (N=123)	FAILURE (N=68)	P-VALUE#
SpO ₂ /FiO ₂	Prior to NHF	180	106	0.005
3pO ₂ /11O ₂	2 h	155	109	0.003
	6 h	160	115	0.001
	12 h	165	113	0.001
	18 h	176	118	0.002
	24 h	194	120	< 0.001
Respiratory	Prior to NHF	28	32	0.778
rate (bpm)	2 h	25	28	0.023
\	6 h	24	26	0.003
	12 h	23	26	< 0.001
	18 h	22	25	0.001
	24 h	21	24	0.004
ROX index	Prior to NHF	5.81	4.06	0.169
	2 h	5.71	4.43	0.001
	6 h	6.55	4.86	< 0.001
	12 h	7.53	4.78	< 0.001
	18 h	8.6	5.10	< 0.001
	24 h	8.68	5.05	< 0.001

In those who had success, a higher ROX and SpO_2/FiO_2 versus lower RR were observed after and throughout the study.



Roca et al. 2018

Results

ROX score margin for failure over time

HOURS OF NHF USE	ROX SCORE FOR FAILURE	AUROC	POSITIVE PREDICTIVE VALUE %
2	< 2.85	0.679	98
6	< 3.47	0.703	98 – 99
12	< 3.85	0.759	99
> 12	< 4.88	NA	80

The prediction accuracy of the ROX index increased over time.



Roca et al. 2018

Results

Hospital and ICU mortality and length of stay (LOS)

	SUCCESS N=123	FAILURE N=68	P-VALUE
ICU mortality	6 (4.9%)	29 (43.3%)	<0.001
Hospital mortality	17 (14.3%)	33 (51.6%)	<0.001
ICU LOS	6 (4 – 9)	20 (10 – 28)	<0.001
Hospital LOS	17 (11 – 29)	35 (23 – 65)	0.002

Those who succeeded on NHF therapy had lower rates of mortality and length of stay.



- 35.6 % (N=68) participants failed, requiring intubation.
- Median duration success: 96h (48 144 hours) (P=0.001).
- Median duration failure: 24h (12 60 hours) (P<0.001).
- The prediction accuracy of the ROX index increased overtime.
- Participants with NHF failure presented a lower increase in the values of the ROX index over 12 hours.
- In those who had success, a higher ROX, SpO₂/FiO₂ and lower RR were observed after and throughout the study.
- ROX index < 2.85 at 2h, < 3.47 at 6h, and <3.85 at 12h post NHF initiation were predictors of failure.
- ROX index ≥ 4.88 for success was validated at three time points, in addition to the one first described (Roca, 2016).



Conclusions

- The ROX index can help to predict the outcome of NHF therapy in patients with AHRF due to pneumonia.
- A ROX index of ≥ 4.88, measured at 2, 6, or 12 hours, is a determinant of NHF success.
- ROX index values were identified and validated for three different time points which can predict NHF failure.
- Among measured components, SpO₂/FiO₂ had a greater predictive weight than RR.
- The index is dynamic and can be measured with ease and repetition at the patient's bedside.



Notes

Roca et al. 2018

- Adult pneumonia patients AHRF only.
- AIRVO[™]2 or MR850 used.



Additional resources

Roca et al. 2018

- PubMed link
 - www.ncbi.nlm.nih.gov/pubmed/30576221
- Full paper not freely available
- Clinical summary
- Clinical trial register
 - www.clinicaltrials.gov, ID: NCT02845128
- ROC explained
 - www.youtube.com/watch?v=21lgj5Pr6u4



Glossary

Roca et al. 2018

AUROC:

Area under receiver operating curve.

ROC:

A receiver operating characteristic curve or ROC curve, is a plot illustrating the diagnosticability of a binary classifier system as its discrimination threshold is varied.

HR:

Hazard ratio is the ratio of (chance of an event occurring in the treatment arm)/(chance of an event occurring in the control arm). The HR has also been defined as the ratio of (risk of outcome in one group)/(risk of outcome in another group), occurring at a given interval of time.

COX PROPORTIONAL HAZARDS MODELING:

Cox regression (or proportional hazards regression) is the method for investigating the effect of several variables upon the time a specified event takes to happen.

RR:

Respiratory rate.

