
SleepStyle[™]200

CPAP SERIES

*Ambient
Tracking[™]
Plus*

Operating Manual


Fisher & Paykel
HEALTHCARE


English	A
Français	B
Español	C


PACKAGE CONTENTS

A

B

SleepStyle™ 200 Series HUMIDIFIED CPAP SYSTEMS

The SleepStyle™ 200 Series are a range of CPAP systems designed for use in the home for the treatment of OSA.

For further assistance, please contact your local Fisher & Paykel Healthcare office – see reverse for addresses. Please keep this manual for future reference.

TABLE OF CONTENTS

1. Symbol Definitions	A - 2
2. Warnings, Cautions and Contraindications	A - 2
3. Oxygen Use Warnings	A - 3
4. Intended Use	A - 3
5. Description of the Device	A - 3
5.1 Important Parts of Your Device - Standard CPAP Range	A - 3
5.2 Important Parts of Your Device - Integrated CPAP Range	A - 3
5.3 Removing Connector	A - 3
5.4 Accessories	A - 3
6. Cleaning and Maintenance	A - 4
7. Setup Instructions	A - 4
8. Controls and Display - Integrated CPAP Range	A - 5
8.1 Display Descriptions and Functions	A - 5
8.2 Patient Menu Summary - Integrated CPAP Range	A - 6
8.3 Patient Menu - Additional Settings Summary - Integrated CPAP Range	A - 6
9. Controls and Display - Standard CPAP Range	A - 7
9.1 Display Descriptions and Functions	A - 7
9.2 Patient Menu Summary - Standard CPAP Range	A - 8
9.3 Patient Menu - Additional Settings Summary - Standard CPAP Range	A - 8
10. Frequently Asked Questions	A - 9
11. Operating Conditions	A - 9
12. Storage and Transport Conditions	A - 9
13. Troubleshooting	A - 9
14. Device and Consumables Disposal Instructions	A - 9
15. Product Specifications	A - 10
16. Appendix	A - 10

PLEASE READ ALL INSTRUCTIONS BEFORE INITIAL USE

Caution: USA Federal Law restricts this device for sale by or on the order of a physician.

1. SYMBOL DEFINITIONS

Type BF
Applied Part

ATTENTION
Consult
accompanying
documents

IPX1
Drip Proof

Alternating
Current

Class II
Double Insulated

Standby or
On/Off

CE 0123

93/42/EEC
Class IIb

2. WARNINGS, CAUTIONS AND CONTRAINDICATIONS

NOTES

- This manual refers to the SleepStyle™ 200 CPAP Series unit as "the device".
- If required to use CPAP by a referring physician, you should use your device every time you sleep. Should your device stop working for any reason, contact your healthcare provider immediately.
- The user of this system shall have sole responsibility and liability for any injury to persons or damage to property resulting from operation of the device which is not in accordance with the operating instructions supplied.
- The device should only be used with ISO 17510-2 compliant masks, connectors and delivery tubes recommended by Fisher & Paykel Healthcare, or your healthcare provider.
- Refer all repairs and maintenance to Fisher & Paykel Healthcare.

WARNINGS

To avoid electric shock from your device, do not:

- Operate the device if the power cord or plug is damaged.
- Operate the device if it has been dropped in water.
- Plug the device into the power socket if it is wet.
- Clean the device while connected to the power socket.
- Store or use the device where it can tilt, fall or be pulled into water. If water has entered the unit enclosure, disconnect the power cord and discontinue use. Seek advice from Fisher & Paykel Healthcare.

To avoid choking or inhalation of a foreign body:

- Never place any object into any opening of the tube.
- Ensure the air filter is fitted during device use.

To ensure optimal therapy, do not:

- Operate the device if dropped or damaged.
- Operate the device if not working properly.
- Adjust the pressure. Pressure must only be adjusted by a healthcare professional.
- Operate the device if the breathing tube has been damaged with holes, tears or kinks.
- Block the exhaust flow on the mask.
- Use the mask if the unit is turned off or not operating properly.

To avoid burns, do not:

- * Fill the chamber with boiling water.
- * Touch the exposed heater plate or chamber base.

To avoid the risk of fire while using oxygen, do not:

- Turn oxygen flow on when the device is not operating; this can lead to accumulation of oxygen within the device.
- Locate the device in a position where ventilation around the device is restricted.
- Use oxygen while smoking or in the presence of an open flame.
- Use any materials which will burn in air or ignite easily at high oxygen concentration.
- Keep any source of ignition near the product. To avoid ignition, it is preferable to keep all sources of ignition out of the room where supplemental oxygen is being used.
- Operate the device without the cover over the serial port adaptor.
- Keep oxygen regulators, cylinder valves, tubing, connections and all other oxygen equipment near oil, grease or greasy substances. Spontaneous and violent ignition may occur if these substances come into contact with oxygen under pressure.

Other:

- * Place the device on a level surface lower than head height to prevent water entering the tubing. If water does enter the tubing, drain excess condensate. Water in the tubing may result in aspiration.
- The device is intended to be used with CPAP masks and connectors that have exhaust flow holes to allow continuous flow of air out of the mask. When the device is turned on and operating properly new air flushes exhaled air out of the mask through the exhaust flow holes. At low CPAP pressures and in the event of power failure or machine malfunction, remove the mask immediately, as flow through the mask may be insufficient to clear all exhaled gas and CO₂ re-breathing may occur which can be hazardous.
- Failure to select the correct altitude level (for any given location) will have an adverse effect on delivered pressure.
- This device is not intended for life support.
- Avoid unnecessary removal of the power cord from the rear of the device. If removal is necessary then hold the connector and remove. Avoid pulling on the power cord (refer to Section 5.3).
- Humidity performance of the device can be compromised when used outside the specified ambient temperature range or humidity range.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65.

CAUTIONS

To prevent water damage to your device:

- * Remove the humidification chamber from the device before filling.
- * Empty water from the chamber before transporting the device. If the device is required to be handled with water in the chamber, avoid tilting the device to prevent water entering its enclosure.

Other:

- To prevent damage to your PC, only operate the device if it is connected to a PC via an isolated serial port adaptor (900HC236).
- To prevent airway irritation, do not use the device when room temperature exceeds 95 °F (35 °C) as this may lead to an elevated airflow temperature.

CONTRAINDICATIONS

- Research indicates the following pre-existing conditions may contraindicate the use of positive pressure for some patients: pneumothorax, bullous lung disease, pneumocephalus, cerebrospinal fluid leak, recent cranial surgery or trauma, abnormalities of the cribriform plate, pathologically low blood pressure or in patients whose upper airways are bypassed.

PRECAUTIONS

- The safety and effectiveness of positive pressure has not been established in patients with respiratory failure or COPD.

ADVERSE EFFECTS

- Nosebleeds, ear and sinus discomfort may occur from the use of positive pressure therapy.

Please contact your physician if you have any questions concerning your therapy.

* Only applicable to integrated CPAP devices

3. OXYGEN USE

If oxygen is required, it is recommended that supplemental oxygen be administered at the mask. Please see instructions specific to your mask type.

NOTE:

- At a fixed flow rate of supplemental oxygen, the inhaled oxygen concentration will vary, depending on the pressure settings, patient's breathing pattern, mask selection and leak rate.


Before using oxygen with the device please see oxygen warnings in Section 2.

4. INTENDED USE


The device is used to assist with patient breathing while sleeping, for the purpose of treating Obstructive Sleep Apnea (OSA). This is done by the delivery of Continuous Positive Airway Pressure (CPAP) in order to prevent airway obstruction. The addition of heated respiratory humidification to the device relieves the drying and irritating effects on the patient airways which usually arises from use of a CPAP system. The CPAP device is for use on adult patients at home or in the sleep lab. Please refer to the SleepStyle™ 200 CPAP series product range listed in the Appendix section of this manual.

5. DESCRIPTION OF THE DEVICE


5.1 Important Parts of Your Device - Standard Range

5.2 Important Parts of Your Device - Integrated Range

5.3 Removing Connector

5.4 Accessories

HC014 Conversion Kit – Heated Humidification

HC385S Standard Humidification Chamber

HC355* Extended Life Humidification Chamber

***Not available in all countries**

900HC010 Connector

900HC221 Breathing Tube to fit 22 mm (.86") conical connector

900HC240 Filter

900HC225 Water Funnel (10-Pack)

6. CLEANING AND MAINTENANCE

Please Note: The device should be cleaned as required.

1. Turn off the device and unplug from the power socket. Wait for the water to cool, then remove chamber.
2. Wipe the exterior of the device with a clean, damp (not wet) cloth and mild dishwashing detergent. Do not use harsh abrasives or solvents, as these may damage the device.

3. DAILY


Clean chamber and tube.

- Remove the breathing tube from the chamber.
- Rinse the tube with soapy water. Hang up hose with hose ends pointing to floor to dry.
- Remove the chamber by pushing down the finger guard and pulling out the chamber.
- Pour out and discard the remaining water.
- For standard (HC385S) chambers, clean and wash with soapy water, then rinse and dry. Extended life chambers (HC355) can be cleaned in a domestic dishwashing machine.

4. WEEKLY

Thoroughly clean the chamber.

- Soak the inside of the chamber for 10 minutes in a solution of one part white vinegar to two parts water.
5. Replace the air filter when it becomes significantly discoloured, at least once every 3 months or after 1000 hours of machine running time.
 - Remove the filter holder from the back of the device and take out the filter.
 - Replace the old filter with a new filter: ensure the vertical black line is facing towards the device.

7. SETUP INSTRUCTIONS

HOW TO SET UP YOUR SYSTEM

1. Remove your device from its packaging.
2. Place the device on a low shelf or on the floor beside your bed, so that the device is positioned below head height.
3. Chamber Setup (only applicable to integrated CPAP devices):
 - a. Remove one water chamber from the packaging.
 - b. Remove the blue caps and discard them (Fig. 1).
 - c. Fill the chamber up to the water filling line with distilled water only (Fig. 2).

Never fill the chamber while it is attached to the device.
When moving your device, ensure that the water chamber is empty. Device failure due to water damage is not covered by warranty.

 - d. To attach the water chamber to the device, press down the finger guard and slide the chamber on (Fig. 3 and 4).
 - e. When the chamber is fitted correctly, the finger guard will click into place.
4. One end of the tube has a white plastic connector. Push this onto the outlet on top of the chamber (Fig. 5) or CPAP outlet (Fig. 5a).
5. Connect the other end of the tube to the mask.
6. Plug the device power cord into your household power socket and into the rear of the device.
7. Switch on the device by pressing the on/off button (See Fig. 6). The current humidity setting is displayed.
8. For changing the setting refer to Section 8 "Controls and Display".


Fig. 1


Fig. 2


Fig. 3


Fig. 4


Fig. 5


Fig. 5a



Fig. 6

Your device is now ready for use.

8. CONTROLS AND DISPLAY - INTEGRATED CPAP RANGE

8.1 Display Description and Functions

KEY	FUNCTION & DESCRIPTION	OPERATION
	ON/OFF Switches the device ON & OFF	TO SWITCH "ON" - Press briefly - The device will default to display the humidity setting TO SWITCH "OFF" - Remove your mask - Press briefly
	RAMP Reduces pressure to a lower level to help you fall asleep. The pressure will gradually return to full operating pressure over 20 minutes	TO START THE RAMP - Press the button NB: to reset ramp during a cycle the device must be switched off first TO SWITCH "OFF" - Press the button
	INCREASE/DECREASE Adjusts settings up or down	Press the and button to raise or lower settings
	MODE To step through options	Press sequentially to step through and view options available


PATIENT MENU	DISPLAY	OPERATION
For viewing of CPAP settings, humidity, compliance data, altitude and pressure settings		TO VIEW ITEMS Press the "" button to sequentially view:
		Humidity Setting (default display) For adjusting humidity levels to minimize upper airway side effects from CPAP use. Press and buttons to raise or lower settings
	hr	Total run hours Displays the total number of hours the device has been switched on (includes hours without mask use)
	hr/C	Average compliance hours per night* Displays average number of hours the device has been used at night with the mask on
		Checksum* Displays information for medical practitioner or dealer use to check compliance
		Altitude setting Displays altitude setting
	cmH ₂ O	Pressure setting Displays pressure setting

PATIENT MENU - ADDITIONAL SETTINGS	DISPLAY	OPERATION
- To access additional settings: From the Default setting "HUMIDITY" Press the "" key for 3 seconds - To view subsequent items: Press the "" key		Altitude Select "m" for metres or "ft" for feet using the and keys
	3000	Altitude level - Set the altitude level by pressing the and keys - Press to exit menu NOTE: Failure to select correct altitude setting (for any given location) will have an adverse effect on delivered pressure.
	cmH ₂ O	Actual Pressure Displays actual pressure NOTE: Only available on devices with Automatic Altitude Adjustment/Leak Compensation feature.

● NOTE: Manual adjustment of altitude is not featured on units with the Automatic Altitude Adjustment/Leak Compensation feature.


● * Only available on devices with Compliance Monitoring.

8.2 Patient Menu Summary - Integrated CPAP Range

* Only available on devices with Compliance Monitoring.


8.3 Patient Menu - Additional Settings Summary - Integrated CPAP Range

• **Warning:** Failure to select correct altitude setting (for any given location) will have an adverse effect on delivered pressure.

9. CONTROLS AND DISPLAY - STANDARD CPAP RANGE

9.1 Display Description and Functions

KEY	FUNCTION & DESCRIPTION	OPERATION
	ON/OFF Switches the device ON & OFF	TO SWITCH "ON" - Press briefly - The device will default to display the humidity setting TO SWITCH "OFF" - Remove your mask - Press briefly
	RAMP Reduces pressure to a lower level to help you fall asleep. The pressure will gradually return to full operating pressure over 20 minutes	TO START THE RAMP - Press the button NB: to reset ramp during a cycle the device must be switched off first TO SWITCH "OFF" - Press the button
 	INCREASE/DECREASE Adjusts settings up or down	Press the and button to raise or lower settings
	MODE To step through options	Press sequentially to step through and view options available


PATIENT MENU	DISPLAY	OPERATION
For viewing of CPAP settings, humidity, compliance data, altitude and pressure settings		TO VIEW ITEMS Press the "" button to sequentially view:
	23x	Model Type (default display) Displays model type as default setting
	hr	Total run hours Displays the total number of hours the device has been switched on (includes hours without mask use)
	hr/☾	Average compliance hours per night* Displays average number of hours the device has been used at night with the mask on
		Checksum* Displays information for medical practitioner or dealer use to check compliance
		Altitude setting Displays altitude setting
	cmH ₂ O	Pressure setting Displays pressure setting

PATIENT MENU - ADDITIONAL SETTINGS	DISPLAY	OPERATION
- To access additional settings: From the default setting Press the "" key for 3 seconds - To view subsequent items: Press the "" key		Altitude Select "m" for metres of "ft" for feet using the and keys
	3000	Altitude level - Set the altitude level by pressing the and keys - Press to exit menu NOTE: Failure to select correct altitude setting (for any given location) will have an adverse effect on delivered pressure.
	cmH ₂ O	Actual Pressure Displays actual pressure NOTE: Only available on devices with Automatic Altitude Adjustment/Leak Compensation feature.

● NOTE: Manual adjustment of altitude is not featured on units with the Automatic Altitude Adjustment/Leak Compensation feature.


● * Only available on devices with Compliance Monitoring.

9.2 Patient Menu Summary - Standard CPAP Range

* Only available on devices with Compliance Monitoring.

9.3 Patient Menu - Additional Settings Summary - Standard CPAP Range

• **Warning:** Failure to select correct altitude setting (for any given location) will have an adverse effect on delivered pressure.

10. FREQUENTLY ASKED QUESTIONS

• For integrated CPAP devices:

a. When I wake up in the morning, my nose and throat feel dry. What can I do?

Try increasing your humidity setting. If this does not help, please contact your clinician for advice.

b. How do I prevent condensation in the tubing?

The humidity setting enables adjustment of humidity, so that there should be fine misting in the six inches (183 cm) of the breathing tube closest to your face. If the humidity setting is too high for the conditions, condensation may occur in the breathing tube.

Note: Using greater than 6' (183 cm) of tubing will increase condensation.

There are several ways of reducing condensation in the tube:

- i Ensure the device is not positioned in a cool draft.
- ii Decrease the humidity setting on your device slightly. Gradually decrease the setting until the condensation no longer occurs.
- iii Warm the air in the room.

c. Do I have to use distilled water with my device?

The use of distilled water will maximize the life of the water chamber. Water from the faucet will often contain minerals which can damage the chamber, causing pitting in the base, even if the water is passed through a filter.

d. When do I replace my water chamber?

Replace the chamber if the plastic walls of the chamber become cracked or cloudy or the chamber base becomes pitted. (Note: If the chamber leaks at all, it should be replaced immediately).

• Can I use my device in other countries?

Yes. Simply use the appropriate electrical socket adaptor and the device will automatically adapt to any voltage supply.

Note: When moving your device, ensure that the water chamber is empty. Machine failure due to water damage is not covered by warranty.

• Can I use supplemental oxygen with my device?

Yes, oxygen can be administered at the mask. Turn the device on before turning on the oxygen. Ensure that the oxygen source is turned off before turning off the device to avoid oxygen accumulating in the machine. Refer to Section 3 for additional information.

• What happens to my device during power failure?

Remove the mask immediately. Upon restoration of the power supply, the device will restart in the same operation mode with the same settings as before the power failure. Refer to Section 2 for additional information.

11. OPERATING CONDITIONS

AMBIENT TEMPERATURE: 41 - 95 °F (5 - 35 °C)

HUMIDITY: 10 - 95% Relative Humidity (RH)

ALTITUDE: 0 - 9000 ft (0 - 3000 m)

Above 4500 ft (1500 m) the maximum operating pressure will be reduced.

12. STORAGE AND TRANSPORT CONDITIONS

The device should be stored and transported in environmental conditions of: 14 - 140 °F (-10 - 60 °C).

13. TROUBLESHOOTING

If you feel that your device is not operating correctly, consult your healthcare provider.

Do not attempt to repair the device yourself.

14. DEVICE AND CONSUMABLES DISPOSAL INSTRUCTIONS

• Unit Disposal Instructions

This unit contains electronics. Please do not discard as regular rubbish. Dispose according to local guidelines for disposing of electronics.

• Consumable Disposal Instructions

Place the mask, breathing tube and water chamber in a rubbish bag at the end of use and discard with normal rubbish.

15. PRODUCT SPECIFICATIONS

DIMENSIONS:	10.8" x 6.7" x 5.5" (275 mm x 170 mm x 140 mm)
WEIGHT:	4.2 lbs (1.9 kg) (3.0 kg packaged in bag incl. accessories)
PERFORMANCE:	
Pressure Range:	4 to 20 cmH ₂ O (In the unlikely event of fault conditions pressure may reach up to 30 cmH ₂ O)
Altitude Range:	0 to 9000 ft or 0 to 3000 m
Maximum Flow Rates	
CPAP Pressure Setting (cmH ₂ O)	4 8 12 16 20
Measured pressure at the patient connection port (cmH ₂ O)	3 7 11 15 19
Maximum flow rate (L/Min) at mask connection port	60 61 62 66 68
Static Pressure Stability (long term)	
-0.05 to 0.18 cmH ₂ O pressure difference, measured at the mask connection port, at the pressure setting of 10 cmH ₂ O	
Dynamic Pressure Stability (short term)	
CPAP Pressure Setting (cmH ₂ O)	4 8 12 16 20
Pressure Difference (cmH ₂ O) at mask connection port	.65 .74 .59 .63 .53
Humidity:	Maximum Humidity = 43.9 mg/L (BTPS), 100% RH at 4 cmH ₂ O, with humidity setting 4 Typical Humidity = 27.2 mg/L (BTPS), 99.2% RH at 10 cmH ₂ O, with humidity setting 4
Gas Temperature:	Max = 105.8 °F (41 °C)
Noise Level:	<30 dBA

STANDARDS COMPLIANCE:	
Complies with:	EN / IEC 60601-1 AS3200.1.0 UL 60601-1 ISO 17510 -1

ELECTRICAL RATINGS:	
Supply Frequency:	50-60 Hz
Supply Voltage and Current:	1.2 A, 1.3 A ~ (100-115 V max), 0.8 A, 1.8 A ~ (220-240 V max)

NOTE: These values represent average current.

DC to AC Inverter Requirements: 115 V - 200 W/300 W surge
230 V - 300 W/500 W surge

Heater Plate:	85 W max
Heater Plate Temperature:	149 °F (65 °C) max

The device complies with the electromagnetic compatibility requirements of IEC 60601-1-2. In certain circumstances the device may affect or be affected by nearby portable mobile radio frequency communication equipment, due to the effects of electromagnetic interference. If this should happen, try moving your device or the location of the equipment causing interference, or alternatively consult Fisher & Paykel Healthcare.


16. APPENDIX

The following table lists the main difference in features across the SleepStyle™ 200 CPAP Series.

Features	Product Code					
	HC231	HC232	HC233	HC234	HC236	HC238
Standard CPAP	✓	✓			✓	
Fully Integrated CPAP with Heated Humidification			✓	✓		✓
Patient Compliance Monitoring (LCD & Full Download)		✓		✓	✓	✓
Auto Altitude Adjustment					✓	✓
Leak Compensation					✓	✓
Ambient Tracking™ PLUS			✓	✓		✓
Proportional Ramp	✓	✓	✓	✓	✓	✓

Note:

- Some models may not be available for sale in all countries

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