

ADULT BREATHING CIRCUIT REUSABLE HEATED WALL 1.5m

900MR810 / 900MR810E (For use with Fisher and Paykel MR810 humidifier only)
Flow: ≥5 LPM to ≤60 LPM

900MR810E


900MR810

Cleaning/Disinfection:

Use only approved disinfection and cleaning methods as shown:


Step 1. Cleaning:


Step 2. High-level Disinfection:

Use one of the disinfection methods shown:


a) Chemical

b) Pasteurization

c) Autoclave (connectors only)

CE 0123 Rx only

Fisher & Paykel
HEALTHCARE

MANUFACTURER Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.co.nz
Web: www.fphcare.com **AUSTRALIA** (Sponsor) Fisher & Paykel Healthcare Pty Limited, 36-40 New Street, PO Box 167, Ringwood, Melbourne, Victoria 3134, Australia Tel: +61 3 9879 5022 Fax: +61 3 9879 5232 **AUSTRIA** Tel: 0800 29 31 23
Fax: 0800 29 31 22 **BENELUX** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **BRAZIL** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo – SP, Brazil Tel: +55 11 2548 8002 **CHINA** 代理人/售后服务机
构: 费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路 31 号 G12 栋 301 号 电话: +86 20 32053486 传真: +86 20 32052132 **FINLAND** Tel: +358 (0)405 406618 Fax: +46 (0)8 36 6310 **FRANCE**
Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 **GERMANY** Tel: +49 7181 98599 0 Fax: +49 7181 98599 66 **INDIA** Tel: +91 80 4284 4000 Fax: +91 80 4123 6044 **IRISH REPUBLIC** Tel: 1800 409 011 **ITALY** Tel: +39 06 7839 2939 Fax: +39 06 7814
7709 **JAPAN** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **KOREA** Tel: +82 2 6205 6900 Fax: +82 2 6309 6901 **NORTHERN IRELAND** Tel: 0800 132 189 **RUSSIA** Tel and Fax: +7 495 782 21 50 **SPAIN** Tel: +34 902 013 346 Fax: +34 902 013
379 **SWEDEN** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **SWITZERLAND** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **TAIWAN** Tel: +886 2 8751 1739 Fax: +886 2 8751 5625 **TURKEY** Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, Alinteri
Bulvari 1161/1 Sokak No. 12-14, P.O. Box 06371 Ostim, Ankara, Turkey Tel: +90 312 354 34 12 Fax: +90 312 354 31 01 **UK** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU,
UK Tel: +44 1628 626 136 Fax: +44 1628 626 146 **USA/CANADA** Tel: 1800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001

English  Refer to MR810 humidifier user instructions for further set up and correct use information.

Intended Use

The Fisher & Paykel Healthcare 900MR810/900MR810E adult breathing circuits are designed for use with the MR810 humidifier in order to provide therapeutic levels of heat and humidity to a patient's inspired respiratory gases, when using a continuous or intermittent ventilator system or a continuous gas flow system.

Addition of heat and humidity to the supply of cold and dry respiratory gases provided through ventilation is beneficial to prevent drying of the patient airways. The MR810 System is designed for use in hospitals, long term care facilities and homes under the prescription of a qualified medical professional.

Symbol:

Cleaning/Disinfection

Use only approved disinfection and cleaning methods as set out in these instructions:

Discard tubes and all circuit components after 20 cleaning cycles, or six months after first use, whichever occurs earlier.
Do not exceed maximum recommended cleaning temperatures as higher cleaning temperatures may damage the circuit.
Do not use automated washer/cleaner as it may damage the circuit.

Cleaning and disinfection instructions:

Step 1 – Cleaning:

Caution: When cleaning, wear appropriate personal protection equipment, such as: disposable gloves, splash visor/safety glasses and protective gown.

- Dis-assemble the circuit
NOTE: Clean circuit as soon as possible after use.
- Rinse circuit components directly after use with warm tap water (<40°C).
- Fill a sink/container with warm water (<40°C) and Sodium Hypochlorite to 0.5% maximum concentration. Then soak tubing and components for 30–40 minutes.
- Clean the inside of the tube removing all visible debris, if needed, with a soft bottle brush. Do not use hard brushes inside the tube.
- Rinse the inside and outside of the tube thoroughly with warm tap water to remove the remaining Sodium Hypochlorite and debris.
- Hang the tubes vertically upright to remove excess water. Tubes can be air dried up to 75°C.

Warning: Do not clean or dry circuit with solvents or alcohol.

- Proceed to high level disinfection, step 2.

Step 2 – High-level Disinfection:

Use one of the disinfection methods shown:

a. Disinfection with 0.3% orthophthalaldehyde (OPA) chemical:

Caution: Before handling OPA
Wear appropriate personal protection materials, such as: disposable gloves, splash visor/safety glasses and protective gown. Ensure disinfection is performed in a fume cupboard or with fume extraction. Read OPA manufacturers instructions for use.

- Fill sink/container with disinfectant OPA 0.3%, such as Johnson and Johnson Cidex OPA (follow manufacturer instructions for correct preparation of OPA). Submerge tube and components fully in disinfectant. Soak for a minimum of 30 minutes. Then remove all parts from the OPA.
- Change gloves.
- Rinse tubes and components with sterile water (follow manufacturer instructions for correct rinsing procedure of OPA).
- Hang the tubes vertically to let them air dry. Tubes can be air dried up to 75°C. Inspect tube and circuit components, check for any cracks or tears or damage. When all components are dry re-assemble the breathing circuit. Store circuit in a dry and dust free environment ready for next use.

b. Disinfection by Pasteurization:

- Ensure pasteurizer is set to control between 75 and 78°C maximum.
NOTE: Temperatures higher than 78°C will damage the breathing tube and shorten its life.
- Remove circuit components, hang the tubes to dry vertically. Tubes can be air dried up to 75°C. Inspect tube and circuit components, check for any cracks or tears or damage. When all components are dry re-assemble the breathing circuit. Store circuit in a dry and dust free environment ready for next use.

c. Autoclave:

Warning: Only the circuit connectors can be autoclaved. Autoclaving will destroy the tubes.

- Ensure the autoclave is set at either 121°C or 136°C. Leave circuit connectors in autoclave for 30 minutes at 121°C or 4 minutes at 136°C. Inspect circuit connectors, check for any cracks or damage. When all connectors are dry, re-assemble the circuit connectors to disinfected tubing. Store circuit in a dry and dust free environment ready for next use.

Technical Specifications:

Type BF

Interface connections: ISO 5356-1 or ISO 5367-1

Ambient temperature and inlet gas operating range: 18°C – 26°C

Maximum delivered temperature: <40°C

	900MR810	900MR810E
Compliance:	4.6 mL/kPa/m	3.4 mL/kPa/m
Compressible volume:	640 mL	1200 mL
Leakage rate:	<25 mL/min @ 6 kPa	<50 mL/min @ 6 kPa
Resistance to flow at rated flow:	0.4 cmH ₂ O @ 60 ± 1 L/min	1.6 cmH ₂ O @ 60 ± 1 L/min
Flow range:	≥5 to ≤60 L/min *	≥5 to ≤60 L/min *

* Tested with MR370 chamber

Warning:

- Clean circuit prior to use and after each patient use, using approved disinfection methods only.
- Use of unapproved cleaning methods may damage the circuit and reduce its useable life.
- Must be set up and used by a trained medical professional.
- Inspect circuit before re-use, do not use if the circuit shows signs of deterioration, such as: cracks, tears, or damage.
- Use only with Fisher & Paykel MR810, humidification chambers and accessories.
- Ensure that the humidifier is always positioned lower than the patient.
- Do not fill the chamber with water in excess of 37°C.
- Set appropriate ventilator alarms.
- Ensure connector with patient symbol is located at patient end of circuit.
- Use a watertrap in the expiratory limb when applicable.
- Check all tube connections are tight before use.
- Perform a pressure and leak test on the breathing system, and check for occlusions, before connecting to a patient.
- Ensure the complete circuit functions correctly with required ventilator settings, before connecting to a patient.
- Check circuit tubing and watertraps for condensate regularly, and drain as required.
- Do not use the heated breathing circuit without gas flow. If gas flow is interrupted turn the humidifier off.
- Do not cover the circuit with materials such as blankets, towels or bed linen.
- Surface temperature may reach 43°C. Prolonged contact of circuit tube with exposed patients' skin may result in pressure injury or skin damage.
- Do not crush the circuit tube.
- Humidity output of this system can be compromised if used outside recommended ambient temperature or flow range.
- No modification to this equipment is allowed.
- Disconnect tube by handling end connectors only, do not pull or twist tubing, as this may cause damage.
- Check the water level in the humidification chamber regularly and top up as required.
- Do not expose tubes to excessive UV radiation.
- Discard tubes and all circuit components after 20 cleaning cycles, or six months after first use, whichever occurs earlier. At end of product life, disinfect the circuit, then dispose of as uncontaminated medical waste.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

Warranty Statement:

Fisher & Paykel Healthcare warrants that the Evatherm™ circuit, when used in accordance with its instructions for use, shall be free from defects in workmanship and materials and will perform in accordance with Fisher & Paykel Healthcare's official published product specifications for a period of 90 days from the date of purchase by the end user. This warranty is subject to the limitations and exceptions set out in detail in <http://www.fphcare.com/evathermwarranty>

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