
**F&P** HumiGard™

Instructions for Use  
- Laparoscopic & Open Procedures

en	English	1
de	Deutsch	7
fr	Francais	13


## ST310 Humidified Insufflation Kit

### Additional components








Barbed connector

Bypass connector

# SETUP FOR LAPAROSCOPIC & OPEN SURGERY

## TO SET UP

- 1 Remove the pre-assembled, non-sterile components from the tray. Depress finger guard on the humidifier and slide on pre-assembled components (chamber and dry line)
- 2 Remove the filter end from the chamber
- 3 Connect filter end to gas source, using the barbed connector if required;
  - a. laparoscopic surgery: insufflator
  - b. open surgery: CO<sub>2</sub> gas supply stand
- 4 Fill chamber with sterile water to level indicator (~180 mL) using funnel
- 5 Ensure heater wire adapter is connected to the humidifier
- 6 Turn on the humidifier
- 7 Remove the sterile insufflation tube from the pouch
- 8 Connect the heater wire adapter to the insufflation tube elbow
- 9 Connect the blue probe of the heater wire adapter to the insufflation tube elbow. Ensure the blue triangle aligns with the notch for a firm connection.
- 10 Push the insufflation tube elbow onto the chamber
- 11 Connect the luer lock to the patient interface;
  - a. laparoscopic surgery: cannula
  - b. open surgery: VITA-diffuser™

---

## INTENDED USE

For use in warming and humidifying surgical insufflation gas.

### Laparoscopic Procedures

Intended for patients of all ages undergoing laparoscopic procedures where carbon dioxide (CO<sub>2</sub>) insufflation gas is used.

### Open Surgery Procedures

Intended for patients of all ages undergoing open abdominal or cardiothoracic surgical procedures where carbon dioxide (CO<sub>2</sub>) insufflation gas is used.

### Operator Profile

To be operated by a healthcare provider with experience in the set up of surgical procedures, including handling of sterile items. The user must read and understand these instructions for use and the instructions for use of the F&P HumiGard™ Surgical Humidification System (MR860AEU) before use.

# Instructions for Use

## WARNINGS

### General Warnings

- Ensure only CO<sub>2</sub> gas is used. Misconnection of flammable gases could result in serious harm.
- Monitor for systematic acidosis. CO<sub>2</sub> insufflation may increase acidosis.
- Intended for single use. Reuse may result in transmission of infectious substances, interruption of treatment, serious harm or death. Incinerate consumables after use.

### Warnings for Laparoscopic Procedures

- Purge the system of air using CO<sub>2</sub> before insufflation, to prevent air embolism.

### Warnings for Open Surgery Procedures

- Insufficient insufflation may result in inefficient de-airing, monitor for gas embolism in cardiac surgery.
- In-line gas monitoring is required to avoid systemic acidosis in cardiac surgery.

## PRECAUTIONS

### General Precautions

For optimal therapy and performance:

- Inspect sterile packaging for damage before use. Discard if damaged.
- Monitor patient temperature throughout the procedure and adjust patient warming protocols to prevent hyperthermia or hypothermia.
- If using implantables, please consult the manufacturer for suitability of use with humidified CO<sub>2</sub>.
- Only use devices and accessories approved by Fisher & Paykel Healthcare with the ST310 Insufflation Kit. For correct use of approved accessories and devices, follow associated instructions for use.
- Monitor to ensure there is always water in the chamber.

### Precautions for Laparoscopic Surgery Procedures

- For safe and effective use do not exceed the maximum set flow rate of  $\leq 45$  L/min.

### Precautions for Open Surgery Procedures

- For safe and effective use do not exceed the maximum set flow rate of  $\leq 15$  L/min.

## APPROVED DEVICES AND ACCESSORIES

- F&P HumiGard™ Surgical Humidifier (MR860AEU)
- 900ST100 Adaptor

**NOTE:** Accessories may not be available in all countries. Contact your local Fisher & Paykel Healthcare office or authorised distributor for more information.

### Open Surgery Accessory only

- VITA-diffuser™ manufactured by Cardia Innovation AB

## SPECIFICATIONS

### General Specifications

- Humidification chamber reservoir: ~180 mL (use sterile water)
- Luer fittings conform to ISO 594-1:1986 and ISO 594-2:1998
- Not made with natural rubber latex

Fisher & Paykel Healthcare has a policy of continuous product improvement and reserves the right to alter specifications without notice.


Specifications for Laparoscopic Procedures	Specifications for Open Procedures
Maximum input set flow rate: $\leq 45$ L/min	Recommended flow rate: 10 L/min
Maximum achieved flow rate: $\leq 15$ L/min	Maximum input set flow rate: $\leq 15$ L/min
Maximum input set pressure: $\leq 30$ mmHg	Maximum achieved flow rate: $\leq 15$ L/min

## TROUBLE SHOOTING: Veress Needle Use



The use of flow rates  $>1$  L/min when using a Veress needle may result in unstable and inaccurate insufflator pressure readings.

To ensure stable and accurate pressure readings on the insufflator when using the Veress needle for flow rates  $>1$  L/min, the following guidance applies:

- Ensure insufflators are regularly serviced and/or calibrated as per the manufacturer's instructions.
- Ensure that the humidification chamber is filled to the recommended water level of  $\sim 180$  mL, or bypass the chamber using the supplied bypass connector while the Veress needle is in use.


## SYMBOLS


	Sterile		Single use
	Non-sterile		Tray is recyclable
	Excess water in the chamber		Use by date
	Correct water level in the chamber		Date of manufacture
	Inspect sterile packaging for damage before use. Discard if damaged.		Contains phthalate
	Peel sterile bag open in the direction indicated.		

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Ticaret Limited Şirketi,
Fisher & Paykel Healthcare Ltd