





Non-Vented Hospital Full Face Mask with Anti-Asphyxiation Valve

Single use only
Ready for Use

CE 0123 Rx only

Fisher & Paykel
HEALTHCARE

MANUFACTURER Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.co.nz Web: www.fphcare.com

AUSTRALIA (Sponsor) Fisher & Paykel Healthcare Pty Limited, 36-40 New Street, PO Box 167, Ringwood, Melbourne, Victoria 3134, Australia Tel: +61 3 9879 5022 Fax: +61 3 9879 5232 **AUSTRIA** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **BENELUX** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **BRAZIL** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraiso, 04004-000, São Paulo – SP, Brazil Tel: +55 11 2548 8002 **CHINA** 代理人/售后服务机构: 费雪派克医疗保健(广州)有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 传真: +86 20 32052132 **FINLAND** Tel: +358 (0)405 406618 Fax: +46 (0)8 36 6310 **FRANCE** Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 **GERMANY** Tel: +49 7181 98599 0 Fax: +49 7181 98599 66 **INDIA** Tel: +91 80 4284 4000 Fax: +91 80 4123 6044 **IRISH REPUBLIC** Tel: 1800 409 011 **ITALY** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **JAPAN** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **KOREA** Tel: +82 2 6205 6900 Fax: +82 2 6309 6900 **NORTHERN IRELAND** Tel: 0800 132 189 **RUSSIA** Tel and Fax: +7 495 782 21 50 **SPAIN** Tel: +34 902 013 346 Fax: +34 902 013 379 **SWEDEN** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **SWITZERLAND** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **TAIWAN** Tel: +866 2 8751 1739 Fax: +866 2 8751 5625 **TURKEY** Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, Alinteri Bulvarı 1161/1 Sokak No. 12-14, P.O. Box 06371 Ostim, Ankara, Turkey Tel: +90 312 354 34 12 Fax: +90 312 354 31 01 **UK** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: +44 1628 626 136 Fax: +44 1628 626 146 **USA/CANADA** Tel: 1800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001

RT043 Non-Vented Hospital Full Face Mask**Indications for Use:**

The Fisher & Paykel Healthcare RT043 non-vented hospital full face mask is an oronasal patient interface for use as an accessory to therapy devices providing non-invasive bi-level or continuous positive airway pressure support therapy. The therapy device must have adequate alarms and safety systems for device failure. This mask is for single patient use with spontaneously breathing adult (> 30 kg) patients with respiratory insufficiency or respiratory failure who are suitable for non-invasive positive pressure support therapy in a hospital or clinical setting.

This mask includes an anti-asphyxiation valve designed to allow a patient to breathe in case of therapy device failure.

Technical Specification:

Operating pressure range: 3 – 25 cmH₂O
Interface connections: ISO 5356-1 Conical connectors
Compatible with MR850 series non-invasive humidification systems.
Mask dead space <255 cc
This product was not made with natural rubber latex

Label – Mask Size Symbols

(S) = Small (M) = Medium (L) = Large

⚠ Warnings

- Verify that the therapy device, including alarms and safety systems, have been validated prior to use.
- Verify operation of the anti-asphyxiation valve before use. With no system flow the external openings in the valve should allow room air to entrain into the mask. With system flow the flaps should close the external openings and allow system air to flow into the mask. Replace the mask if the anti-asphyxiation valve does not operate correctly or becomes fouled with secretions.
- Only wear mask when therapy is delivered.
- For Bi-Level/CPAP use with a single limb circuit, an exhalation port must be fitted.
- A minimum pressure of 3 cmH₂O is required.
- Set MR850 humidifier to non-invasive mode when administering humidity.
- Remove all packaging before use.
- This mask may only be used in a hospital or clinical setting where the patient is adequately monitored by trained medical staff. Failure to monitor the patient may result in loss of therapy, serious injury or death.
- This mask is fitted with an anti-asphyxiation valve and is not for use with dual limb ventilation systems.
- Do not block or seal the vent hole(s) on the exhalation port.
- In the case of therapy device failure, the use of this mask requires the same level of attention and assistance as in the use of a tracheal tube.
- Reuse may result in transmission of infectious substances, interruption to treatment, serious harm or death.
- Contact the physician if the patient experiences unusual chest discomfort, shortness of breath, stomach distension, belching or severe headache upon waking or when receiving positive airway pressure.
- To minimise the risk of vomiting during sleep, the patient should avoid eating or drinking for three hours before using the mask.
- If the patient experiences skin irritations, consult physician.
- This product is intended for use for a maximum of 7 days. Discard according to hospital protocol.
- California residents please be advised of the following, pursuant to Proposition 65:
This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

Contraindications

- This mask may not be suitable for persons with the following conditions: Impaired cardiac sphincter function, excessive reflux, impaired cough reflex, and hiatal hernia. It should not be used if the patient is vomiting or feeling nauseous, or if the patient is unconscious, uncooperative, unresponsive or unable to remove the mask.
- Do not use on patients with hypotension, secondary to hypovolemia. If symptoms of these conditions occur discontinue treatment immediately.

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