
REF OPT1042 (S)

REF OPT1044 (M)

REF OPT1046 (L)

⊗ MD C €0123 Rx Only

Fisher & Paykel
HEALTHCARE

Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Email: info@fphcare.co.nz Web: www.fphcare.com

Australia (AU) (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 **Brazil (BR)** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo – SP, Brazil Tel: +55 11 2548 7002 **China (CN)** 代理人/售后服务机构:费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12 栋301号 电话: +86 20 32053486 **France (FR)** Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Email: c.s@fphcare.fr **Germany (DE)** Fisher & Paykel Healthcare GmbH, Wiesenstrasse 49, 73614 Schorndorf, Germany Tel: +49 7181 98599 0 **Hong Kong (HK)** Tel: +852 2116 0032 **India (IN)** Tel: +91 80 2309 6400 **Japan (JP)** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea (KR)** Tel: +82 2 6205 6900 **Mexico (MX)** Tel: +52 55 9130 1626 **Poland (PL)** Fisher & Paykel Healthcare Poland Sp. z o.o., Pl. Andersa 7, 61-894 Poznań, Poland Tel: +48 664 846 464 **Russia (RU)** Tel: +7 495 782 21 50 **Switzerland (CH)** Fisher & Paykel Healthcare GmbH, Säntisstrasse 2, 9501 Wil / SG, Switzerland Tel: 0800 83 47 63 **Taiwan (TW)** Tel: +886 2 8751 1739 **Turkey (TR)** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 **UK (GB)** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 **USA (US)/Canada (CA)** Tel: 1800 446 3908 or +1 949 453 4000

Austria (AT) Tel: 0800 29 31 23 **Benelux (BE NL LU)** Tel: +31 40 216 3555 **Denmark (DK)** Tel: +45 70 26 37 70 **Finland (FI)** Tel: +358 9 251 66 123 **Ireland (IE)** Tel: 1800 409 011 **Italy (IT)** Tel: +39 06 7839 2939 **Norway (NO)** Tel: +47 21 60 13 53 **Spain (ES)** Tel: +34 902 013 346 **Sweden (SE)** Tel: +46 8 564 76 680

OPTIFLOW™ 3S NASAL CANNULA

en

INTENDED USE / INDICATIONS FOR USE

The F&P Optiflow 3S nasal cannula is a 14-day single use nasal cannula interface for use with specified respiratory gas humidifiers to deliver nasal high flow (NHF) therapy to spontaneously breathing adult patients. This product is designed to be used by appropriately qualified healthcare professionals in a hospital/institutional environment.

WARNINGS

- Appropriate patient monitoring must be used at all times. Failure to monitor the patient may result in loss of therapy, serious injury or death.
- This device is not intended for life support. Do not use on patients who cannot tolerate a brief interruption of therapy, to avoid serious injury or death.
- Do not use system near a naked flame or any ignition source, including electrosurgery, electrocautery, or laser surgery instruments. Exposure to oxygen increases the risk of fire that may result in patient injury or death.
- Nasal delivery of respiratory gases generates flow-dependent positive airway pressure (PAP). This must be taken into account where PAP could have adverse effects on a patient, to avoid serious injury.
- Failure to set up and use this device as per these instructions can compromise performance and affect safety.

CAUTIONS

- For single patient use only. Reuse may result in transmission of infectious substances.
- Do not over-tighten headstrap, to avoid skin damage.
- Using the cannula may cause skin reddening. Monitor for signs of skin damage or irritation and adjust the cannula if necessary.
- Do not use if all packaging label perforations are torn or if packaging is open, to avoid contamination.
- Achievable flow rates may depend on individual patient and/or flow source.

Failure to comply with the following may result in loss of therapy:

- Do not reprocess (e.g. soak, wash or sterilize) any part of the device as this will cause degradation of materials.
- Do not crush or stretch tube, to avoid damaging the device.
- Properly reinsert the tubing connection after swapping tube entry sides, to prevent leakage.
- Check for condensate regularly. Drain condensate away from patient as required, to prevent blockages or condensate from reaching the patient.
- Do not use with an air entrainer. This will reduce the amount of flow and/or oxygen reaching the patient.

SIDE-EFFECTS

The following side-effects may occur during nasal high flow therapy:

- Skin reddening

SETUP

Airvo™ Series Humidifier with an AirSpiral™ heated breathing tube and chamber kit (e.g. 90OPT561).

Flow range: **OPTI042** 10 – 50 L/min
 OPTI044, OPTI046 10 – 60 L/min

OR

MR850 Respiratory Humidifier in invasive mode, RT series kit with 22 mm heated inspiratory tube and chamber (e.g. RT232).

Flow range: 10 – 60 L/min

OR

F&P 950™ Humidifier in Optiflow mode with F&P 950 adult heated circuit kit (e.g. 950A40).

Flow range: 10 – 70 L/min

Note: The F&P 950 Respiratory Humidifier is currently not cleared for sale in the USA.

Note: The F&P 950 Respiratory Humidifier may not be available in all countries.

SPECIFICATIONS

- Ambient operating temperature: 18 °C to 28 °C (64 °F to 82 °F)
- Permissible transport and storage conditions: –10 °C to 50 °C (14 °F to 122 °F)
- This product is intended to be used for a maximum of 14 days.

BEFORE USE

Before connecting the interface, check for adequate gas flow and ensure that the system has warmed up, to prevent loss of therapy.

Please refer to the WARNINGS and CAUTIONS within these instructions.

FITTING INSTRUCTIONS

- 1

Select an appropriate size. The prongs must not create a seal in the nares, to avoid over pressure. A clear gap must be visible around each prong.
- 2

Use the headgear buckle to apply and remove the cannula from the patient. Place the headgear so that it is not touching the patient's ears, to avoid skin injury.
- 2a

Adjust the headgear to fit.
- 3

The headgear straps can be separated for extra stability.
- 4

To swap tube entry side:

4a

Unclip the cannula tube clip from the side of the cannula.

4b

Push out the tubing connection and reinsert it on the opposite side.

4c

Ensure the tubing connection is properly inserted.

4d

Reconnect the cannula tube clip to the side of the cannula before use.
- 5

Attach the breathing tube clip to a secure location to prevent the cannula from pulling on the ears or off the face.

5a

If using an MR850, ensure the probe cable is not crushed by the breathing tube clip.

DURING USE

- As appropriate for the patient's condition, monitor for any disruptions to the patient receiving flow.
- To prevent condensate, remove/minimize the impact of anything that may cool the tubing, e.g. a fan, air conditioning.
- Notice to User: If a serious incident has occurred while using this device please inform your local Fisher & Paykel Healthcare representative and Competent Authority in your country.

DISPOSAL

- Dispose of products and packaging according to local guidelines. Hospitals should discard cannula according to their standard method for disposing of contaminated product.

SYMBOL DEFINITIONS

	Catalogue number		Prescription only		Correct FITTING INSTRUCTION method
	Manufacturer		Single use		Incorrect FITTING INSTRUCTION method
	Authorised European Representative		Refer to instructions		Polyethylene Terephthalate
	Batch code		Medical device		Small nasal cannula size
	Date and country of manufacture		Storage and transport temperature range		Medium nasal cannula size
	Use-by date		Conforms with Regulation (EU) 2017/745		Large nasal cannula size
	Unique device identifier		Do not open with sharp object		Recommended opening method
	Importer		Distributor		

INTELLECTUAL PROPERTY INFORMATION

F&P, Optiflow, 950, AirSpiral and Airvo are trademarks of Fisher & Paykel Healthcare Limited.

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F&P, Optiflow, 950, AirSpiral, dan Airvo merupakan merek dagang dari Fisher & Paykel Healthcare Limited.