

Configuration for Adult Invasive Ventilation

**F&P™ 850 System:  
Adaptable Humidification**


- Designed for use with adult, pediatric and neonatal patients to humidify invasive ventilation, Optiflow™ high flow therapy and noninvasive ventilation

**Reliable Performance**

- One-touch therapy selection combined with dual temperature feedback systems
- Flow-sensor technology for improved safety and humidity delivery
- Built-in algorithms designed to reduce circuit condensate

**Evaqua™ 2 Circuits Minimize Condensate**

- Insulating Microcell™ technology in the inspiratory limb and permeable Evaqua technology in the expiratory limb reduce condensate
- Saves clinician time and reduces ventilation and expiratory filter issues caused by mobile-circuit condensate
- 14-day duration of use



Some products are unavailable in certain markets