

Configuration for Adult Optiflow™ High Flow Therapy

**F&P™ 850 System:
Adaptable humidification**

- Designed for use with adult, pediatric and neonatal patients to humidify invasive ventilation, Optiflow high flow therapy and noninvasive ventilation.

A Range of F&P Circuit Options


- F&P circuits designed to minimize condensation and maintain humidity delivery from the humidifier to the patient interface

Reliable Performance

- One-touch therapy selection combined with dual temperature feedback systems
- Flow-sensor technology for improved safety and humidity delivery
- Built-in algorithms designed to reduce circuit condensate

Optiflow+ High Flow Interfaces

- The three-way connector fits all F&P circuits
- Customizable head strap provides a secure, soft fit
- Head-strap clip removes the need for a lanyard
- Soft cheek pads designed to add stability and reduce pressure on the face
- Contoured prongs for a softer, streamlined fit



Some products unavailable in certain markets