

**F&P™ 850 System:
Reliable Performance**


- One-touch therapy selection combined with dual temperature feedback systems
- Flow-sensor technology for improved safety and humidity delivery
- Built-in algorithms designed to reduce circuit condensate

**F&P NIV Mask with
Standard Elbow**

- **F&P Nivairo™** Full Face Hospital NIV masks are designed to reduce pressure on the bridge of the nose while maintaining a good seal.
- **F&P Visairo™** Hospital Under Nose masks feature a high-performance seal that eliminates nasal bridge skin breakdown while enabling an open field of vision.

**Evaqua™ 2 Circuits:
Minimize Condensate**

- Insulating Microcell™ technology in the inspiratory limb and permeable Evaqua technology in the expiratory limb reduce condensate
- Saves clinician time and reduces ventilation and expiratory filter issues caused by mobile-circuit condensate
- 14-day duration of use



Some products are unavailable in certain markets