

900PT563 Airvo™ Neb

Nebulize without interruption
to patient therapy.


NEW





Designed for maximum delivery efficiency.

Independent testing* of the Airvo Neb adapter demonstrated an emitted dose as high as 69%†



Better together. The best of Airvo 2 technology in one kit.

This kit comes as standard with our AirSpiral™ heated breathing tube technology, which produces up to 93% less condensate than previous Airvo breathing tubes‡




Built with advanced nebulizer technology in mind.

Designed for use with the latest technology; vibrating mesh nebulizers and the smooth inner surfaces of AirSpiral tubing



Engineered with convenience in mind.

As the Airvo Neb chamber adapter is part of the water chamber there's no need to hold it in place during treatment



Product Code	900PT563	AirSpiral™ tube and chamber kit including Airvo Neb adapter (10 pack) (Aerogen® Solo nebulizer not included)
---------------------	----------	--

For more information, contact your local Fisher & Paykel Healthcare representative at **+1 (800) 792 3912**

This product is only available in the USA.

For further information about the 900PT563, including warnings and precautions, please read the user instructions.

F&P, Airvo, Optiflow and AirSpiral are all trademarks of Fisher & Paykel Healthcare Ltd.

For patent information, see www.fphcare.com/ip

The vibrating mesh nebulizer depicted in this document is the Aerogen Solo. Contact your local Aerogen distributor for information.

Aerogen Solo is a registered trademark of Aerogen Ltd.

*Testing was undertaken in co-operation with Aerogen™.

†Based on the average emitted dose at 10 L/min. Testing carried out with Airvo 2, 900PT563, OPT970 and Aerogen Solo nebulizer. At 20 L/min and 30 L/min the emitted dose was 50% and 37% respectively.

‡When compared with Airvo tube 900PT501 in internal F&P testing.