




Visit our patient resources page

[Learn more](#)

**1 Getting to know your mask**



F&P Eson™ 2 brochure




[Download](#)



Introduction to the F&P Eson 2

[Watch](#)


**2 Getting started**



F&P myMask™ App

With Easy access to information and support, as you need it, F&P myMask App is the must-have mask companion. Access step-by-step videos of common tasks such as mask fitting, fine-tuning, cleaning and assembly.

Scan to Download




**3 Additional support material**



Fitting and cleaning your mask

[Watch](#)



User Instructions

[Download](#)

Replacement schedule

Cushion



Replace every 15 days or if the cushion begins to leak constantly, becomes discolored or is torn.

Frame



Replace every 3 months if the frame becomes discolored or cracked.

Headgear



Replace every 6 months or if the headgear becomes discolored, stretches beyond adjustment or the hook and fasteners no longer latch.