



Visit our patient resources page

Learn more

1

Getting to know your mask



F&P Evora™ Nasal brochure

Download




Introduction to the F&P Evora Nasal

Watch

2

Getting started



F&P myMask™ App

With Easy access to information and support, as you need it, F&P myMask App is the must-have mask companion. Access step-by-step videos of common tasks such as mask fitting, fine-tuning, cleaning and assembly.

Scan to Download



3

Additional support material



Fitting and cleaning your mask

Watch



User Instructions

Download

Replacement schedule

14 Days

Cushion

Replace every 15 days or if the cushion begins to leak constantly, becomes discolored or is torn.



3 Months

Frame

Replace every 3 months if the frame becomes discolored or cracked.

6 Months

Headgear

Replace every 6 months or if the headgear becomes discolored, stretches beyond adjustment or the hook and fasteners no longer latch.