

F&P Solo

AutoFit Nasal Mask
AutoFit Nasal Pillows Mask

Simply AutoFit™



www.fphcare.com

Fisher & Paykel
HEALTHCARE

World's first AutoFit mask

F&P Solo™ combines years of relentless innovation and cutting-edge components to usher in a new era of automatic masks. The world's first AutoFit mask, F&P Solo, simplifies mask setup like no other, using AutoLock™ Technology.



Scan QR code to see online how F&P Solo works



Fit

Users simply stretch the mask over their head and onto their face.

Touch

With one touch, users adjust for their preferred fit by sealing the mask under their nose with a gentle push.

Sleep


AutoLock technology comfortably holds the mask in place, dynamically responding to sleep movement and CPAP-therapy pressures.

One revolutionary headgear design



Reduces Overtightening

Without conventional straps to undo and pull, the AutoFit headgear reduces overtightening.



AutoLock technology

AutoLock technology holds the seal in place and requires only a single pull to release the headgear.

(x-ray view)

Two cushion options

F&P Solo is sold with nasal or nasal pillows cushions. These can also be bought individually and interchanged.



Nasal Cushion




Nasal Pillows Cushion

Simple for you and your patients




Mask Selection

One headgear system with two cushion options to simplify mask choice for you and your patient.




Setup with the F&P myMask™ App

F&P myMask supports your patients with mask setup and fitting.



User-friendly packaging

Step-by-step pictorial instructions provides confidence to set up the mask independently.



F&P myMask support

F&P myMask helps with fitting, cleaning and troubleshooting.

F&P Solo Nasal	
Code	F&P Solo Nasal Description
SLN1SMLU	Nasal Mask Fit Pack
SLN1SU, SLN1SSU*	Nasal Mask Small
SLN1MU, SLN1MMU*	Nasal Mask Medium
SLN1LU, SLN1LLU*	Nasal Mask Large
SLN1WU	Nasal Mask Wide
400SLN111	Nasal Cushion Small
400SLN112	Nasal Cushion Medium
400SLN113	Nasal Cushion Large
400SLN114	Nasal Cushion Wide

F&P Solo Pillows	
Code	F&P Solo Pillows Description
SLP1SMLU	Pillows Mask Fit Pack
SLP1SU, SLP1SSU*	Pillows Mask Small
SLP1MU, SLP1MMU*	Pillows Mask Medium
SLP1LU, SLP1LLU*	Pillows Mask Large
400SLP111	Pillows Cushion Small
400SLP112	Pillows Cushion Medium
400SLP113	Pillows Cushion Large

* Double cushion models. Not available in all countries.

F&P Solo	
Code	F&P Solo Description
400SL121	Headgear Spare
400SL141	Tube, Frame and Swivel Spare

References: 1.85 out of 91 participants rated their experience with the F&P Solo Nasal or F&P Solo Pillows masks as simple or very simple. Internal validation study conducted with 91 participants in USA in 2022. F&P, AutoFit, AutoLock, F&P Solo and myMask are trademarks of Fisher & Paykel Healthcare Limited. For patent information, please see www.fphcare.com/ip. 632012 REV A © 2024 Fisher & Paykel Healthcare Limited.