

REF OJR416HM

REF OJR418HM

Home Nasal Interface




F&P Optiflow™ Junior 2 Home has been designed to deliver Nasal High Flow therapy at home.





F&P Optiflow Junior 2 is a nasal interface specifically designed for the anatomy and flow requirements of infants and children on Nasal High Flow therapy.

It's the next generation of care, providing your small patients with respiratory support.

Like our original Optiflow Junior, Optiflow Junior 2 has soft, anatomically shaped prong design, minimizes condensate and kinking with the FlexiTube and offers easy removal and maintenance with the Wigglepads™ 2.

It uses Waveflex™ technology to allow for natural facial movement when patient's cheeks are compressed. This improves prong stability in the nostrils.



PRODUCT		OJR416HM 	OJR418HM 																				
Nasal cannula size		L	XL																				
Weight guide		 <table><tr><td>Weight (Kg)*</td><td>3</td><td>3.5</td><td>18</td><td>20</td></tr><tr><td>Correlated age**</td><td>37.5 wkGA</td><td>40 wkGA</td><td>4.9 yr</td><td>5.6 yr</td></tr></table>	Weight (Kg)*	3	3.5	18	20	Correlated age**	37.5 wkGA	40 wkGA	4.9 yr	5.6 yr	 <table><tr><td>Weight (Kg)*</td><td>5</td><td>7</td><td>25</td><td>30</td></tr><tr><td>Correlated age**</td><td>47.5 wkGA</td><td>4.7 mo</td><td>7.6 yr</td><td>12 yr</td></tr></table>	Weight (Kg)*	5	7	25	30	Correlated age**	47.5 wkGA	4.7 mo	7.6 yr	12 yr
Weight (Kg)*	3	3.5	18	20																			
Correlated age**	37.5 wkGA	40 wkGA	4.9 yr	5.6 yr																			
Weight (Kg)*	5	7	25	30																			
Correlated age**	47.5 wkGA	4.7 mo	7.6 yr	12 yr																			
myAIRVO 2 flow rate (L/min)		2 - 20	2 - 25																				
PRODUCT SPECIFICATIONS																							
Quantity		Box of 5																					
Usage		Single patient use; Maximum 7 days																					
Materials not present		Not manufactured with natural rubber latex, PVC or Phthalates (DEHP, DBP, BBP)																					
PRODUCT CODES																							
F&P Wigglepad2		WJR112; 2/pack (20 in each pack)																					
Compatible circuits		900PT560; 900PT560E; MYAIRVOKIT1																					
Compatible humidifier		myAIRVO 2 (in Junior mode) PT100xx																					

wkGA = weeks of gestation; mo = months; yr = years

*Weight data is based on F&P product validation studies. **Age data is a correlation to weight data based on a combination of Fenton, WHO and CDC growth charts.

Find out more at
www.fphcare.com/myairvo

Fisher & Paykel
 HEALTHCARE