


## Masks at a glance





**Next generation of Dynamic Support Technology™**

Designed to optimize performance throughout the night.

**Stability wings**

Offer support to the seal.




**Floating seal**

Decoupled from the frame. Soft, pliable silicone seal.

**VentiCool™ Technology**

Breathable headgear designed to allow heat and moisture to dissipate.



**Removable Frame Clip**

Enables the headgear settings to stay in place during fitting and cleaning.

**Exhaust Holes**

Designed for minimal noise.



**Compact Full Face Mask**


Might be right for you if you

Frequently have a stuffy nose or congestion

Want a smaller mask

Might not be right for you if you

Have a sensitive nasal septum




**Dynamic Support Technology**

Designed to keep the seal in place during the night.



**VisiBlue**

**Ease-of-use features:**  
Blue color cues help patients more easily fit, disassemble, assemble and orient the mask.




**VentiCool™ Technology**

Breathable headgear designed to allow heat and moisture to dissipate.

**RollFit™ XT Seal**

Seal dynamically adapts to minimize leaks and pressure on the bridge of the nose.



**Full Face Mask**

Might be right for you if you

Frequently have a stuffy nose or congestion

Might not be right for you if you

Feel claustrophobic

**Soft-knit fabric**

Comfortable headgear designed with AirEdges to avoid leaving marks on the face.



**CapFit™ headgear**

Minimal, one-size headgear for a simple and intuitive fit.



**Floating seal**

Engages to wrap around the nose, providing a flexible and comfortable seal.



**Stability wings**

Designed to keep the floating seal in place throughout the night.



**Compact Nasal Mask**

Might be right for you if you

Don't want something in your nose

Need a mask that is easy to seal and fit

Might not be right for you if you

Frequently have a stuffy nose or congestion



**VisiBlue**

Blue color cues help patients more easily fit, disassemble, assemble and orient the mask.

**Intuitive Headgear**

Headgear – incorporating VisiBlue – can be stretched directly over the head for an easy fit, and then adjusted as required.



**Washable diffuser**

Removable, washable diffuser to muffle the noise and reduce draft.



**RollFit Seal**

Soft silicone cushion, featuring RollFit Technology, automatically adjusts to the bridge of the nose.



**Nasal Mask**

Might be right for you if you

Don't want something in your nose

Need a mask that is easy to seal and fit

Might not be right for you if you

Frequently have a stuffy nose or congestion


**VisiBlue**

Blue color cues help patients more easily fit, disassemble, assemble and orient the mask.



**Washable Diffuser**

The diffuser, which reduces noise and air draft, is washable and durable.



**Adjustable Headgear**

Helps keep the headgear in place as you move.



**AirPillow™ Seal**

The innovative AirPillow inflates to form a gentle and effective seal.

**Nasal Pillows Mask**

Might be right for you if you	Feel claustrophobic
	Want a clear line of sight
	Want a smaller mask
Might not be right for you if you	Have very sensitive nostrils
	Frequently have a stuffy nose or congestion