


Masks at a glance





Compact Full-Face Mask

Might be right for you if you	Frequently have a stuffy nose or congestion
Might be right for you if you	Want a clear line of sight
Might be right for you if you	Want a smaller full-face mask
Might not be right for you if you	Have a sensitive nasal septum

Full-Face Mask


Might be right for you if you	Frequently have a stuffy nose or congestion
Might not be right for you if you	Feel claustrophobic

AutoFit headgear
Simplifies mask setup like no other - simply stretch to fit and touch to adjust.

One touch to adjust
Users fine-tune the mask with one touch for an optimal fit.

AutoLock™ technology
AutoLock technology holds the seal in place and requires only a single pull to release the headgear.

Interchangeable cushions
The AutoFit headgear is interchangeable with F&P Solo Pillows and F&P Solo Nasal cushions.



Nasal Cushion

Nasal Pillows Cushion


*F&P Solo is sold as either a nasal or nasal pillows mask. These cushions can be purchased separately.



Blue color cues help patients more easily fit, disassemble, assemble and orient the mask.



Washable diffuser
Removable, washable diffuser to muffle the noise and reduce draft.



Intuitive Headgear
Headgear - incorporating VisiBlue - can be stretched directly over the head for an easy fit, and then adjusted as required.

RollFit Seal
Soft silicone cushion, featuring RollFit Technology, automatically adjusts to the bridge of the nose.



AutoFit Nasal and AutoFit Nasal Pillows Mask

Might be right for you if you	Want a clear line of sight
Might be right for you if you	Want a smaller mask
Might be right for you if you	Want a mask that is easy to fit and seal
Might not be right for you if you	Frequently have a stuffy nose or congestion

Nasal Mask

Might be right for you if you	Don't want something in your nose
Might not be right for you if you	Frequently have a stuffy nose or congestion



Soft-knit fabric
Comfortable headgear designed with AirEdges to avoid leaving marks on the face.



Stability wings
Designed to keep the floating seal in place throughout the night.



Floating seal
Engages to wrap around the nose, providing a flexible and comfortable seal.



CapFit™ headgear
Minimal, one-size headgear for a simple and intuitive fit.



VisiBlue
Blue color cues help patients more easily fit, disassemble, assemble and orient the mask.



Washable Diffuser
The diffuser, which reduces noise and air draft, is washable and durable.



Adjustable Headgear
Helps keep the headgear in place as you move.



AirPillow™ Seal
The innovative AirPillow inflates to form a gentle and effective seal.

Compact Nasal Mask

Might be right for you if you	Don't want something in your nose
	Want a clear line of sight
	Want a smaller nasal mask
Might not be right for you if you	Need a mask that is easy to fit and seal
	Frequently have a stuffy nose or congestion

Nasal Pillows Mask

Might be right for you if you	Feel claustrophobic
	Want a clear line of sight
	Want a smaller mask
Might not be right for you if you	Have very sensitive nostrils
	Frequently have a stuffy nose or congestion

F&P, AirEdge, AutoFit, AutoLock, Brevida, CapFit, Dynamic Support Technology, Eson, Evora, Solo, F&P Solo, RollFit, VentiCool, VisiBlue and Vitera are trademarks of Fisher & Paykel Healthcare Limited. For patent information, see www.fphcare.com/ip. Prescription only. Use only as directed. Always follow the instructions for use. Your healthcare professional will advise you whether this product is suitable for you/your condition.
632004 REV B © 2024 Fisher & Paykel Healthcare Limited