

myAirvo™ 2 humidifier with integrated flow generator



Introducing the F&P myAirvo 2.

With its integrated flow generator and innovative oxygen delivery system, the myAirvo 2 provides an effective, versatile solution for delivering the benefits of humidification to patients in the home and in long-term care.

High Performance

- Designed to exceed the ISO standard for respiratory humidification
- Quiet operation day and night
- Up to 15 L/min of supplemental oxygen

Versatile

- Variety of patient interfaces for upper and bypassed airways
- Wide flow range: 2-60 L/min













Set up and use

- Graphical display menu for settings
- Integrated O₂ mixing from concentrator or bottle
- Inbuilt O₂ sensor
- No probes or external air supply required

Specifications

myAirvo 2 Humidifier with Integrated Flow Generator	
Intended use	The myAirvo 2 humidifier is for the treatment of spontaneously breathing patients who would benefit from receiving high flow warmed and humidified respiratory gases. This includes patients who have had upper airways bypassed. The flow may be from 2 - 60 L/min depending on the patient interface. The myAirvo 2 is for patients in homes and long-term care facilities.
Humidity output	Complies with ISO 80601-2-74:2017 standard
Flow settings	10 to 60 L/min (default) 2 to 25 L/min (Junior mode)
Patient interfaces	F&P Optiflow+ Series Interfaces F&P Optiflow Junior 2 Home Nasal Cannulas F&P Optiflow Junior Nasal Cannulas
Models	PT100AZ / PT100EE / PT100EW / PT100UK / PT100US / PT100AN / PT100JP / PT100KR
Dimensions	295 mm x 170 mm x 175 mm (11.6" x 6.7" x 6.9")
Weight	2.2 kg (4.8 lb) - unit only 3.4 kg (7.5 lb) - packaged in bag incl. accessories
Electrical ratings	50-60 Hz 100-115 V - 2.2 A (2.4 A max) 220-240 V - 1.8 A (2.0 A max)

Consumables

Product codes	
	Optiflow™+ Series Nasal Cannulas MYOPT9SMALL (2-pack) MYOPT9MEDIUM (2-pack) MYOPT9LARGE (2-pack) OPT942E (1-pack) OPT944E (1-pack) OPT946E (1-pack)
	Optiflow Junior 2 Home Nasal Cannulas OJR416HM (5-pack) OJR418HM (5-pack) Optiflow Junior Nasal Cannulas OPT316 (20-pack) OPT318 (20-pack)
	Optiflow+ Tracheostomy Interfaces MYOPT9TRACHE (2-pack) OPT970E (1-pack)
	Optiflow+ Mask Adapter Interfaces MYOPT9MASK (2-pack) OPT980E (1-pack)
	AirSpiral™ Tube MYAIRSPIRAL (1-pack) 900PT560E (1-pack) 900PT560 (10-pack)
	Water Bag 900PT401 (2-pack)
	AirSpiral Tube and Auto-fill Chamber Kit MYAIRVOKIT1 (1-pack)
	Reusable Chamber MYAIRVOCHAMBER1 (1-pack)
	Compact Stand 900PT400 (1-pack)
	Air Filter 900PT913 (2-pack)
	Oxygen Inlet Kit 900PT422 (1-pack)
	Disinfection Kit 900PT600 (1-pack)

Find out more at

www.fphcare.com/myairvo

Some products may not be available in your country. Please contact your Fisher and Paykel Healthcare representative.