

F&P myMask

Select.
Size.
Assist.

Fisher & Paykel Healthcare's
obstructive sleep apnea mask app

www.fphcare.com

Fisher & Paykel
HEALTHCARE





We understand the ever-increasing demands on you and your resources when it comes to delivering the highest standard of personalized patient care. That's why we have updated F&P myMask App to include **mask selection and sizing capabilities alongside its existing education and help.**

Now you can easily Select, Size and Assist your patients in person or remotely, using proactive support tools to help increase patient competency and mask compliance.


The Selection and Sizing tools are invitation-only services. To register, or for more information, ask your Fisher & Paykel Healthcare representative.

Select.

So many masks to choose from, how do you decide?

Answer a simple series of questions to better understand your patient's needs. Then receive a personalized recommendation of masks along with their features.


In person or remotely, have confidence in knowing that your patient has answered specific questions to identify the recommended masks for you or your patient to select.



Size.

Is sizing a challenge for you or your team?


Scan your patient's face in real-time to discover their mask size. F&P myMask accurately sizes Fisher & Paykel Healthcare PAP masks*, in person or remotely on the phone. Save time, build confidence and have peace of mind knowing that only the patient's mask type and size are being shared.†



Assist.

Are callbacks a costly, repetitive issue for your organization?

Step-by-step personalized education and help are in your hands. Have confidence with access to knowledge and resources on how to effectively fit, fine-tune and clean Fisher and Paykel Healthcare PAP masks, along with troubleshooting and support in the Help menu. All this information is in one place, available at your convenience.



* The F&P myMask App sizing accuracy was tested on 73 participants in 2023. When compared to the use of calipers, the F&P myMask App was found to have a median accuracy of > 91% (Confidence Interval 90 - 95%) across three measures: Face height, nose height and nose width.

† We take patient data and privacy seriously and respect anonymity. There is no need for a patient to create an account to use any of the features on the myMask App including gaining access to selection and sizing or for mask help and support. No patient images are transmitted during the myMask selection and sizing process and no image or video is stored on the device after processing. To manage data security, we are ISO 27001 certified internationally and ISO 13485 certified in New Zealand.

Introducing F&P Services

This user-friendly website enables you to invite clinicians and patients to access and utilize our new Selection and Sizing tools.

Patient mode access enables:

Your personal Dashboard

Personalize your contact details for sharing with your patients.

Track patients who have used your app invite to determine what mask size is suitable.

Patient-text Management

Invite clinicians and patients to utilize the latest technology for mask recommendations and sizing, at their convenience.

To register for these new services, contact your Fisher & Paykel Healthcare representative.



Patient mode remote



Patients receive a text-message invitation to download app



Download app



Mask recommendations



Size recommendation



Clinician mode in person



Clinician invitation to download F&P myMask App



Mask recommendations



Size recommendation



F&P myMask App supports these PAP masks



F&P Evora™ Full
Compact Full-Face Mask



F&P Evora™ Nasal
Compact Nasal Mask



F&P Vitera™
Full-Face Mask



F&P Eson™ 2
Nasal Mask



F&P Simplus™
Full-Face Mask



F&P Brevida™
Nasal Pillows Mask



For a full product set-up description, cautions, hazards and warnings, refer to user instructions provided with the product. Use device as directed. Always follow the instructions for use, e.g., the user instructions. F&P, Brevida, Eson, Evora, myMask, Simplus and Vitera are trademarks of Fisher & Paykel Healthcare Limited. For patent information, please see www.fphcare.com/ip.
PM-630326 REV B © 2023 Fisher & Paykel Healthcare Limited