

Infant Resuscitation Mask Sizing Guide



	RD803-10 Micro-Premature (XS)	RD804-10 Premature (S)	RD805-10 Neonatal (M)	RD806-10 Infant (L)	RD807-10 Pediatric (XL)
Diameter	35 mm (1.38")	42 mm (1.65")	50 mm (1.97")	60 mm (2.36")	72 mm (2.83")
Suggested Weight Range*	400 g - 1 kg (0.88 - 2.2 lb)	800 g - 1.5 kg (1.8 - 3.3 lb)	1.5 - 2.5 kg (3.3 - 5.5 lb)	2.5 - 5 kg (5.5 - 11 lb)	5 - 10 kg (11 - 22 lb)
	Extremely Low Birth Weight (ELBW)	Very Low Birth Weight (VLBW)	Low Birth Weight (LBW)	Term	Pediatric
Material	Thermoplastic Elastomer (The product does not contain any latex or DEHP plasticizer).				
Cleaning and Disposal	Single patient use - supplied clean, not sterile. Discard according to hospital protocol.				

* This is only to be used as a guide. Actual measurements should be taken for an ideal fit and optimum therapy delivery.