

**SwingFit™.**  
Performance defined.



# Discover the definition of performance

We know that for a sound sleep, mask performance is crucial. This is why our industry-leading engineers examined mask design from every angle to create a mask that you and your patients can trust. The result: F&P Nova™ Nasal - it's performance defined.

**1**  
Easy fitting

Performance  
you can trust

98%

of trial users trust  
F&P Nova Nasal  
to work<sup>1</sup>

**2**  
Stable seal

**3**  
Sound sleep

# Combining innovative technologies

Packed with features working in synergy with the intuitive SwingFit headgear, F&P Nova Nasal delivers all-round performance. Discover a mask that's easy to fit and provides a stable seal, with a clear field of view and minimal noise and draft.

**1**  
SwingFit headgear

Simply anchor the cushion on the nose and swing the headgear around the head to secure with one click.

**AirEdge**  
Soft fabric edge helps avoid marks on the face.

**3**  
Diffuser  
Minimizes noise and draft.

**Tube clip**  
Optional over-the-head tube placement.

**Clear field of view**  
Structured side arms for visual freedom.

**2**  
RollFit™ cushion with pivoting clips  
Adapts to movement for a stable seal.



# F&P Nova Nasal

## Product Codes

Code	Description
NVN1SZ	Nova Nasal Mask - Small
NVN1MZ	Nova Nasal Mask - Medium
NVN1LZ	Nova Nasal Mask - Large
NVN1SMLZ	Nova Nasal Mask - Fit pack - Small, Medium, Large


## Spare Parts

Code	Description
400NVN111	Nova Nasal Spare Cushion - Small
400NVN112	Nova Nasal Spare Cushion - Medium
400NVN113	Nova Nasal Spare Cushion - Large
400NVN121	Nova Nasal Headgear - Small
400NVN122	Nova Nasal Headgear - Medium-Large
400NVN123	Nova Nasal Back Panel - Small*
400NVN124	Nova Nasal Back Panel - Medium-Large*
400NVN141	Nova Nasal Frame Assembly Spare
400NVN151	Nova Nasal Diffuser Spare
400NVN161	Nova Nasal Tube Clip Spare
400NVN162	Nova Nasal Headgear Pad Spare



## F&P myMask™ App

Instructs your patients on mask use with step-by-step fitting, fine-tuning and cleaning videos. Download the F&P myMask App. Scan the QR code or download directly.



\* Headgear includes back panel. Back panel can also be purchased separately.

**Always read the label and follow the directions for use.** Prescription only. This product may not be right for you. Your healthcare provider will advise you whether this product is suitable for you/your condition. If symptoms persist, talk to your healthcare provider.

1.42 out of 43 users agreed or strongly agreed with the following statement: "I trust the F&P Nova Nasal mask to work." Internal validation trial conducted by Fisher & Paykel Healthcare on 43 participants in the United States, 2023.

F&P, F&P Nova, Nova, SwingFit, RollFit and myMask are trademarks of Fisher & Paykel Healthcare Limited. For patent information, please see [www.fphcare.com/ip](http://www.fphcare.com/ip). App Store is a service mark of Apple Inc., registered in the US and other countries and regions. Google Play and the Google Play logo are trademarks of Google Inc. 633606 REV A © 2025 Fisher & Paykel Healthcare Limited.

**Fisher & Paykel**  
HEALTHCARE