

Stable. Simple. Soft.
Optimized therapy delivery.



Introducing **F&P Optiflow™ 3S**

New nasal cannula interface range for Optiflow Nasal High Flow therapy.

The Optiflow 3S aims to provide a positive patient and clinician experience.



Stable

Enhanced Patient Compliance

A secure fit ensures consistent delivery of therapy. The Optiflow 3S:

- Offers stability, resulting in **reduced interface displacement**
- Allows the patient to **eat, drink, and communicate** without interruption to therapy, therefore improving compliance



Simple

Patient-Centered Care

An interface that is user-focused enhances patient care efficiency and enables a positive patient experience. The Optiflow 3S:

- Provides **flexibility** in setup
- Simplifies the **fitting process**
- Is compatible with a range of Fisher & Paykel Healthcare devices and circuits, streamlining the patient journey



Soft

Uncompromised Patient Comfort

A comfortable patient is a compliant patient, and patient compliance ultimately improves outcomes. The Optiflow 3S:

- Provides a **comfortable fit**
- Is designed to **minimize pressure** on the face and head

Features

Tube:

Consists of durable, **breathable** Evaqua™ material that reduces mobile condensate

Headstrap management:

- **Soft**, breathable fabric wicks moisture away from the head
- Tangle-free straight out of the packet
- Headstrap designed for **stability**
- Quick-release buckles ease the **fitting process**
- Hook tabs help to keep straps off the face



Cannula design:

Cannula arms are encased in soft material and are designed to hug the natural curves of the face to provide **comfort**

Structured yet flexible arms provide a **secure** fit without the need to over-tighten

Pillow feature designed to **minimize pressure** on the nose and upper lip

Connector:

Fisher & Paykel Healthcare's universal connector is color coded to provide swift sizing identification to help ease the fitting process

Tube clip:

Supports the weight of the circuit to keep prongs in place, which results in **less interface displacement**

Neatly positions the tube to the side of the patient's face

Side-swapping:

Provides setup **flexibility** for clinicians by allowing delivery of therapy from either side of the patient

Specifications

F&P Optiflow 3S	Product Order Codes
F&P Optiflow Small	OPT1042
F&P Optiflow Medium	OPT1044
F&P Optiflow Large	OPT1046
Box Quantity	20
Usage	Single-patient use; maximum 14 days
Materials	DEHP free

www.fphcare.com/optiflow3S

This product may not be available in your country. Please contact your local Fisher & Paykel Healthcare representative. For further information about the Optiflow 3S, including warnings and precautions, please read the user instructions. F&P, Optiflow, and Evaqua are all trademarks of Fisher & Paykel Healthcare Ltd. For patent information, see www.fphcare.com/ip

