


The Optiflow™ Junior 2+ nasal cannula has been designed to fit a wider range of patients, bridging the gap between the XL Optiflow Junior 2 cannula and the Small Optiflow+ cannula. It retains all the existing benefits of the Optiflow Junior 2 cannula range including Waveflex™ technology to assist with prong stability and Wigglepads™ 2 to allow for easy application, adjustment and maintenance of the cannula.



**Increased Flow**

Larger tubing to provide more flow from 1-50 L/min.\*



**Wider Patient Range**

Designed to fit patients between 6 months and 12 years old.



\* Depending on platform

## OJR520 Nasal Cannula Specifications



### OJR520 XXL

### MR850

### AIRVO™ 2

#### COMPATIBLE ACCESSORIES

F&P Wigglepad 2

WJR114



#### MAXIMUM FLOW RATES\* (L/min)

	1 – 36	10 – 50
--	--------	---------

#### PRODUCT SPECIFICATIONS

Compatible circuit	RT330, OPT014	900PT531, 900PT561
Cannula weight	28.2 g	
Connection	Easy-click connector	
Quantity	Box of 10	
Box components	F&P Optiflow Junior 2+ nasal cannula, User Instructions	

#### PERFORMANCE

Ambient temperature range	18 – 26 °C / 64 – 79 °F	
Humidifier mode	Invasive mode	Default mode
Usage	Single patient use; maximum 7 days	
Duration of use	7 days	
Recommended gas source	Medical gas	

#### COMPONENTS AND COMPOSITIONS

Predominant materials	Thermoplastic elastomer; hydrocolloid; ABS	
Materials not present	Not manufactured with natural rubber latex, PVC or phthalates (DEHP, DBP, BBP)	
Manufacturing mode	Noninvasive device; produced in a controlled working environment	
Disposal	Incineration or according to hospital protocol for cannula; clamshell packaging and label recyclable PET	

#### REGULATORY

Classification	AU-IIa, EU-IIa, Canada-II, US-II	
Country of origin	New Zealand	
Notified body	TÜV SÜD Product Service GmbH CE0123	

\* Flow rates above describe technical capability of the product when used at sea level. Ensure clinical judgement is used when prescribing flow rates.

Please note that the information in this specifications sheet (including product information and images) is summarized and provided for illustrative purposes only. Please refer to the relevant User Instructions for more information and confirm details with your local Fisher & Paykel Healthcare representative prior to placing an order. Information subject to change without notice.

F&P, AIRVO, Optiflow, Wigglepads, and Dolphin are trademarks of Fisher & Paykel Healthcare Limited.  
For patent information, see [www.fphcare.com/ip](http://www.fphcare.com/ip). PM-619584 REV C © 2022 Fisher & Paykel Healthcare Limited