

Choosing the right


F&P Nasal Mask



F&P Evora™



F&P Eson™ 2



F&P Eson 2

Eson 2 is a stable, over-the-nose mask.
It's secure on the face with a narrow,
transparent stability bar.

The RollFit seal works in conjunction with
the stability bar to create dynamic support
during movement.

Available in sizes: Small (ESN2SA),
Medium (ESN2MA), and Large (ESN2LA).



Compact Nasal Mask



F&P Evora

Evora sits under the nose with minimal headgear on the face. It's simple and intuitive due to the CapFit headgear that can be taken on and off in one quick movement.

The stability wings work in conjunction with the floating seal to create our next generation of Dynamic Support Technology™.

Available in sizes: Small (EVO1SA), Medium (EVO1MA), Large (EVO1LA), and Wide (EVO1WA).