


# F&P Full Face Mask

Stability for active sleepers

F&P Vitera™

## F&P myMask App

Designed to support mask setup, the F&P myMask™ app helps you to effectively fit, fine-tune and clean your mask.\* Support when it matters most.



Scan to download or download directly  
\* F&P myMask is available in selected countries only

## F&P Evora Full

### Product Codes

Code	Description
EVFIXA	Evora Full Mask - Extra Small
EVFIMA	Evora Full Mask - Small to Medium
EVFILA	Evora Full Mask - Large
EVFIXMLA	Evora Full Mask Fitpack - Extra Small, Small to Medium, Large

### Spare Parts

Universal frame fits all seal sizes

Code	Description
400EVF114	Evora Full Seal - Extra Small Spare
400EVF115	Evora Full Seal - Small to Medium Spare
400EVF116	Evora Full Seal - Large Spare
400EVF121	Evora Full Headgear - Standard Spare
400EVF122	Evora Full Headgear - Extra Large Spare
400EVF131	Evora Full Mask - No Headgear, Extra Small Seal
400EVF132	Evora Full Mask - No Headgear, Small to Medium Seal
400EVF133	Evora Full Mask - No Headgear, Large Seal
400EVF141	Evora Full Headgear Clips & Frame Clip Spare
400EVF151	Evora Full Tube Spare
400EVF161	Evora Full Frame Spare
400FPH341	F&P Swivel White Spare

## F&P Vitera

### Product Codes

Code	Description
VITISA	Vitera Mask - Small
VITIMA	Vitera Mask - Medium
VITILA	Vitera Mask - Large
VITISMLA	Vitera Mask Fitpack - Small, Medium, Large

### Spare Parts

Universal frame fits all seal sizes

Code	Description
400VIT114	Vitera Seal - Small Spare
400VIT115	Vitera Seal - Medium Spare
400VIT116	Vitera Seal - Large Spare
400VIT121	Vitera Headgear - Small Spare
400VIT122	Vitera Headgear - Medium/Large Spare
400VIT131	Vitera Mask - No Headgear, Small Seal
400VIT132	Vitera Mask - No Headgear, Medium Seal
400VIT133	Vitera Mask - No Headgear, Large Seal
400FPH241	F&P Swivel
400VIT151	Vitera Clips and Forehead Clip
400VIT171	Vitera Elbow

F&P, Dynamic Support Technology, Evora, myMask, RollFit, VentiCool, VisiBlue and Vitera are trademarks of Fisher & Paykel Healthcare Limited. For patent information, see www.fphcare.com/ip. Apple, the Apple logo and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Google Play and Google Play logo are trademarks of Google LLC. 627655 REV A © 2022 Fisher & Paykel Healthcare Limited

# Choosing the right

Minimal contact and clear line of sight

F&P Evora™ Full

**Universal Full Face Mask**



**F&P Vitera**  
**Redefining Mask Performance**


F&P Vitera incorporates Dynamic Support Technology™ to provide stability and comfort throughout the night: the low profile stability bar works in conjunction with the RollFit XT cushion allowing the cushion to adjust dynamically to patient movements.

The RollFit XT cushion also adapts to minimize leaks and pressure on the bridge of the nose.

Available in sizes: Small (VIT1SA), Medium (VIT1MA), Large (VIT1LA), Fitpack (VIT1SMLA).



**Compact Full Face Mask**



**F&P Evora Full**  
**Performance Made Compact**

F&P Evora Full is a compact, under-the-nose full face mask. It provides minimal contact with full performance. The mask features the next generation of Dynamic Support Technology; a combination of a 'Floating Seal', supported by Stability Wings. These technologies work together to allow for freedom of movement while keeping the mask comfortably in place.

Extra Small (EVF1XA), Small to Medium (EVF1MA), Large (EVF1LA) and Fitpack (EVF1XMLA).

