

## RT330 Optiflow Tubing Kit

### RT330 CIRCUIT, CHAMBER AND PRESSURE MANIFOLD

#### PERFORMANCE SPECS

<b>Max Flow Rate (L/Min)</b>	0.5 – 25
<b>Length of Tubing</b>	1.75m/5.7ft
<b>Ambient Range °C</b>	18 – 26 °C
<b>Connection</b>	Unique Swivel connection to Optiflow Cannula (OPT31X)
<b>Compatible Humidifier</b>	MR850
<b>Humidifier Mode</b>	Invasive
<b>Duration of Use</b>	7 days
<b>Carton Size</b>	Box of 10
<b>Recommended Gas Source</b>	Medical Gas
<b>Pressure Manifold</b>	40cm H <sub>2</sub> O pressure pop off

#### COMPONENTS & COMPOSITION

<b>Predominant materials</b>	High Density Polyethylene; Polypropylene; PTFE
<b>Materials Not Present</b>	Contains no natural rubber Latex, PVC or Phthalates (DEHP, DBP, BBP)
<b>Manufacturing Mode</b>	Noninvasive device; Produced in CWE
<b>Compatible With</b>	OPT31X Range
<b>Destruction</b>	Incineration, or according to hospital protocol

#### REGULATORY

<b>Classification</b>	AU-IIa, EU-IIa, Canada-II, USA-II
<b>Country of Origin</b>	New Zealand
<b>Intended Use</b>	Single patient use; Max 7 days
<b>EU Notified Body</b>	TÜV SÜD Product Service GmbH CE0123

