

Simply AutoFit™



World's first AutoFit mask

F&P Solo™ combines years of relentless innovation and cutting-edge components to usher in a new era of automatic masks. The world's first AutoFit mask, F&P Solo, simplifies mask setup like no other, using AutoLock™ Technology.



Scan QR code to see online
how F&P Solo works



Fit

Users simply stretch the mask over their head and onto their face.

Touch

With one touch, users adjust for their preferred fit by sealing the mask under their nose with a gentle push.

Sleep

AutoLock technology comfortably holds the mask in place, dynamically responding to sleep movement and CPAP-therapy pressures.




One revolutionary headgear design



Reduces Overtightening

Without conventional straps to undo and pull, the AutoFit headgear reduces overtightening.



AutoLock technology

AutoLock technology holds the seal in place and requires only a single pull to release the headgear.


(x-ray view)

Two cushion options

F&P Solo is sold with nasal or nasal pillows cushions. These can also be bought individually and interchanged.



Nasal Cushion



Nasal Pillows Cushion




Simple for you and your patients



Mask Selection

One headgear system with two cushion options to simplify mask choice for you and your patient.




Setup with the F&P myMask™ App

F&P myMask supports your patients with mask setup and fitting.



User-friendly packaging

Step-by-step pictorial instructions provides confidence to set up the mask independently.



F&P myMask support

F&P myMask helps with fitting, cleaning and troubleshooting.

| F&P Solo Nasal | |
|----------------|----------------------------|
| Code | F&P Solo Nasal Description |
| SLN1SMLU | Nasal Mask Fit Pack |
| SLN1SU | Nasal Mask Small |
| SLN1MU | Nasal Mask Medium |
| SLN1LU | Nasal Mask Large |
| SLN1WU | Nasal Mask Wide |
| 400SLN115 | Nasal Cushion Small |
| 400SLN116 | Nasal Cushion Medium |
| 400SLN117 | Nasal Cushion Large |
| 400SLN118 | Nasal Cushion Wide |

| F&P Solo Pillows | |
|------------------|------------------------------|
| Code | F&P Solo Pillows Description |
| SLP1SMLU | Pillows Mask Fit Pack |
| SLP1SU | Pillows Mask Small |
| SLP1MU | Pillows Mask Medium |
| SLP1LU | Pillows Mask Large |
| 400SLP114 | Pillows Cushion Small |
| 400SLP115 | Pillows Cushion Medium |
| 400SLP116 | Pillows Cushion Large |

| F&P Solo | |
|----------|------------------------------|
| Code | F&P Solo Description |
| 400SL121 | Headgear Spare |
| 400SL141 | Tube, Frame and Swivel Spare |

References: 1. 85 out of 91 participants rated their experience with the F&P Solo Nasal or F&P Solo Pillows masks as simple or very simple. Internal validation study conducted with 91 participants in USA in 2022. F&P, AutoFit, AutoLock, F&P Solo and myMask are trademarks of Fisher & Paykel Healthcare Limited. For patent information, please see www.fphcare.com/ip. 632007 REV A © 2024 Fisher & Paykel Healthcare Limited.