

World's first AutoFit™ mask



F&P Solo Nasal



F&P Solo Pillows



93%

of trial users rated their experience with either the F&P Solo™ Nasal or F&P Solo Pillows mask as simple or very simple.¹



F&P Solo

World's first AutoFit mask

Fisher & Paykel Healthcare's world-first AutoFit technology is available in nasal and nasal pillows. F&P Solo Nasal and F&P Solo Pillows simplify setup like no other, using AutoLock™ technology – stretch to fit, touch to adjust.



Scan QR code for the online F&P Solo experience.

Fisher & Paykel
HEALTHCARE


Fit

Users simply stretch the mask over their head and onto their face.



Touch

With one touch, users adjust for their preferred fit by sealing the mask under their nose with a gentle push.



AutoLock technology holds the seal in place and requires only a single pull to release the headgear.


Sleep

AutoLock technology comfortably holds the mask in place, dynamically responding to sleep movement and CPAP-therapy pressures.



One revolutionary headgear design

Simple,
every night



Users can easily fine-tune the mask during nightly use for an optimal fit.



X-ray view

Quality
components



AutoLock technology,
manufactured with carbon fiber
material for seamless movement.

Reduces
overtightening



Without conventional straps
to undo and pull, the AutoFit
headgear reduces overtightening.

Two comfortable cushion options, one headgear

Nasal Cushion

F&P Solo Nasal features Dynamic Support Technology in the nasal cushion, allowing freedom of movement while keeping the cushion comfortably in place.


Pillows Cushion

F&P Solo Pillows offers the next generation of AirPillow™ cushion, which comfortably sits in and around the nose for a gentle, stable seal.

The AutoFit headgear is interchangeable with either the nasal or nasal pillows cushions. Cushion types sold separately.




Simplify the mask experience across the patient's journey




Mask Selection

The AutoFit headgear is interchangeable with F&P Solo Pillows and F&P Solo Nasal cushions.




Setup with the F&P myMask™ App

For an enhanced experience, F&P myMask supports your patients with mask setup and fitting.




User-friendly packaging

Step-by-step pictorial instructions provides confidence to set up the mask independently.




F&P myMask support

F&P myMask helps with fitting, cleaning and troubleshooting. Support when it's needed most.



Simple for you and your patients

F&P Solo has done the hard work so you and your patients don't have to.



F&P Solo Nasal Product Codes

Code	Description
SLNISA	Nasal Mask Small A Model
SLNIMA	Nasal Mask Medium A Model
SLNILA	Nasal Mask Large A Model
SLNIWA	Nasal Mask Wide A Model
SLNISMLA	Nasal Mask Fitpack (S, M, L) A Model
400SLN115	Nasal Cushion Small Spare
400SLN116	Nasal Cushion Medium Spare
400SLN117	Nasal Cushion Large Spare
400SLN118	Nasal Cushion Wide Spare
400SLN131	Nasal Mask No Headgear, Small
400SLN132	Nasal Mask No Headgear, Medium
400SLN133	Nasal Mask No Headgear, Large
400SLN134	Nasal Mask No Headgear, Wide

F&P Solo Spare Parts

(Shared between F&P Solo Nasal and F&P Solo Pillows)

Code	F&P Solo Description
400SL122	Headgear Spare
400SL141	Tube, Frame and Swivel Spare

F&P Solo Pillows Product Codes

Code	Description
SLPISA	Pillows Mask Small A Model
SLPIMA	Pillows Mask Medium A Model
SLPILA	Pillows Mask Large A Model
SLPISMLA	Pillows Mask Fitpack (S, M, L) A Model
400SLP114	Pillows Cushion Small Spare
400SLP115	Pillows Cushion Medium Spare
400SLP116	Pillows Cushion Large Spare
400SLP131	Pillows Mask No Headgear, Small
400SLP132	Pillows Mask No Headgear, Medium
400SLP133	Pillows Mask No Headgear, Large

References: 1. 85 out of 91 participants rated their experience with the F&P Solo Nasal or F&P Solo Pillows mask as "simple" or "very simple". Internal validation trial conducted by Fisher & Paykel Healthcare on 91 participants in the United States, 2022. 2. 79 out of 83 participants were able to set themselves up on the F&P Solo Nasal or F&P Solo Pillows mask without assistance. Internal validation trial conducted by Fisher & Paykel Healthcare on 83 participants (naïve to PAP therapy) in New Zealand, 2022. F&P, AirPillow, AutoFit, AutoLock, F&P Solo, Solo and myMask are trademarks of Fisher & Paykel Healthcare Limited. For patent information, please see www.fphcare.com/ip. 632993 REV A © 2024 Fisher & Paykel Healthcare Limited.