

F&P Solo Nasal

AutoFit Nasal Mask

Simply AutoFit™



www.fphcare.com

Fisher & Paykel
HEALTHCARE

World's first AutoFit mask

F&P Solo™ combines years of relentless innovation and cutting-edge components to usher in a new era of automatic masks. The world's first AutoFit mask, F&P Solo Nasal, simplifies mask setup like no other, using AutoLock™ Technology.



Fit

Users simply stretch the mask over their head and onto their face.

Touch

With one touch, users adjust for their preferred fit by sealing the mask under their nose with a gentle push.

Sleep

AutoLock technology comfortably holds the mask in place, dynamically responding to sleep movement and CPAP-therapy pressures.

A revolutionary headgear design




Reduces overtightening

Without conventional straps to undo and pull, the AutoFit headgear reduces overtightening.



Comfortable nasal cushion


Dynamic Support Technology allows freedom of movement while keeping the cushion comfortably in place.



AutoLock technology

AutoLock technology holds the seal in place and requires only a single pull to release the headgear.

Simple for you and your patients




Setup with the F&P myMask™ App

F&P myMask supports your patients with mask setup and fitting.



User-friendly packaging

Step-by-step pictorial instructions provides confidence to set up the mask independently.



F&P myMask support

F&P myMask helps with fitting, cleaning and troubleshooting.

Product Codes

Code	Description
SLNISA	Solo Nasal Mask Small A Model
SLN1MA	Solo Nasal Mask Medium A Model
SLN1LA	Solo Nasal Mask Large A Model
SLN1WA	Solo Nasal Mask Wide A Model
SLN1SMLA	Solo Nasal Mask Fit Pack (S, M, L) A Model

Spare Parts

Code	Description
400SL122	Solo Headgear Spare
400SL141	Solo Tube, Frame and Swivel Spare
400SLN115	Solo Nasal Cushion - Small Spare
400SLN116	Solo Nasal Cushion - Medium Spare
400SLN117	Solo Nasal Cushion - Large Spare
400SLN118	Solo Nasal Cushion - Wide Spare
400SLN131	Solo Nasal Mask - No Headgear, Small
400SLN132	Solo Nasal Mask - No Headgear, Medium
400SLN133	Solo Nasal Mask - No Headgear, Large
400SLN134	Solo Nasal Mask - No Headgear, Wide

References: 1. 42 out of 45 participants rated their experience with the F&P Solo Nasal mask as simple or very simple. Internal validation study conducted with 45 participants in USA in 2022. F&P, AutoFit, AutoLock, F&P Solo and myMask are trademarks of Fisher & Paykel Healthcare Limited. For patent information, please see www.fphcare.com/ip. 632011 REV A © 2024 Fisher & Paykel Healthcare Limited.