

F&P Solo

AutoFit Nasal Mask

Simply AutoFit™



www.fphcare.com/Solo

Fisher & Paykel
HEALTHCARE

World's first AutoFit mask

F&P Solo™ combines years of relentless innovation and cutting-edge components to usher in a new era of automatic masks.

The world's first AutoFit mask, F&P Solo, simplifies mask setup like no other, using AutoLock™ Technology.



Fit

Users simply stretch the mask over their head and onto their face.



Touch

With one touch, users adjust for their preferred fit by sealing the mask under their nose with a gentle push.



AutoLock technology holds the seal in place and requires only a single pull to release the headgear.


Sleep

AutoLock technology comfortably holds the mask in place, dynamically responding to sleep movement and CPAP-therapy pressures.



One revolutionary headgear design

Simple,
every night



Users can easily fine-tune the mask during nightly use for an optimal fit.



Quality
components



AutoLock technology,
manufactured with carbon fiber
material for seamless movement.

Reduces
overtightening



Without conventional straps
to undo and pull, the AutoFit
headgear reduces overtightening.

Comfortable nasal cushion



Nasal Cushion

Sits under the nose allowing for freedom of movement while keeping the cushion comfortably in place.



S



M



L



W

Simplify the mask experience across the user's journey




Setup with the F&P myMask™ App

For an enhanced experience, F&P myMask supports your patients with mask setup and fitting.




User-friendly packaging

Step-by-step pictorial instructions provides confidence to set up the mask independently.



F&P myMask App support

F&P myMask helps with fitting, cleaning, and troubleshooting. Support when it's needed most.



Simple for you and your patients

F&P Solo has done the hard work so you and your patients don't have to.



Headgear (one size)

Nasal cushion

Frame + Breathing Tube + Swivel



Under our F&P SureFit™ 30-Day Mask Guarantee if you are not completely satisfied with your F&P Solo mask, simply return it to the store you purchased it from – within 30 days of purchase – and we will replace the original mask at no charge (conditions apply).

Mask Codes	
Code	Description
SLN1SA	Nasal Mask Small
SLN1MA	Nasal Mask Medium
SLN1LA	Nasal Mask Large
SLN1WA	Nasal Mask Wide
SLN1SMLA	Nasal Mask Fit Pack/SML

Spare Parts	
Code	Description
400SLN115	Nasal Cushion Small
400SLN116	Nasal Cushion Medium
400SLN117	Nasal Cushion Large
400SLN118	Nasal Cushion Wide
400SLN131	Nasal Mask No H/G Small
400SLN132	Nasal Mask No H/G Medium
400SLN133	Nasal Mask No H/G Large
400SLN134	Nasal Mask No H/G Wide
400SL122	Headgear Spare
400SL141	Tube, Frame & Swivel Spare

References: 1. 42 out of 45 participants rated their experience with the F&P Solo Nasal mask as simple or very simple. Internal validation study conducted with 45 participants in USA in 2022.
 2. 40 out of 42 participants were able to set themselves up on the F&P Solo Nasal mask without assistance. Defined as orienting the mask on the head/face as intended and achieving an acceptable seal without human interaction. Internal self-setup trial conducted with 42 users, naïve to CPAP therapy, in New Zealand in 2022. F&P, AutoFit, AutoLock, F&P Solo and myMask are trademarks of Fisher & Paykel Healthcare Limited. For patent information, please see www.fphcare.com/ip. 630497 Rev A © 2023 Fisher & Paykel Healthcare Limited.